

HIGHLY ADJUSTABLE TO ELEVATE YOUR MOBILITY

The Quickie® 5^R is designed to deliver superior rigid wheelchair performance day after day. As you and your environment evolve, or as you refine your riding capabilities, the 5^R's frame and components adjust to elevate your mobility.





OPEN FRAME DESIGN

The Quickie 5th's open frame takes the ride quality and geometry that you love about the GT and improves upon it with a swept-back hydro-formed caster design for added style and a smooth ride.

The 5th offers two frame options. The **Classic** features a traditional frame length for maximum stability. The shortened **Performance** frame length offers maximum responsiveness.



NEW AXLE PLATE DESIGN

The 5th's swept-back axle plate design isn't just for looks. It is also incredibly quick to adjust and is available in two lengths to ensure you are getting the adjustments you need without adding extra weight.



Quickie Xtender Power Assist

Developed by Quickie and Yamaha, the Xtender features motorized rear wheels that increase the force you apply to the handrims by *four times*.



XTR Suspension

The XTR Suspension Frame helps you fearlessly conquer curbs, bumps, and tough terrain. Add Mountain wheels, Knobby tires, wide casters, and Frog Legs suspension caster forks for the ultimate off-road wheelchair experience.



Topolino Carbon Core Wheels

One of the lightest wheelsets on the market today at only 2.5 lbs. per pair. Topolino Carbon Core Wheels enhance your overall independence and mobility while making transport and storing much easier.



Custom-Colored JAY® J3™ Backs

Add a pop of matching or contrasting color to your 5th with a painted JAY J3 back. Choose from any of the standard or matte frame colors available on Quickie wheelchairs, and we will paint your J3 back at no additional charge.

THE 5^R
GIVES YOU
ACCESS TO THE
PREMIUM
OPTIONS
YOU WANT TO
ELEVATE
YOUR MOBILITY

COMPLETE OPTION MIX FOR FUNCTION & STYLE

The 5^R offers our greatest mix of rigid wheelchair options, including handrim options, HCPCS-coded tension-adjustable upholstery, single-sided forks, and colored components.



built4me

Our Built-4-Me Team is dedicated to creating product customizations beyond what is listed on our product order forms. Experts in the engineering and manufacturing of custom rehab products will listen to your specific requests and build a chair as individual as your fingerprint!

Customizations include:

- Custom paint
- Spinergy wheels with colored spokes
- Custom upholstery
(Shown on Quickie Q7 to right)

**Conforms with ANSI/RESNA WC-19 standards. Please refer to Quickie® & Zippie® WC-19 Transit Securement Supplement (MK-100072) for additional information. Supplement is available for download at www.SunriseMedical.com.*

5^R SPECS

HCPCS Code: K0005

Overall Weight: 21.5 lbs

Transport Weight: 12.9 lbs

Weight Capacity: 265 lbs

Seat Width: 12" to 20"

Seat Depth: 12" to 20"

Center of Gravity Adjustment:

-1.5" to 4"

Back Angle Adjustment:

Fore (-) 17° to Aft (+) 8°

Front Seat Height:

16" to 21"

Rear Seat Height:

13" to 20"

Transit Approved: Yes*

Warranty:

Limited/Lifetime



© 06.2015 Sunrise Medical (US) LLC
MK-100136 Rev. A

800.333.4000 www.SunriseMedical.com

