



Squiggles Seat

February 2023

Mark For:

Date:

Dealer Acct #:

Dealer:

Dealer Contact:

Dealer Address:

Dealer City: ST: ZIP:

Dealer Phone: Fax:

Confirmation Email:

Confirm Via: ☐ Fax ☐ Email

Submitting for:



Quote



Order

PO#:

ADDITIONAL SHIPPING INFORMATION

Ship To:

Attention:

Address:

Address:

Ship To City: ST: ZIP:

Ship To Phone: Fax:

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

Squiggles Seating System

Weight Capacity: 48.4 lbs

Seat Shell

Selection Required

	Item Code	Description	HCPC	Retail
☐	120-610S	Squiggles Seat Shell		
	120-610-1	Seat Base with Depth Adjustment	E2294	\$1,072
	120-610-2	Seat Back - Fixed	E2293	\$1,072
	120-610-3	Seat Back Reclining Mechanism	E1014	\$468
	120-610-4	Seat Base Hip Guides / Support	E0956	\$130

Seat Cover

Selection Required. Includes seat base and backrest cushion covers.

	Item Code	Description	HCPC	Retail
☐	120-710-03	Squiggles Seat Covers - Blue		N/C
☐	120-710-01	Squiggles Seat Covers - Green		N/C
☐	120-710-02	Squiggles Seat Covers - Orange		N/C
☐	120-710-04	Squiggles Seat Covers - Pink		N/C

Base

	Item Code	Description	HCPC	Retail
☐	120-800	Size 1 Hi Low Chassis (Foot Pedal Operated.)	E1399	\$1,919
	122-710	Hi Low Chassis Push Handle		\$279
☐	102-2753	Size 2 Hi Low Chassis (Foot Pedal Operated. Includes Push Handle)	E1399	\$2,199
	120-818	Adapter to Mount Squiggles Seat on Size 2 Hi Low Chassis		\$214
☐	194-600	Mobility Base Interface Plate (Allows mounting Squiggles Seat onto Sunrise Wheelchairs)	K0108	\$504

Foot Positioning

Sandals and Footplate Raisers require Squiggles Seat Adjustable Footrest

	Item Code	Description	HCPC	Retail
☐	120-639	Squiggles Seat Adjustable Footrest	K0040	\$237
☐	152-1600	Small Sandals w/ Straps (Pair)	E0951, E0952	\$267
☐	152-2600	Medium Sandals w/ Straps (Pair)	E0951, E0952	\$291
☐	120-816	Footplate Raiser 1" (Each)	E0951, E0952	\$90

Quantity: (1 to 3)

Harness

	Item Code	Description	HCPC	Retail
☐	120-844-03	Pelvic Harness Small - Blue	E0978	\$161
☐	120-844-01	Pelvic Harness Small - Green	E0978	\$161
☐	120-844-02	Pelvic Harness Small - Orange	E0978	\$161
☐	120-844-04	Pelvic Harness Small - Pink	E0978	\$161
☐	117-848-03	Pelvic Harness Medium - Blue	E0978	\$161
☐	117-848-01	Pelvic Harness Medium - Green	E0978	\$161
☐	117-848-02	Pelvic Harness Medium - Orange	E0978	\$161
☐	117-848-04	Pelvic Harness Medium - Pink	E0978	\$161
☐	120-789	Pelvic Harness Spacer Pads - Black (Pair)	K0108	\$90
☐	120-799	Small Trunk Harness	E0960	\$289
☐	120-757	Medium Trunk Harness	E0960	\$289

Lateral Support

Cover selection required with Lateral support hardware.

	Item Code	Description	HCPC	Retail
☐	120-855	Flip-Away Lateral Support Hardware (Pair)		\$337
☐	120-838-03	Lateral Support Covers - Blue (Pair)	E0956	\$154
☐	120-838-01	Lateral Support Covers - Green (Pair)	E0956	\$154
☐	120-838-02	Lateral Support Covers - Orange (Pair)	E0956	\$154
☐	120-838-04	Lateral Support Covers - Pink (Pair)	E0956	\$154
☐	120-656	Rigid Lateral Support Hardware (Pair)	E1028	\$205
☐	120-838-03	Lateral Support Covers - Blue (Pair)	E0956	\$154
☐	120-838-01	Lateral Support Covers - Green (Pair)	E0956	\$154
☐	120-838-02	Lateral Support Covers - Orange (Pair)	E0956	\$154
☐	120-838-04	Lateral Support Covers - Pink (Pair)	E0956	\$154
☐	120-954	Laterals w/ Protection Pads - Black (Pair)	E0956	\$433
☐	137-808-03	Chest Support - Blue <i>(Requires Laterals)</i>	E0960	\$189
☐	137-808-01	Chest Support - Green <i>(Requires Laterals)</i>	E0960	\$189
☐	137-808-02	Chest Support - Orange <i>(Requires Laterals)</i>	E0960	\$189
☐	137-808-04	Chest Support - Pink <i>(Requires Laterals)</i>	E0960	\$189

Tray

	Item Code	Description	HCPC	Retail
☐	120-672	Tray - Gray	E0950	\$401
☐	120-688	Extra High Tray +2" Adjustment	E0950	\$507
☐	120-270	Padded Tray Insert - Gray	K0108	\$332
☐	117-769	Grab Rail / Toy Bar for Tray	K0108	\$298

Pommel

	Item Code	Description	HCPC	Retail
☐	120-641	Pommel	E0957	\$113

Headrest

Cushion selection required with Headrest Support hardware

	Item Code	Description	HCPC	Retail
☐	120-859	Contoured Headrest Support Hardware	E1028	\$376
☐	137-800-03	Contoured Headrest Cushion - Blue	E0955	\$172
☐	137-800-01	Contoured Headrest Cushion - Green	E0955	\$172
☐	137-800-02	Contoured Headrest Cushion - Orange	E0955	\$172
☐	137-800-04	Contoured Headrest Cushion - Pink	E0955	\$172
☐	117-697	Flat Headrest Support Hardware	E1028	\$231
☐	117-850-03	Flat Headrest Cushion - Blue	E0955	\$118
☐	117-850-01	Flat Headrest Cushion - Green	E0955	\$118
☐	117-850-02	Flat Headrest Cushion - Orange	E0955	\$118
☐	117-850-04	Flat Headrest Cushion - Pink	E0955	\$118
☐	118-667	Lateral Bracket Assembly Hardware (Pair) <i>Requires Flat Headrest Support Hardware and Cushion</i>	E1028	\$212
☐	137-807-03	Flat Headrest Lateral Cushions - Blue (Pair)	E0956	\$192
☐	137-807-01	Flat Headrest Lateral Cushions - Green (Pair)	E0956	\$192
☐	137-807-02	Flat Headrest Lateral Cushions - Orange (Pair)	E0956	\$192
☐	137-807-04	Flat Headrest Lateral Cushions - Pink (Pair)	E0956	\$192
☐	117-745	Removable Head Support Bracket for 1/2" Box Section	E1028	\$71

Accessories

Cover selection required with Shoulder Support Hardware

	Item Code	Description	HCPC	Retail
☐	120-779-03	1" Back Spacer Pad and Flip-Away Lateral Supports - Blue	E0956, E1297	\$813
☐	120-779-01	1" Back Spacer Pad and Flip-Away Lateral Supports - Green	E0956, E1297	\$813
☐	120-779-02	1" Back Spacer Pad and Flip-Away Lateral Supports - Orange	E0956, E1297	\$813
☐	120-779-04	1" Back Spacer Pad and Flip-Away Lateral Supports - Pink	E0956, E1297	\$813
☐	120-702	Shoulder Support Hardware	K0108	\$275
☐	120-828-03	Shoulder Support Cover - Blue	K0108	\$156
☐	120-828-01	Shoulder Support Cover - Green	K0108	\$156
☐	120-828-02	Shoulder Support Cover - Orange	K0108	\$156
☐	120-828-04	Shoulder Support Cover - Pink	K0108	\$156

Technical Specifications		
--------------------------	--	--

Specification	Imperial	Metric
Age (Approximate)	1 - 5 years	1 - 5 years
Stature	29.5" - 43.3"	750 mm - 1100 mm
Maximum User Weight	48.4 lbs	22 kg
Seat Depth	7.5" - 11.4"	190 mm - 290 mm
Seat Width	6.3" - 10.2"	160 mm - 260 mm
Hip Width	6.3" - 8.7"	160 mm - 220 mm
Backrest Height	12.6" - 15"	320 mm - 380 mm
Backrest Height with Shoulder Support	16.5" - 18.9"	420 mm - 480 mm
Minimum Chest Width	6.3"	160 mm
Maximum Knee Width	3.5"	90 mm
Minimum Distance Between Pommel and Side of Seat	7.5"	190 mm
Minimum Seat to Footplate	4.8"	122.5 mm
Seat to Floor	11.4" - 26.4"	290 mm - 670 mm
Backrest Angle (Prone / Recline)	10° / 25°	10° / 25°
Footplate Angle (Plantar / Dorsiflexion)	10° / 10°	10° / 10°
Tray Size	18.1" x 22"	460 mm x 560 mm
Seat Weight	13.2 lbs	6 kgs
Tray Height	4.5" - 6.7"	115 mm x 170 mm
Squiggles Seat Hi-Lo Chassis Foot Operated (Excluding Push Bar)	Imperial	Metric
Weight	19.8 lbs	9 kgs
Length	25.2"	640 mm
Width	21.3"	540 mm
Minimum Dimensions for Storage	12.4"	315 mm



Sunrise Medical (US) LLC2842 Business Park Ave. • Fresno, CA 93727 • USA

Customer Service: 800-333-4000 Fax: 877-237-4214 Email: customseating@sunmed.com www.sunrisemedical.com

© 2/2023 MK-129856 Rev F

Page 4