


# ROHO<sup>®</sup> MOSAIC<sup>®</sup>

Comfort For You



**permobil**



## Experience the Value

The re-engineered, cost-effective ROHO MOSAIC cushion is now even better. With a more comfortable fit, you'll enjoy more sitting time while still protecting your skin and soft tissues.

## Designed For You

By always matching your shape and following your movements throughout the day, the MOSAIC cushion is designed for you. Now with more anatomically-designed cell heights and a more effective, controlled air flow, you can sit longer, with less pain, knowing you have proven ROHO air-cell-based technology supporting you.

## Proven Protection

The innovative design of the MOSAIC cushion allows you to sink comfortably into it, shaping to the unique contours of your body. Available in an expanded range of sizes to better fit your seating surface, MOSAIC provides basic preventative skin and soft tissue protection at an affordable price.

## Specifications

Cover Type	2-Way Stretch Cover Heavy Duty Cover (optional)
Sizes	16" x 16" (41cm x 41cm) 16" x 18" (41cm x 46cm) 17" x 17" (44cm x 44cm) 18" x 16" (46cm x 41cm) 18" x 18" (46cm x 46cm) 20" x 18" (51cm x 46cm)
Product Weight	1.35lbs (0.6kg)*
Weight Limit	315lbs (143kg)
Cushion Warranty	12-Month Limited Warranty
Cover Warranty	6-Month Limited Warranty
U.S. Medicare Code	E2601 (pending approval)

\*Product weight varies depending on cushion size.

Since U.S. Medicare coding and coverage criteria is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process. ©2017 Permobil. The following are registered trademarks of ROHO Inc: MOSAIC®, ROHO®. Permobil has a policy of continual product improvement and reserves the right to amend this document. The current version of this document is available at permobilus.com.