

Invacare® MVP® and MVP® jr. Folding Wheelchair



Yes, you can.®

MVP® Wheelchair



Anything else is a compromise.™

Invacare® MVP® and MVP® jr. Folding Wheelchair

Versatile, growable, and designed with performance in mind, the Invacare® MVP® and MVP® jr. are the premier folding chairs in the Invacare Custom Manual wheelchair line. With a refined sense of style and three rear frame options, these chairs are easy to match the consumers needs – responsiveness (curved frame), stability (square frame), or growth and positioning (extended frame). The three-piece modular frame design easily accommodates growth by replacing the front frame section. Backed by the quality and reliability that is built into every Invacare product, the MVP and MVP jr. give consumers the durability and performance they need and set the standard for every folding chair in the industry. Think MVP when consumers want maximum versatility and performance with the convenience of a folder.



The MVP with curved rear frame, tapered swingaway front riggings, 24" high performance spoke wheels, 5" aluminum casters, and Teal frame.

Features



Adjustable angle back option on the square frame provides four angle positions for postural requirements: 90°, 95°, 100°, 105°.



Front head tubes offer a cleaner look and easy adjustment to the front caster to compensate for rear seat-to-floor height changes.



Hideaway Undermount Wheel Locks provide great action and remain out of the way while propelling and transferring. Available with fixed front frame only.

Features

- Lightweight aluminum folding frame
- Adjustable angle back option on the square frame
- Adjustable angle footplate and clamp
- New front head tubes style
- Hideaway Undermount Wheel Locks available
- Improved wheels include black rim and black hub
- Three rear-frame options – Curved (responsive), Square (stability) and Extended (growth and positioning)
- Standard axle plate on the Square and Extended frames allows center-of-gravity adjustment under the back canes for maximum stability
- Virtually unlimited seat-to-floor height/wheel combinations
- Three piece modular frame easily accommodates growth by replacing the front frame
- Heavy-duty package increases weight limit from 250 lb. to 350 lb. (MVP only)
- Less-Seat-Rail option provides lower seat-to-floor heights
- Camber choices from 0° to 9°
- Choice of fixed or swing-away front end
- Thumb Saver Wheel Locks are standard
- Tapered front riggings available



Available in frame sizes as small as 12" x 12", the MVP jr. is a great solution for your pediatric needs.

Rear-Frame Options

The MVP® and MVP® jr. are available with three rear frame options – curved, square and extended frames.



The curved rear frame option offers great maneuverability and quick-response performance, while providing a sporty appearance and the shortest overall length.



The square rear frame option performs with outstanding stability. The axle plate allows the rear wheels to be placed directly below the back canes.



The extended rear frame option allows adjustment for the positioning needs of both adolescents and individuals with special seating requirements.

Adjustable Angle Footplate and Clamp



Adjustable angle footplates offer softer edges.



The footplate clamp is double-sided for increased durability and superior performance.

MVP® and MVP® jr. Extended Frame with Less-Seat-Rail (LSR)



To accommodate the consumer's growth, the frame is depth adjustable and features an adjustable angle back.



The Less-Seat-Rail option allows easy placement of seating systems.



Yes, you can.®

Specifications

	MVP	MVP jr.
Seat Width	15" - 22"	12" - 16"
Seat Depth	9" - 22"	9" - 19"
Minimum overall width*	Seat width plus 10"	Seat width plus 10"
Adjustable angle back**	90°, 95°, 100°, 105°	90°, 95°, 100°, 105°
Adjustable height back	10" - 18"	10" - 18"
Camber angles	0°, 3°, 6°, 9°	0°, 3°, 6°, 9°
Tubing outside diameter O. D.	1" aluminum	1" aluminum
Weight of chair	17 lb. (16"x16" without wheels and accessories)	
Weight capacity	250 lb.	150 lb.
Heavy-duty package	350 lb.	n/a
Warranty on frame and crossbraces	Limited Lifetime	Limited Lifetime

*Measurement taken from outside of handrims with 3° camber.
 **Not available with Heavy-duty package.

Rear Seat-to-Floor

(range adjusts in 7/8" increments)	Curved	Square/Extended
16" wheel	14.75" - 16.5"	14.75" - 17.75"
18" wheel	14.75" - 17.5"	14.75" - 18.75"
20" wheel	14.75" - 19"	14.75" - 20"
22" wheel	15.25" - 20"	15.75" - 21"
24" wheel	15.75" - 20.5"	16.5" - 21.25"
26" wheel*	16.75" - 21.25"	17.5" - 21.25"

*Not available on MVP jr.

Front Seat-to-Floor

	Swingaway		Fixed	
	Min.	Max.	Min.	Max.
3" caster	15.25"	18.75"	14.75"	18"
5" caster	16"	19.75"	15.5"	19"
6" caster	17"	20.25"	16.5"	19.25"
6" caster	18"	20.25"	17.5"	19.75"
8" caster	18.75"	21.25"	18.25"	20.5"

Knee-to-heel*

Front Rigging Type	Min.	Max.**
70° tapered	15.75"	19.75"
70°	13.50"	18.50"
70° MFX	8.75"	14.75"
Elevating	15.25"	21.25"
90°	5.25"	10.25"
60°	14.75"	18.50"
60° MFX	11"	15"

*Ranges adjust in 1" increments.
 **Measurement taken from top of seat rail to back of footrest.

Arm Heights*

Swingaway padded	8" - 10" (in 1" increments)
Cantilever (bolt-on)	5" - 10" (in 1" increments)
Adjustable height Dual point	10" - 14" (in 1" increments) (conventional or spacesaver)
T-arm low range	6" - 10" (in 1/2" increments) (full and desk range)
T-arm high range	8" - 14" (in 1" increments) (full and desk range)

*Measurement taken from top of seat rail to top of armrest.

Invacare Corporation
 www.invacare.com

USA

One Invacare Way
 Elyria, Ohio
 44035-4190
 (800) 333-6900

Canada

570 Matheson Blvd. E.
 Unit 8
 Mississauga, Ontario
 L4Z 4G4 Canada
 (800) 668-5324

Invacare, the Blue Medallion Design, Yes, you can. and MVP are registered trademarks and Anything else is a compromise is a trademark of Invacare Corporation. Specifications are subject to change without notification.

© 2007
 Invacare Corporation
 Form No. 01-350
 Rev. 03/07