

MedRISER™ Series



600Lb Patient Lift

- **Electric**— powered electric lift (height and base) will support up to 600lbs and can be used in a variety of lifting situations.
- **6 point spreader** — added width helps enhance comfort to end user while transferring. Pair with a variety of slings for transfers, rehab/gait training, turning and repositioning.
- **Lift Height** — built to have an extremely long lift arm achieving a high lift height which allows it to be utilized for a variety of lifting needs and care settings.
- **GaitTraining**—achieving an extremely high lift height will allow for gait training, ambulation, and various other rehabilitation exercises.
- **OptionalScale**— have the ability to weight patients anywhere and at anytime with a bolt on scale.
- **Ergonomics** -enhanced thicker push handles with bottom "kick" support bar for additional leverage and staff safety.
- **Made in the USA** — an affordable quality lift solution with a commitment to keeping inventory and parts on hand.(no customs or long lead times)

Med-Riser Lift SPECIFICATIONS

Safe Working Load	600 LBS.
Front Casters.....	4"
Rear Casters.....	4"
Max Height.....	83" (pivot pin)
Low Height.....	19.625" (spreader bar)
Legs Open.....	44" (center of caster)
Legs Closed.....	24" (center of caster)
Overall Length.....	50.2"
Warranty	3 Year Limited (Parts Only) 1 Year (Casters, Hand Control) 90 Days (Battery)



Accessory Options



600lb Scale


PROUDLY MADE IN THE U.S.A.

(877) 867-7365
MED-MIZER.COM



Gait Training Bars

Sling Options



HighBack Sling (mesh)



HighBack Sling (padded)



Hygiene Sling



Ambulation Sling

MED MIZER®

is at the forefront of innovating medical technology for post acute care. We are focused on partnering with you to help address healthcare challenges in the following key areas—




FALLS
PREVENTION




BARIATRIC
CARE



MOBILITY



REHABILITATION



WOUND CARE