

# MEDRISER™

## FlexSeries

Padded "Flex" Knee/Shin Plate improves patient comfort

### 500Lb Sit to Stand Lift



- Electric—powered electric lift (height and base) will support up to 500lbs and can be used in various transfer situations.
- "Flex" Kneeadapted from our specialty seating line the pressure relief strapping and pad aids in providing superior comfort to leg and knee during transfer.
- Optional scale to easily weigh patients.
- Multiple Sling Bar Options — Easily change the sling bar to help enhance resident comfort and support based on clinical preference.
- Removable "Padded" Foot Rest — enhanced comfort with soft foam removable foot rest which aids in safety due to textured design.
- Dual controls - controls on the pendant and control box allow the lift to be operable at any time.
- Ergonomics— enhanced thicker push handles with bottom "kick" support bar for additional leverage and staff safety.
- Made in the USA — an affordable quality lift solution committed to keeping inventory and parts on hand. (no customs or long lead times).

## Med-Riser Lift SPECIFICATIONS

|                         |                                                                                    |
|-------------------------|------------------------------------------------------------------------------------|
| Safe Working Load ..... | 500 LBS.                                                                           |
| Front Casters.....      | 4"                                                                                 |
| Rear Casters.....       | 4"                                                                                 |
| Max Height.....         | 60.5"                                                                              |
| Low Height.....         | 31" (spreader bar)                                                                 |
| Legs Open.....          | 37.25" (center of caster)                                                          |
| Legs Closed.....        | 22.5" (center of caster)                                                           |
| Overall Length.....     | 40.5"                                                                              |
| Warranty .....          | 3 Year Limited (Parts Only)<br>1 Year (Casters, Hand Control) 90<br>Days (Battery) |



Shown with  
optional  
comfort-wide  
slingbar

## Accessories



Calf Strap



Comfort-Wide Slingbar



Scale Option



Handlebar (standard only)

PROUDLY MADE IN THE U.S.A.

(877)867-7365

[MED-MIZER.COM](http://MED-MIZER.COM)



Seat Strap Mount



Seat Strap

**MED MIZER®**

is at the forefront of innovating  
medical technology for health care.  
We are focused on partnering with you  
to help address healthcare challenges  
in the

## MedRiser Sling




FALLS  
PREVENTION



BARIATRIC  
CARE



MOBILITY



REHABILITATION



WOUND CARE