



MDS806300DLX

SPECIFICATIONS

Weight Capacity	300 lbs	Warranty on Upholstery Parts & Components	6 months
Seat Width	16", 18"	Footplates	composite
Seat Depth	16"	Dual Axle	yes + transport
Warranty on Frame & Crossbar	limited lifetime	HCPCS Code	K0001, K0002

Item #	Description
MDS806150NDLX	16" Permanent full-length arms, swing-away detachable footrests
MDS806200NDLX	16" Permanent full-length arms, swing-away detachable elevating legrests
MDS806250NDLX	16" Removable desk-length arms, swing-away detachable footrests
MDS806300NDLX	16" Removable desk-length arms, swing-away detachable elevating legrests
MDS806150DLX	18" Permanent full-length arms, swing-away detachable footrests
MDS806200DLX	18" Permanent full-length arms, swing-away detachable elevating legrests
MDS806250DLX	18" Removable desk-length arms, swing-away detachable footrests
MDS806300DLX	18" Removable desk-length arms, swing-away detachable elevating legrests
COLOR: Black	

The Most Features in Any Standard Wheelchair

Exclusively from Medline:

optional Quick Release Transport Wheel Accessory. (WCAWHEELKIT)



Quick release axles:

push-button design allows tool-free removal of the wheels.



Scissor frame style:

reduces weight up to 5 lbs over slide tube style frames.




Three axle positions:

adjustable for Standard, Hemi and Transport chair.



Strong, Cast Aluminum Seat Cradles:

not plastic as typically seen on chairs.



Other Features

- Nylon upholstery is breathable and light weight.
- Threaded seat screw inserts reduce stripping from upholstery changes.
- Padded nylon armrests.
- No spacers on front casters, making wheel changes easier.
- Carbon steel frame
- Sealed bearings in fork housing, rear wheels and front wheels keep dirt out and are maintenance free.
- Mag-style composite wheels are extremely durable and easy to clean.
- Heavy-duty fork and welded fork stems minimize wobbling.
- Quick-release wheel lock allows easy relocation of wheel lock for setting in transport mode.