



Ultra-Lightweight Folding Electric Wheelchair

The NEW Modo Breeze is the perfect power chair for people who need first-time or occasional mobility support.

Only
\$ 1,999
MAP



Ultra-light & compact



Sleek carbon fiber frame



Color LCD screen



Reinforced footrest and suspension

LED Headlight + USB Charging Port

2" Contoured Pincore Cushion

Calf Strap

Airline-Friendly 10Ah Li-Ion Battery

Reinforced Footrest Support

Front Suspension Springs



Rear Storage Pocket

Adjustable Back-Support Tension Straps

Adjustable Lap Belt

Rear Wheel Fenders

Anti-Tip Safety Bars

Spare Battery Bracket

HIGHLIGHTED FEATURES

- Unit weighs 28.6 lbs. without battery and cushions
- 300 lbs. weight capacity
- Reinforced footrest support adds strength and stability for long-term durability
- 2" contoured pincore seat cushion enhances comfort and air circulation for all-day comfort
- Airline-friendly 10Ah Li-ion battery
- Programmable quick disconnect joystick features a color LCD screen with built-in headlight + USB charging port
- Collapsible push handles and flip-up footrest for easy storage
- X-frame design allows the unit to fold nearly flat in a vertical position
- Flip-up armrests and footplates for easy front and side entrance

SPECIFICATIONS

Frame:	Carbon Fiber
Battery:	10Ah / 24V / 240 Wh, lithium-ion battery
Battery weight:	3.9 lbs. (1.8 kg)
Driving range:	Up to 9.5 miles
Motor:	180 W (x2); Brushless DC
Brake:	Electromagnetic Brake
Seat cushion:	2" thick
Max speed:	3.7 mph
Front wheel:	6.5" solid tire
Rear wheel:	11" solid tire
Size (unfolded):	38.97" x 20.86" x 34.25" (99*53*87 cm)
Size (folded):	29.92" x 8.66" x 27.16" (76*22*69 cm)
Max incline:	12 degrees



WARRANTY:

Frame: 5 years/ Seat and back cushion: 2 years/ Electrical and mechanical components, batteries, and tires: 1 year. For complete warranty terms and conditions, visit whill.inc/us/warranty-policy