



Night Rider® HD

BARIATRIC SLEEP SYSTEM

With the Night Rider® HD Adjustable Bed System, ***Dreams Made Easy™*** isn't just a saying — it's how you sleep every night. This bariatric bed is designed for the way you live, giving you the **comfort, independence, and support** you deserve.

With a wide range of bed positioning options — including **full tilt, Trendelenburg, and Reverse Trendelenburg** — you can find the angle that best fits your comfort while also addressing specific health needs.



Experience a bed that adapts to you — delivering **comfort, independence, and restorative sleep** so you wake refreshed and ready for the day.

Dreams, Made Easy™



Comfort Designed With Purpose

The Night Rider® HD combines the care you need with the luxury you want helping you stay active while aging with more comfort and independence.

- ✔ **Effortless adjustment:**
Head, foot, height (Hi-Low) to help you find that perfect bed positioning
- ✔ **Available in:** Full & Queen sizes
- ✔ **Advanced Positioning:**
Help manage health conditions both at home and in the facility
 - Tilt into Trendelenburg to support circulation
 - Reverse Trendelenburg to make breathing easier or help if you suffer from conditions such as acid reflux
- ✔ **Customized height: Super low height of 12.75":**
(for bed transfers) up to 24.75" (ideal for caregivers to provide care) making it easier to get in and out of bed
- ✔ **Seamless Sleep Deck:** 3DL material that makes it water-resistant, scratch and easy to clean – built to last and designed for everyday life
- ✔ **Sturdy design:** 750 lb (340 kg) weight capacity
- ✔ **Rolling casters:** let you move the bed around when cleaning or adjusting its location in the bedroom
- ✔ **Soothing massage:** helps you unwind and reduce body tension
- ✔ **Emergency power-down:** function for added peace of mind
- ✔ **Backlit wireless remote:** makes night time adjustments easy



Optional Add-Ons:

- ✔ **PressureGuard EasyAir or APM2 mattresses:**
Offer superior pressure relief and support
- ✔ **Support Rails:** Added assistance for getting in and out of bed
- ✔ **Headboard:** Optional laminate headboard to go against the wall
- ✔ **Bamboo Rail Covers:** Soft, stylish protection for extra comfort
- ✔ **Battery Backup:** Keeps you powered even during outages

GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 678.291.3190

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.transfermaster.com



© 2025, GF Health Products, Inc. All Rights Reserved.
TransferMaster and Supernal are registered
trademarks of GF Health Products, Inc.



Comfortable In Every Setting

Created with your comfort in mind, the Night Rider® HD bed brings you the benefits of an **adjustable bariatric hospital bed alternative with advanced positioning** in a beautifully discreet style that blends seamlessly into your home.

Designed for the Way You Live

Backed by more than 100 years of Graham-Field's Made in USA expertise, our sleep solutions are designed for the way you live — delivering **Dreams, Made Easy™** for every stage of active aging.

Night Rider® HD

BARIATRIC SLEEP SYSTEM



For more information visit us online at
www.transfermaster.com

TransferMaster
by graham-field