

PIGLEO-T



PEDIATRIC
COLLAPSIBLE
LIGHTWEIGHT
TILT-IN SPACE STROLLER

Nissin Wheelchairs USA
3551 Voyager St., Unit 102
Torrance, CA 90503
(877) NISSIN7 (310) 370-8484 www.nissinusa.net

 **NISSIN**
WHEELCHAIRS
...YOUR CHAIR IS HERE!

PIGLEO-T



**LIGHTWEIGHT
COLLAPSIBLE
PEDITRIC
TILT-IN SPACE STROLLER**

SPECIFICATIONS & OPTIONS

Standard

- Adjustable Height Stroller Handles
- 4-Point Chest Harness
- Angle Adjustable Back
- Anti-Tippers
- Foot Straps

Two Color Choices



Options

Ventilator Box




Sun Visor



Head Pillow



Seat Wedge



Hip Guides



Pigleo Patch



Available in two sizes Standard and Large

Specifications

Specifications	Open
Weight Capacity	80 lbs.(Large 100lbs.)
Adjustable Handle Height	39.8" - 49.6"
Frame Height	33"
Frame Length	14.2"
Back Width	15.75"
Tilt Range	-5° - 45°
Seat Width	13" (Large 15")
Seat Depth	9.8" - 12.2"
Seat Height	18.5"
Seat Back Height	24.4"
Back Angle	90°, 95°, 100°
Leg Length	7.1" - 11.8"
Product Weight	20.5 lbs.
Folded Dimensions(both)	H29" x W 19.7" x L17.7"

