

# TiLite 2GX Order Form

Effective March 5, 2018 (New Items in Red)

FOLDING FRAME

Page 1 of 6

DEALER INFORMATION		RTS INFORMATION	
Dealer:	Account #:	RTS:	Phone:
PO:	Ordered By:	Email:	
Email:		<b>SHIP TO INFORMATION</b>	
Mark For:		Name:	Ship Via:
TiLite Sales Rep:		Address:	
<b>TiLite Information:</b>	Promo Code:	City:	State: Zip:
Ship Date:	Order #:	Phone:	Fax:

**Please Email or Fax to the TiLite Customer Service Team at orders@tilite.com or 866.586.2413 / 509.586.2413**

## FRAME STYLE

- G2FS1 ☐ **2GX2 Titanium Frame** (Full Adjustability. 265 Pound Weight Limit. Light Weight, High Strength and Vibration Dampening. Not Available with Frame Transit Tie-Down.).....**\$4,265**
- XFS1 ☐ **Swing Away** .....**STD** K0005
- XFS3 ☐ **Hemi Height Frame** (Available with a Front Seat Height of 13.5" - 20" (34-51cm). Not Available with Depth Adjustable Back/Frame, 12" or 13" Seat Depth or Heavy Duty Frame.) .....**NCO**
- XFS2 ☐ **Fixed Front** .....**NCO** K0005
- XHDR1 ☐ **Heavy Duty Frame** (350 Pound Weight Limit. Includes Tension Adjustable by Straps Back Upholstery. Available with 18" to 22" Seat Widths and 16" to 20" Seat Depths. Not Available with 90° front end).....**\$435**
- XPWR1 ☐ **Power Adaptable Frame** (No charge when 20" Rear Seat Width **A** and 20" Seat Depth **B** are selected. Requires LiteSpeed Billet Aluminum Front Wheels and 0 degree Camber. Not available with Single-Sided Forks).....**\$340**
- XRRF1 ☐ **Reinforced Frame** (No charge when 20" Rear Seat Width and 20" Seat Depth are selected.) .....**\$340**
- XMAXSD1 ☐ **Permobil Smart Drive MX2+** .....**\$7,131** E0986
- XMAXSD2 PushTracker Size ☐ Small ☐ Medium ☐ Large.....**STD**
- XMAXSD4 ☐ Carrying Bag .....**\$78**
- XMAXSD5 ☐ Thumb Throttle Buttons .....**\$437**

## TRANSIT OPTIONS

- XWTB ☐ **Transport Bracket Option** (Consists of four Tie-Down Brackets. Unoccupied Use Only.) (Not Available when Swing Away frame and Heavy Duty Frame are both selected. Not available with Swing Away Tubular Armrest or Swing Away -Flip Back Tubular Armrests. Requires a minimum seat depth of 14" or 15" seat Depth with Hemi Height Frame.).....**\$350**

## FRAME FINISHES

### Natural Finishes

- XFF1 ☐ **Satin** (Available with Titanium Frame only.).....**STD**
- XSCBT1 ☐ **Satin Finish Buffing Pads**.....**\$20**
- XFF2 ☐ **Polished** (Available with Titanium Frame only.).....**\$795**
- XFF3 ☐ **Tattoos** (Tattoos on front frame tubes. Titanium frame finished in Satin . Aluminum frame tattooed in blasted aluminum with clear coat. Tattoos will be trimmed and adjusted to fit your frame.).....**\$335**

XFF45 ☐ **Need For Speed**



XFF46 ☐ **Neon Dreams**



XFF44 ☐ **Flower Garden**



XFF43 ☐ **Mosaic**



### High Gloss Paint

### Metallic Paint

### Pearlescent Paint

### Matte Paint

### Candy Paint

### Matte Textured Paint

XFF8 ☐ **Super Wet Black**

XFF14 ☐ **Ocean Blue Metallic**

XFF17 ☐ **Midnight Blue Pearl**

XFF19 ☐ **Black Lava**

XFF22 ☐ **Ruby Red**

XFF40 ☐ **Hyper Blue**

XFF10 ☐ **TiLite Red**

XFF15 ☐ **Tangerine Metallic**

XFF18 ☐ **White Pearl**

XFF27 ☐ **Titanium Grey**

XFF23 ☐ **Azure Blue**

XFF41 ☐ **Flame Orange**

XFF29 ☐ **Neon Pink**

XFF30 ☐ **Acid Green Metallic**

XFF25 ☐ **Electric Plum**

XFF42 ☐ **OD Green**.....**STD**

## SEAT WIDTH **A** Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

- XSW1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19".....**STD**
- XSW2 ☐ 20" ☐ 21" ☐ 22" (21" and 22" Seat Widths **A** are only available with Heavy Duty Frame.).....**\$495** E2201

## SEAT DEPTH **B** Measure from 1/2" in front of back post to the front edge of seat sling. (12"-14" Seat Depths **B** not available with 25" (559) and 26" (560)

Rear Wheels.)

- XSD1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19".....**STD**
- XSD2 ☐ 20".....**\$525** E2203
- XSD3 ☐ **Custom Frame Depth:** \_\_\_\_\_ " Frame Depth Longer Than Seat Depth (Use 1" increments with a maximum of 2". A Custom Frame Depth with a maximum of 20" is \$475.).....**NCO**



tilite.com  
2701 West Court Street | Pasco, Washington 99301  
509.586.6117 Tel | 800.545.2266 Toll-Free  
509.586.2413 Fax | 866.586.2413 Toll-Free Fax

TiLite® OF0031\_2GX\_2018RevA



# TiLite 2GX Order Form

Effective March 5, 2018 (New Items in Red)

FOLDING FRAME

Page 2 of 6

## FRAME MEASUREMENTS



**FRONT SEAT HEIGHT ③** Measure from floor to top of seat tube at front edge of seat sling. Adjustable in 1/2" increments. (13.5" - 14" Front Seat Height ③ requires Hemi Height Frame Style. 20.5" - 21.5" Front Seat Height ③ not available with Hemi Height Frame Style.)

### FRONT SEAT HEIGHT RANGES

Front Wheel Size	Front Seat Heights ③	Hemi Front Seat Height ③
3" Front Wheel	14.5" - 18.5"	13.5" - 17.5"
3" Front Wheel w/Glide™	16.5" - 17.5"	15" - 17.0"
4" Front Wheel	15.5" - 19"	14.5" - 18.0"
4" Front Wheel w/Glide™	17.0" - 18."	15.5" - 17"
5" Front Wheel	16.5" - 19.5"	16" - 18.5"
5" Front Wheel w/Glide™	17.5" - 18.5"	16.5" - 18.0"
6" Front Wheel	17.5" - 20"	16.5" - 19.0"
6" Front Wheel w/Glide™	19"	17.5" - 18"
8" Front Wheel	19.5" - 21"	18.5" - 19.5"

XFSH1	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"
	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"
	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"
	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"	<input type="checkbox"/> 21"
				.....STD

**REAR SEAT HEIGHT ④** Measure from floor to top of seat tube at back post. Adjustable in 3/8" increments.

### REAR SEAT HEIGHT RANGES

Rear Wheel Size	Rear Seat Height ④
22"	13" - 18"
24"	14" - 19"
25"	14" - 19"
26"	15" - 20"

XRSH1	<input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"
	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"
	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"
	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"
	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20".....STD

**SEAT BACK HEIGHT ⑤/SEAT BACK/PUSH HANDLES** Measure from top of back post to top of seat tube at rear of frame.

### SEAT BACK TYPE

- XSB13 ☐ Fixed Angle, Adjustable-Height Seat Back.....STD
- XLB1 ☐ Folding Stabilizer Bar (Required with Transit Tie-Down Option unless Depth Adjustable Backrest is selected. Includes Tension Adjustable by Strap Back Upholstery.).....\$390
- XSB14 ☐ Fixed Angle, Adjustable Height Seat Back w/8-Degree Bend. (Not available with Heavy Duty Frame.).....NCO
- XSB18 ☐ Folding Adjustable-Angle and Height Seat Back (Not available with Heavy Duty Frame or Ride Designs® Corbac™.).....\$230
- XSB16 ☐ Depth Adjustable Back/Frame (Not Available with Hemi Height Frame, Heavy Duty Frame, 19" or 20" Seat Depths ⑤, Swing Away Tubular or Swing Away Flip-Back Armrest. Seat Depths ⑤ of 12"-17" allow 3" of rearward adjustment. Seat Depth ⑤ of 18" allow 2" of rearward adjustment. Requires Tension Adjustable Seat Upholstery).....NCO

### PUSH HANDLES

- XSB6 ☐ Integrated Push Handles.....NCO
- XSB7 ☐ Integrated Fold-Down Push Handles (Requires Tension Adjustable by Straps Back Upholstery. Not available with Heavy Duty Frame.).....\$115
- XSB9 ☐ Bolt-On Height Adjustable Push Handles (Requires Tension Adjustable by Straps Back Upholstery. Not available with Heavy Duty Frame, Folding Adjustable-Angle and Height Seat Back, Swing Away Tubular Armrests, Swing Away - Flip Back Tubular Armrests or Ride Designs® Corbac™.).....\$335

XSB20/XSB14 **FIXED ANGLE ADJUSTABLE SEAT BACK HEIGHT ⑤** (11"- 15" Fixed Seat Back Height not available with Heavy Duty Frame, Swing Away Tubular Arms 12"-15", Bolt-On Height Adjustable Push Handles or Folding Stabilizer Bar.)

Push Handles (Selected on Pg. 2)	Short	Tall	Set at: (Use 1" increments)
None	<input type="checkbox"/> 11"- 15"	<input type="checkbox"/> 16"- 20"	Set at: _____
Integrated	<input type="checkbox"/> 11"- 15"	<input type="checkbox"/> 16"- 20"	Set at: _____
Integrated Fold-Down	<input type="checkbox"/> 11"- 15"	<input type="checkbox"/> 16"- 20"	Set at: _____
Bolt-On Height Adj.	-----	<input type="checkbox"/> 16"- 20"	Set at: _____

XSB18 **FOLDING ADJUSTABLE- ANGLE AND HEIGHT SEAT BACK ⑤** (Not available with a Folding Stabilizer Bar or Heavy Duty Frame)

Push Handles (Selected on Pg. 2)	Short	Medium	Tall	Set at: (Use 1/2" increments)
None	<input type="checkbox"/> 11"- 13"	<input type="checkbox"/> 13.5"- 16.5"	<input type="checkbox"/> 16.5"- 19.5"	Set at: _____
Integrated	<input type="checkbox"/> 12.5"- 14.5"	<input type="checkbox"/> 13.5"- 16.5"	<input type="checkbox"/> 17"- 19.5"	Set at: _____
Integrated Fold-Down	<input type="checkbox"/> 15"- 17"	<input type="checkbox"/> 16"- 19"	<input type="checkbox"/> 18.5"- 20.5"	Set at: _____

# TiLite 2GX Order Form

Effective March 5, 2018 (New Items in Red)

FOLDING FRAME

Page 3 of 6

## FRAME MEASUREMENTS

### CAMBER ❶

XCBR8 ☐ 0° ☐ 2° ☐ 4° ☐ 6°.....STD

**FRONT ANGLE ❷** Front Angle ❷ is measured from inside of front frame tube to top of bottom frame tube.

#### Swing Away Front End

XFA3 ☐ 70° ☐ 80° ☐ 85° (Not available with 6" Front Wheels.) ☐ 90° (Requires aluminum footplate. Not available with Heavy Duty Frame, 5" or 6" Front Wheels.).....STD

XFA5 Elevating Leg Rest ☐ L **OR** ☐ R (When choosing a left or right Elevating Leg Rest, please select the opposite Front Angle (XFA3), above. Not available with 14" or narrower Seat Width ❶ or Neoprene Impact Guards.).....\$170 E0990

XFA7 ☐ Elevating Leg Rest Pair (Not available with Heavy Duty Frame, 14" or narrower Seat Width A or Neoprene Impact Guards.).....\$335 E0990

XFA0 ☐ Omit Hangers.....NCO

#### Fixed Front End

XFA6 ☐ 70° ☐ 80° ☐ 85° (85° Not available with 6" Front Wheels).....STD

**SEAT TO FOOTREST ❸/FOOTRESTS** Measure from front edge of seat sling to top rear of footrest.

If your Front Angle ❷ is **70°, 80°, OR 85°** select one Seat to Footrest ❸ length and one Footrest type:

**Standard Footrests** Measurement must be at least 2" shorter than Front Seat Height ❸. 70° Fixed Front has a minimum seat to footrest of 14.5"

☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20"

XFTR5 ☐ Flip-up Footrest.....STD

XFTR6 1-Piece Flip-up Footrest - Pivot on ☐ L **OR** ☐ R (Not available with Heel Loops.).....\$175 K0045

XFTR7 ☐ 1-Piece, Over-Center Folding Footrest (Fixed Front Only. Not available with Heel Loops.).....\$175 K0045

XFTR9 ☐ Depth Adjustable Flip-up Footrest.....\$110 K0045

XFTR10 ☐ Depth & Angle Adjustable Flip-up Footrest (Available for 15"- 20" Seat to Footrest ❸ only.).....\$270 K0040, K0045

XFTR11 ☐ Angle Adjustable Flip-up Footrest (Available for 15"- 20" Seat to Footrest ❸ only.).....\$170 K0040, K0045

**High Mount Footrests** (Not Available with Luggage Carriers, or Impact Guards.)

☐ 7" ☐ 7.5" ☐ 8" ☐ 8.5" ☐ 9" ☐ 9.5" ☐ 10" ☐ 10.5" ☐ 11" ☐ 11.5" ☐ 12" ☐ 12.5"

☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5"

XFTR20 ☐ High Mount Flip-up.....\$190 K0037

XFTR15 ☐ Angle Adjustable High Mount Flip-up Footrest.....\$325 K0040, K0037

XFTR16 ☐ Depth & Angle Adjustable High Mount Flip-up Footrest.....\$375 K0040, K0037

#### Elevating Leg Rest

☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" ☐ 20" ☐ 21"

XFTR17 ☐ Depth Adjustable Flip-up Footrest.....STD

If your Front Angle ❷ is **90°**, select one Seat to Footrest ❸ length and one Footrest type:

☐ 7" ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13"

XFTR5 ☐ Flip-up Footrest.....STD

XFTR9 ☐ Depth Adjustable Flip-up Footrest.....\$110 K0037

XFTR10 ☐ Depth & Angle Adjustable Flip-up Footrest.....\$270 K0040, K0037

XHLP1 ☐ **Heel Loops** (Not available with 13" or narrower Seat Width ❶ or 1-Piece Footrests.).....\$50 E0951

## CUSTOM CONFIGURATION OPTIONS

**FRONT WHEELS** (See Front Seat Height Ranges Chart on Page 2 for Front Seat Height ❸ limitations.)

XFW24 ☐ 3"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub **OR** ☐ Black Hub.....\$190

XFW14 ☐ 4"x.75" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or Glide™ Single-Sided Suspension Fork.).....NCO

XFW3 ☐ 4"x.75" Light-Up Micro-Caster (Not available with Power Adaptable Frame, Heavy Duty Frame or Glide™ Single-Sided Suspension Fork.).....\$80

XFW25 ☐ 4"x.75" LiteSpeed Billet Aluminum Wheel w/Poly Tire ☐ Silver Hub **OR** ☐ Black Hub.....\$150

XFW26 ☐ 4"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Assist, Power Add-on System or Heavy Duty Frame.).....\$95

XFW27 ☐ 4"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub **OR** ☐ Black Hub.....\$200

XFW6 ☐ 5"x1" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or with 90° Front Angle ❷ or Glide™ Single-Sided Suspension Fork.).....NCO

XFW28 ☐ 5"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or 90° Front Angle ❷.).....\$95

XFW29 ☐ 5"x.75" LiteSpeed Billet Aluminum Wheel w/Poly Tire ☐ Silver Hub **OR** ☐ Black Hub (Not available with 90° Front Angle ❷.).....\$170

XFW30 ☐ 5"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub **OR** ☐ Black Hub (Not available with 90° Front Angle ❷.).....\$215

XFW9 ☐ 6"x1" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or 85° or 90° Front Angle ❷ or Glide™ Single-Sided Suspension Fork.).....NCO

XFW34 ☐ 6"x1" LiteSpeed Billet Aluminum Wheel w/Poly Tire ☐ Silver Hub **OR** ☐ Black Hub (Not available with 85° or 90° Front Angle ❷ or Glide™ Single-Sided Suspension Fork.).....\$190

XFW32 ☐ 6"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Adaptable Frame, Heavy Duty Frame, 85° or 90° Front Angle ❷, or Glide™ Single-Sided Fork.).....\$100

XFW33 ☐ 6"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub **OR** ☐ Black Hub (Not available with 85° or 90° Front Angle ❷ or Glide™ Single-Sided Suspension Fork.).....\$225

XFW16 ☐ 8"x1" Plastic Wheel w/Poly Tire (Not available with Power Assist, Power Add-on, Heavy Duty Frame or Glide™ Single-Sided Available with 70° Front Angle ❷ only.).....(Not available with 12"-14" Seat Depths).....NCO

XFW17 ☐ 8"x1.25" Plastic Wheel w/Pneumatic Tire (Not available with Power Assist or Power Add-on, Heavy Duty Frame or Glide™ Single-Sided Suspension Fork. Available with 70° Front Angle ❷ only.).....(Not available with 12"-14" Seat Depths).....\$100

XGL1 ☐ Out-Front Glide™ Suspension Fork (Available in Black anodize only. Not available with 8" Front Wheels or Heavy Duty Frame).....\$485 E1015

Rider's Weight: ☐ 125 lbs. ☐ 125-180 lbs. ☐ 180-265 lbs.

# TiLite 2GX Order Form

Effective March 5, 2018 (New Items in Red)

FOLDING FRAME

Page 4 of 6

## CUSTOM CONFIGURATION OPTIONS

XGL2 ☐ Out-Front Glide™ Single-Sided Suspension Fork (Available in Black anodize only. Not available with 6" or 8" Front Wheels or Heavy Duty Frame.).....\$625 E1015  
Rider's Weight: ☐ 125 lbs. ☐ 125-180 lbs. ☐ 180-265 lbs.

### STANDARD REAR WHEELS (See Rear Seat Height Ranges Chart on Page 2 for Rear Seat Height limitations.)

XRW14 TiLite SHADOW™.....STD  
☐ 22"(501)  
☐ 24"(540)  
☐ 25"(559) (Requires 15" minimum Seat Depth.)  
☐ 26"(590) (Requires 15" minimum Seat Depth.)

### PERFORMANCE REAR WHEELS (See Rear Seat Height Ranges Chart on Page 2 for Rear Seat Height limitations.)

XRW7 Spinergy SPOX ☐ Black ☐ Red ☐ Yellow ☐ Blue ☐ White .....\$845  
☐ 22"(501) - Std. Hub  
☐ 24"(540) - Std. Hub  
☐ 25"(559) - Std. Hub (Requires 15" minimum Seat Depth.)  
☐ 26"(590) - Sport Hub (Requires 15" minimum Seat Depth. Not available with Quad Release Axles.)  
XRW8 Spinergy LX (Available with Heavy Duty Frame only when 0 or 2 degree camber is selected.) ☐ Black ☐ Red ☐ Yellow ☐ Blue  
☐ White ☐ Go Green ☐ Bright Orange ☐ Flamingo Pink.....\$1,060  
☐ 22"(501) - Std. Hub  
☐ 24"(540) - Std. Hub  
☐ 25"(559) - Std. Hub (Requires 15" minimum Seat Depth.)  
☐ 26"(590) - Sport Hub (Requires 15" minimum Seat Depth. Not available with Quad Release Axles.)  
XRW16 Spinergy Carbon Blade (Not available with Heavy Duty Frame.).....\$1,480  
☐ 24"(540)  
☐ 25"(559) (Requires 15" minimum Seat Depth.)  
XRW17 X-Core Mag (Not available with Heavy Duty Frame, Short Tabs, Quad Release Axles, Spoke Guards or Omit Handrims Option. Requires Aluminum Sliver Anodized Handrims.).....NCO  
☐ 22"(501)  
☐ 24"(540)

### AXLES

XAXL1 ☐ Stainless Quick Release.....STD  
XAXL2 ☐ Titanium Quick Release.....\$255  
XAXL3 ☐ Quad Release (Not available with Mag or 26"(590) Spinergy Spox or Spinergy LX Rear Wheels.).....\$155

### PNEUMATIC REAR TIRES

XRTR1 ☐ Treaded (Not available with 25"(559) Rear Wheels.).....STD  
XRTR2 ☐ Airless Inserts.....\$90 E2213  
XRTR3 ☐ 25"(559) Primo Treaded (Standard with 25"(559) Rear Wheels.).....STD  
XRTR4 ☐ Primo Slick.....NCO  
XRTR30 ☐ 24"(540), 25"(559) and 26"(590) Primo Sentinel .....\$125  
XRTR8 ☐ 24"(540), 25"(559) and 26"(590) Schwalbe® Marathon Plus Evolution.....\$140  
XRTR9 ☐ 24"(540) and 25"(559) Schwalbe® Slick.....\$145  
XRTR24 ☐ 24"(540) Kenda Kobra Knobby.....\$95

### SOLID REAR TIRES

XRTR22 ☐ 22"(501), 24"(540) and 25"(559) Primo Express (Not available with Spinergy Carbon Blade Rear Wheels.) .....\$110  
XRTR23 ☐ 24"(540) and 25"(559) Primo Active .....\$135

### HANDRIMS

#### Handrims by OUT-FRONT®

XRIM1 ☐ Aluminum – Silver Anodized.....STD  
XRIM21 ☐ Aluminum – Black Anodized.....NCO  
XRIM12 ☐ Titanium – Satin Finish (Not available with 20"(451) or 22"(501) Rear Wheels, or 24"(540) Knobby tires.).....\$470  
The Natural-Fit® (Original, Larger Profile)  
XRIM14 ☐ Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....\$395  
XRIM15 ☐ Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....\$430  
XRIM16 ☐ Without Thumb Grip.....\$360  
The Natural-Fit® LT (Reduced Weight, Smaller Profile)  
XRIM17 ☐ Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....\$415  
XRIM18 ☐ Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....\$475  
XRIM19 ☐ Without Thumb Grip.....\$390  
XRIM10 ☐ The Surge (All-In-One Oval with Gription™ Strip, Larger Profile) .....\$400  
XRIM11 ☐ The Surge LT (All-In-One Oval with Gription™ Strip, Smaller Profile).....\$430  
XRIM20 ☐ Q-Grip (Exclusive High-Friction, Anti-Slip Coating) .....\$330

#### Additional Handrim Options

XRIM5 ☐ Plastic Coated – Black.....\$160  
XRIM6 ☐ Projections - 8 Vertical Projections (Not available with 22"(501), 26"(590), Carbon Blade Rear Wheels.).....\$370  
XRIM22 ☐ Omit Handrims.....NCO

#### Tab Length

XRIM7 ☐ Long Tabs.....STD  
XRIM8 ☐ Short Tabs (Not available with Mag Rear Wheels or The Natural-Fit® Handrims with Thumb Grip.).....\$25



# TiLite 2GX Order Form


Effective March 5, 2018 (New Items in Red)

FOLDING FRAME

Page 5 of 6

## CUSTOM CONFIGURATION OPTIONS

### WHEEL LOCKS

Wheel Locks by OUT-FRONT®  (Composite wheel locks are not always compatible with 12"-13" Seat Depth. Uni-locks will be substituted at no charge.)

XWLK7 <input type="checkbox"/> Compact Composite Push to Lock.....	NCO	
XWLK8 <input type="checkbox"/> Compact Composite Pull to Lock .....	NCO	
XWLK10 <input type="checkbox"/> Extension Handles for Compact Composite Push/Pull to Lock.....	\$85	E0961
XWLK9 <input type="checkbox"/> Composite Scissor Lock.....	\$145	

#### Additional Wheel Locks

XWLK3 Uni-Lock Set At: <input type="checkbox"/> Push OR <input type="checkbox"/> Pull.....	\$165	
XWLK4 <input type="checkbox"/> Extension Handles for Uni-Lock.....	\$85	E0961

### BACK UPHOLSTERY (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

XBUP1 <input type="checkbox"/> Padded Nylon (Black Only. Not available with Heavy Duty Frame, Folding Back, Folding Stabilizer Bar, Depth Adjustable Frame, Fold-Down Push Handles or Ultrasuede Backrest Stripe.).....	STD	
Tension Adjustable by Straps (Included with Heavy Duty Frame. Required with Aluminum Fendered Side Guards and Fold-Down Push Handles.)		
XBUP3 <input type="checkbox"/> Padded Nylon.....	\$340	E2611
XBUP5 <input type="checkbox"/> AIR Breathable Upholstery (Not available with Ultrasuede or NiteLite Backrest Stripe.).....	\$340	E2611
XUPC3 Ultrasuede Backrest Stripe <input type="checkbox"/> Black <input type="checkbox"/> Sapphire <input type="checkbox"/> Red Pepper <input type="checkbox"/> Burnt Orange <input type="checkbox"/> Ultraviolet		
<input type="checkbox"/> NiteLite Reflective (Not available with Padded Nylon or AIR Breathable Back Upholstery.).....	\$105	
XBUP0 <input type="checkbox"/> Omit Back Upholstery.....	NCO	

### SEAT UPHOLSTERY

XSUP4 <input type="checkbox"/> Tension Adjustable (Black only).....	STD	
XSUP0 <input type="checkbox"/> Omit Seat Upholstery.....	NCO	



### SEAT CUSHION (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

XCSH1 Foam Cushion <input type="checkbox"/> 2" <input type="checkbox"/> 3" <input type="checkbox"/> 4".....	\$120	
XCSH2 <input type="checkbox"/> Comfort Curve Cushion .....	\$120	E2601



### SIDE GUARDS

XSDG1 <input type="checkbox"/> Cloth (Not available with Ride Designs® Corbac™ Back Support.).....	\$100	
<b>Aluminum</b>		
XSDG2 <input type="checkbox"/> Rigid Removable <input type="checkbox"/> Small OR <input type="checkbox"/> Std.....	\$200	
XSDG3 <input type="checkbox"/> Quad Removable w/Hole <input type="checkbox"/> Small OR <input type="checkbox"/> Std.....	\$215	
XSDG6 <input type="checkbox"/> Fendered Adjustable (Requires Tension Adjustable by Straps Back Upholstery. Not available with Angle-Adjustable Fold-Down Seat Backs.).....	\$335	
<b>Carbon Fiber</b>		
XSDG7 <input type="checkbox"/> Rigid Removable <input type="checkbox"/> Small OR <input type="checkbox"/> Std.....	\$310	
XSDG8 <input type="checkbox"/> Quad Removable w/Hole <input type="checkbox"/> Small OR <input type="checkbox"/> Std.....	\$320	
XSDG11 <input type="checkbox"/> Fendered Adjustable (Requires Tension Adjustable by Straps Back Upholstery. Not available with Angle-Adjustable Fold-Down Seat Backs.).....	\$565	

### ARMRESTS (Desk Arms come with built-in Side Guards.)

XARM1 Swing Away Tubular (Not available with Heavy Duty Frame, Folding Stabilizer Bar with Short or Medium Seat Back Height ranges or Depth Adjustable Back.).....	\$225	E0973
<input type="checkbox"/> 9" - 12" OR <input type="checkbox"/> 12" - 15"		
XARM4 <input type="checkbox"/> Swing Away - Flip Back Tubular (Not available with Heavy Duty Frame, Folding Stabilizer Bar or Depth Adjustable Back.).....	\$315	E0973
XARM2 Removable Height Adjustable Desk Arm with Rigid Side Guard (Requires 14" minimum Seat Depth  .).....	\$335	E0973
<input type="checkbox"/> 9" to 12" OR <input type="checkbox"/> 12" to 15"		
<input type="checkbox"/> Standard Length Pad (10") OR <input type="checkbox"/> Full Length Pad (14")		
XARM3 Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop (Requires 16" minimum Seat Depth  .).....	\$375	E0973
<input type="checkbox"/> 9" to 12" OR <input type="checkbox"/> 12" to 15"		
<input type="checkbox"/> Standard Length Pad (10") OR <input type="checkbox"/> Full Length Pad (14")		
<b>Available Armrest Option</b>		
XARM5 <input type="checkbox"/> Covered in Black Ultrasuede (Not available with Desk Arms).....	\$110	

### SEAT BELTS

XBLT1 <input type="checkbox"/> Velcro Adjustable, Black Nylon Webbing.....	\$55	E0978
XBLT2 <input type="checkbox"/> Bodypoint® 2" Push-Button Buckle (Auto-Style.).....	\$95	E0978
XBLT4 <input type="checkbox"/> Bodypoint® 1" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" - 22" Rear Seat Widths  .).....	\$180	E0978
XBLT5 <input type="checkbox"/> Bodypoint® 1.5" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" - 22" Rear Seat Widths  .).....	\$170	E0978
XBL10 <input type="checkbox"/> Bodypoint® Optional Band Clamp Mounting Hardware.....	\$20	K0108
XBLT15 <input type="checkbox"/> Bodypoint® EvoFlex® Pelvic Stabilizer Swivel Buckl.....	\$185	E0978
XBLT16 <input type="checkbox"/> Bodypoint® EvoFlex® Pelvic Stabilizer Push Button .....	\$185	E0978
XBLT17 <input type="checkbox"/> Bodypoint® EvoFlex® Pelvic Stabilizer Rehab Latch .....	\$185	E0978

# TiLite 2GX Order Form

Effective March 5, 2018 (New Items in Red)

FOLDING FRAME

Page 6 of 6

## CUSTOM CONFIGURATION OPTIONS

### CALF STRAP

XCLF2 Velcro Adjustable.....	\$40	K0038
<input type="checkbox"/> Small (Recommended for 9"- 12" Footrest Widths)		
<input type="checkbox"/> Medium (Recommended for 11"- 15" Footrest Widths)		
<input type="checkbox"/> Large (Recommended for 14"- 17" Footrest Widths)		
<input type="checkbox"/> Xtra Large (Recommended for 16"- 19" Footrest Widths)		
XCLF3 Bodypoint® Padded Velcro Adjustable.....	\$100	K0038
<input type="checkbox"/> Medium (Recommended for 12"- 15" Footrest Widths)		
<input type="checkbox"/> Large (Recommended for 16"- 20" Footrest Widths)		
XCLF1 <input type="checkbox"/> Buckle Adjustable, Black Nylon Webbing.....	\$30	K0038

### BACK PACK

XBCK3 <input type="checkbox"/> Black Permobil Back Pack.....	\$95	
--	------	--

### SEAT POUCH

XPCH3 <input type="checkbox"/> Black Permobil Seat Pouch.....	\$60	
XPCH4 <input type="checkbox"/> Bodypoint® Mobility Bag.....	\$100	

### SPOKE GUARDS

XGRD1 <input type="checkbox"/> Clear with TiMotion Logo (Not available with Mag Rear Wheels.).....	\$135	K0065
--	-------	-------

### ANTI-TIPS

XTIP1 <input type="checkbox"/> Rear Aluminum - Black.....	\$160	E0971
XTIP2 <input type="checkbox"/> User Friendly Flip-up.....	\$335	E0971
XTIP9 <input type="checkbox"/> Tip Assist .....	\$105	E0971

### ACCESSORIES

XUEST1 Upper Extremity Support Tray <input type="checkbox"/> 12" to 14" <b>OR</b> <input type="checkbox"/> 13" to 15" Seat Widths <b>A</b> (Requires Armrests.).....	\$290	E0950
XLUG1 <input type="checkbox"/> Luggage Carriers (Not available with High Mount Swing Away Footrests or 1-Piece Flip-up Footrest.).....	\$145	
XGS1 <input type="checkbox"/> GRIP Solutions Lap Board.....	\$55	
XMPCT2 Neoprene Impact Guards with Ultrasuede Stripe (Not available with High Mount Footrests, 90° Hangers or Elevating Leg Rest.)		
<input type="checkbox"/> Black <input type="checkbox"/> Sapphire <input type="checkbox"/> Red Pepper <input type="checkbox"/> Burnt Orange <input type="checkbox"/> Ultraviolet <input type="checkbox"/> TiLite NiteLite Reflective.....	\$105	
XCH1 <input type="checkbox"/> Crutch Holder.....	\$225	
XAFP1 4" Logo Patch <input type="checkbox"/> Army <input type="checkbox"/> Navy <input type="checkbox"/> Air Force <input type="checkbox"/> Marines .....	\$105	

## GENERAL

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, The Titanium Experts, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademarks of Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. Ride Designs® and Corbac® are registered trademarks of Aspen Seating, LLC. ROHO®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. Supracor®, Stimulite® and Corbee® are registered trademarks of Supracor, Inc. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. VariLite® is a registered trademark of Cascade Designs, Inc.

## QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to Permobil dealers and distributors). A Quote or an Order Confirmation will be sent to you within 24 hours.

## ORDER THIS CHAIR