

# Invacare® Top End® Little Excelerator™ Handcycle



PROVIDER / CUSTOMER INFORMATION:			
Provider Information:		End User Information:	
Date:	Name:	Height:	
IVC Acct #:	Ship to Address:	Weight (lbs.):	
Contact:	Address Line 2:	Leg Length (in.):	
Phone #:	City:	Hip Width (in.):	
Fax/Email:	State:	Zip:	Phone #:
<input type="radio"/> Ship to Provider	<input type="radio"/> Ship to End User		Email:

ORDER REQUEST TYPE:			
<input type="radio"/> Quote	<input type="radio"/> Order	PO #	<input type="text"/>

BASE MODEL:		
<input checked="" type="radio"/> 10008	Top End Little Excelerator Handcycle	\$4,225

Assembly: <i>(see definitions)</i>		
<input checked="" type="radio"/> 10048	Assembled	\$0
Shipping:		
Shipping Fee <i>(Continental US Only - See definitions)</i>		\$150

SEAT WIDTH:		
<input checked="" type="radio"/> 15013	14" Seat Width	\$0

CAMBER ANGLE:		
<input checked="" type="radio"/> 45004	9 Degree Camber	\$0

WHEELBASE: <i>(Front to Rear Axle)</i>		
<input checked="" type="radio"/> 75074	56" Wheelbase	\$0

WHEEL SETS:		
<input checked="" type="radio"/> 50120	20" Front Wheel & 20" Rear Wheels with Quick Release Axles	\$0

TIRES:		
<input checked="" type="radio"/> 50212	20" Treaded Tires	\$0

FOOTREST:		
<input checked="" type="radio"/> 60002	Dual, Height/Angle Adjustable	\$0



DRIVE CONTROL:		
CRANK TYPE:		
<input checked="" type="radio"/> 65074	Crank Arms: 155mm Long by 12" Wide	\$0
HAND PEDAL TYPE: <i>(see definitions)</i>		
<input type="radio"/> 65127	Medium 1-3/4" Medium Ovalized	\$0
<input type="radio"/> 65128	Horizontal Hand Pedals	\$127
<input type="radio"/> 65134	Tri-Pin Hand Pedals	\$361
<input type="radio"/> 65135	Quad Cuff Hand Pedals	\$481
<input type="radio"/> 65133	Quad-Grips Hand Pedals with Wrist Wraps & 3" Crank Width Extensions	\$962
SHIFTING:		
<input checked="" type="radio"/> 65060	Twist Activated Shifter with Internal Gearing	\$0
<input type="radio"/> 65061	Quad Twist Shifting Adaptor <i>(see definitions)</i>	\$181
BRAKING:		
<input checked="" type="radio"/> 65080	Reverse Activated Drum Brake	\$0
<input checked="" type="radio"/> 65090	Locking Parking Brake	\$0
CHAINRING:		
<input checked="" type="radio"/> 65100	Single 38 Tooth Chainring	\$0

# Invacare® Top End® Little Excelerator™ Handcycle



UPHOLSTERY:		
MATERIAL & COLOR:		
<input checked="" type="radio"/>	30001	Nylon, Black Upholstery \$0
SEAT UPHOLSTERY:		
<input checked="" type="radio"/>	30014	Adjustable Tension Seat Upholstery \$0
<input checked="" type="radio"/>	25060	Hook and Loop Positioning Strap \$0
BACK UPHOLSTERY:		
<input checked="" type="radio"/>	30029	Fold - Over Back Upholstery \$0
SEAT CUSHION:		
<input checked="" type="radio"/>	25063	Seat Cushion \$0

FRAME COLOR:			
<input type="radio"/>	55018	Black	Wet Black \$0
<input type="radio"/>	55019	Black	Black Prism \$0
<input type="radio"/>	55017	Black	Stealth Black \$0
<input type="radio"/>	55010	Blue	Midnight Blue \$0
<input type="radio"/>	55011	Blue	Royal Blue \$0
<input type="radio"/>	55012	Blue	Electric Blueberry \$0
<input type="radio"/>	55023	Grey	Stealth Grey \$0
<input type="radio"/>	55053	Green	Agriculture Green \$0
<input type="radio"/>	55008	Green	Grasshopper Green \$0
<input type="radio"/>	55054	Green	Neon Green \$0
<input type="radio"/>	55006	Orange	Sunkissed Orange \$0
<input type="radio"/>	55005	Orange	Simple Orange \$0
<input type="radio"/>	55051	Orange	Neon Orange \$0
<input type="radio"/>	55015	Pink	Bubblegum Pink \$0
<input type="radio"/>	55057	Pink	Neon Pink \$0
<input type="radio"/>	55055	Purple	Lilac Purple \$0
<input type="radio"/>	55001	Red	Simple Red \$0
<input type="radio"/>	55002	Red	Tangerine Red \$0
<input type="radio"/>	55003	Red	Raspberry Red \$0
<input type="radio"/>	55050	Red	Neon Red \$0
<input type="radio"/>	55020	Silver	Silver Metallic \$0
<input type="radio"/>	55016	White	White Snow \$0
<input type="radio"/>	55007	Yellow	Safety Yellow \$0
<input type="radio"/>	55052	Yellow	Neon Yellow \$0

ACCESSORIES:		
STANDARD:		
<input checked="" type="radio"/>	25061	Safety Flag \$0
<input checked="" type="radio"/>	25062	Chain Guard \$0
UPGRADABLE:		
SAFETY:		
<input type="radio"/>	25015	Light, Front (White) \$18
<input type="radio"/>	25016	Light, Rear (Red) \$18
<input type="radio"/>	25097	Mirror, Rear View \$37
HYDRATION:		
<input type="radio"/>	25021	Water Bottle and Cage \$60
POSITIONING:		
<input type="radio"/>	25001	Click Strap, 8" Small \$90
<input type="radio"/>	25002	Click Strap, 10" Medium \$97
<input type="radio"/>	25003	Click Strap, 12" Large \$109
PERFORMANCE ENHANCEMENTS:		
<input type="radio"/>	25023	Crank Width Extensions (+1.5" of crank width per set) \$24
<input type="radio"/>	25022	Wheel Alignment Gauge \$259
TRANSPORTATION:		
<input type="radio"/>	25006	Crutch Holder \$181
APPAREL:		
<input type="radio"/>	25084	Lemon Yellow Top End T-Shirt \$24
<input type="radio"/>	25085	Grey Top End T-Shirt \$24
<input type="radio"/>	25086	Grey Top End Racer Back Tank Top \$24
<i>Please select a size for your shirt selection:</i>		
<input type="radio"/>	25087	Extra Small
<input type="radio"/>	25088	Small
<input type="radio"/>	25089	Medium
<input type="radio"/>	25090	Large
<input type="radio"/>	25091	Extra Large
<input type="radio"/>	25092	XX Large
<input type="radio"/>	25093	XXX Large

E-mail to Top End Order Box

**Invacare Corporation Top End Division**

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.

©2015 Invacare Corporation

1181421 REV F KPIPARO 02-2022

[Reset Form](#) [Clear All](#)

# Invacare® Top End® Little Excelerator™ Handcycle



HANDCYCLE SPECIFICATIONS:	
<b>Weight:</b>	45 pounds
<b>Turning Radius:</b>	10 feet for 360 degree circle, 7 feet on a 3-point turn
<b>Transfer Height:</b>	17.5" from ground with 2" cushion
<b>Frame Style:</b>	Single main tube, upright position, 1-3/4" X .065" ovalized carbon steel tubing
<b>Seat:</b>	Seat adjust forward and rearward with sliding mechanism
<b>Fork:</b>	Carbon Steel
<b>Propulsion / Shifting:</b>	Internal front hub, twist activated shifter mounted on crank support
<b>Cranks:</b>	155mm adjustable height
<b>Back:</b>	Adjustable height
<b>Footrest:</b>	Dual adjustable height and angle, includes straps
<b>Hand Pedals:</b>	Ovalized ergonomic with frictionless bearings
<b>Brake:</b>	Reverse activated drum and Locking Parking Brake mounted on fork
<b>Wheelbase:</b>	56" axle to axle
<b>Camber:</b>	9 degree
<b>Wheels:</b>	20" metal spoked wheels
<b>Axles:</b>	Quick Release
<b>Tires / Tubes:</b>	20" X 1.75" tires (40 PSI)
<b>Positioning Straps:</b>	Hook and loop strap included
<b>Cushions:</b>	Seat cushion included
<b>Safety / Visibility:</b>	Flag and wheel reflectors included
<b>Chain Guard:</b>	Full length, shock resistant composite material

DEFINITIONS:	
<b>Assembly and Shipping:</b>	Arrives factory assembled, rear wheels and some accessories will need to be installed. Shipping rate is \$150 for all shipments in the continental US. Freight quotes for shipments outside the continental US are available upon request.
<b>Quad Twist Shifter Adaptor:</b>	Three pronged adaptor which allows user to tap for shifting.
<b>Tri-pin Hand Pedals:</b>	3-Point contact hand pedal similar to hand controls used for vehicles. Allows for the hands to be kept in place for people with limited grip strength.
<b>Quad Grips Hand Pedals:</b>	Combination of Tri-Pin hand pedal and Quad Cuff hand pedal. Allows for the hands to be kept in place for people with limited grip strength.
<b>Quad Cuff Hand Pedals:</b>	Includes Power Plates on the top and the bottom of the pedal to assist through both up and down stroke. Includes an adjustable metal cuff with a hook and look strap to secure hands.

### Warranty Info:

Invacare warrants the frames when purchased new and unused to be free from defects in materials and workmanship for a period of three (3) years from the date of purchase from Invacare or a dealer, with a copy of the seller's invoice required for coverage under this warranty. Invacare warrants the upholstered materials (seat and back) and remaining components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of thirteen (13) months from date of purchase from Invacare or a dealer, with a copy of the seller's invoice required for coverage under this warranty.