

Invacare® Top End® XLT™ Handcycle



PROVIDER / CUSTOMER INFORMATION:			
Provider Information:		End User Information:	
Date:	Name:	Height:	
IVC Acct #:	Ship to Address:	Weight (lbs.):	
Contact:	Address Line 2:	Leg Length (in.):	
Phone #:	City:	Hip Width (in.):	
Fax/Email:	State:	Zip:	Phone #:
<input type="radio"/> Ship to Provider		<input type="radio"/> Ship to End User	
		Email:	

ORDER REQUEST TYPE:		
<input type="radio"/> Quote	<input type="radio"/> Order	PO # <input type="text"/>

BASE MODEL:		
<input checked="" type="radio"/> 10004	Top End XLT Handcycle	\$4,490

Assembly: <i>(see definitions)</i>		
<input checked="" type="radio"/> 10048	Assembled	\$0
Shipping:		
Shipping Fee <i>(Continental US Only - See definitions)</i>		\$150

SEAT WIDTH:		
<input type="radio"/> 15015	15" Seat Width	\$0
<input type="radio"/> 15019	17" Seat Width	\$0

WHEELBASE: <i>(Front to Rear Axle)</i>		
<input checked="" type="radio"/> 75024	53" Wheelbase	\$0

CAMBER ANGLE:		
<input checked="" type="radio"/> 45008	15 Degree Camber	\$0

WHEEL SETS:		
<input checked="" type="radio"/> 50122	26" (571) High Performance Front and Rear Wheels with Quick Release Axles	\$0

TIRES:		
<input checked="" type="radio"/> 50214	26" X 1" High Pressure Clincher Tires	\$0

FOOTREST:		
<input checked="" type="radio"/> 60004	Dual, 4.5" Wide, Adjustable Length	\$0

HAND PEDALS:		
<input type="radio"/> 65127	Medium 1-3/4" Ovalized	\$0
<input type="radio"/> 65128	Horizontal Hand Pedals	\$127
<input type="radio"/> 65134	Tri-Pin Hand Pedals	\$361
<input type="radio"/> 65135	Quad Cuff Hand Pedals	\$481
<input type="radio"/> 65133	Quad Grips Hand Pedals with Wrist Wraps & 3" Crank Width Extensions	\$962



DRIVE CONTROL:		
SHIFTING:		
<input checked="" type="radio"/> 65060	Twist Activated Shifter with Internal Gearing	\$0
<input type="radio"/> 65061	Quad Twist Shifting Adaptor <i>(see definitions)</i>	\$181
BRAKING:		
<input checked="" type="radio"/> 65080	Reverse Activated Drum Brake	\$0
<input checked="" type="radio"/> 65090	Locking Parking Brake	\$0
CHAINRING:		
<input checked="" type="radio"/> 65100	Single 38 tooth Chainring	\$0
CRANK TYPE: <i>(see definitions)</i>		
<input type="radio"/> 65054	Standard: 175mm long by 12" wide	\$0
<input type="radio"/> 65141	V - Crankset	\$360
<input type="radio"/> 65140	Mountain Drive	\$1,204
<input type="radio"/> 65142	Mountain Drive with V Cranks	\$1,565

UPHOLSTERY:		
MATERIAL & COLOR:		
<input checked="" type="radio"/> 30001	Nylon, Black Upholstery	\$0
SEAT UPHOLSTERY:		
<input checked="" type="radio"/> 30014	Adjustable Tension Seat Upholstery	\$0
<input checked="" type="radio"/> 25060	Hook and Loop Positioning Strap	\$0
BACK UPHOLSTERY:		
<input checked="" type="radio"/> 25064	Back Cushion	\$0
SEAT CUSHION:		
<input checked="" type="radio"/> 25063	Seat Cushion	\$0

Invacare Corporation Top End Division

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.

©2015 Invacare Corporation

1181429 REV E KPIPARO 02-2022

Invacare® Top End® XLT™ Handcycle



ACCESSORIES:		
STANDARD:		
<input checked="" type="radio"/>	25061 Safety Flag	\$0
<input checked="" type="radio"/>	25062 Chain Guard	\$0
<input checked="" type="radio"/>	25099 Mirror, Footrest Mounted	\$0
UPGRADABLE:		
SAFETY:		
<input type="radio"/>	25015 Light, Front <i>(White)</i>	\$18
<input type="radio"/>	25016 Light, Rear <i>(Red)</i>	\$18
<input type="radio"/>	25017 Leg Guard	\$60
HYDRATION:		
<input type="radio"/>	25021 Water Bottle and Cage	\$60
POSITIONING:		
<input type="radio"/>	25001 Click Strap, 8" Small	\$90
<input type="radio"/>	25002 Click Strap, 10" Medium	\$97
<input type="radio"/>	25003 Click Strap, 12" Large	\$109
PERFORMANCE ENHANCEMENTS:		
<input type="radio"/>	25023 Crank Width Extensions <i>(+1.5" of crank width per set)</i>	\$24
<input type="radio"/>	25024 Cordless Computer <i>(speed, distance, etc.)</i>	\$120
<input type="radio"/>	25022 Wheel Alignment Gauge	\$259
<input type="radio"/>	25036 Indoor Stationary Training Roller	\$301
TRANSPORTATION:		
<input type="radio"/>	25006 Crutch Holder	\$181
<input type="radio"/>	25025 Tow Bar for Wheelchair	\$271
APPAREL:		
<input type="radio"/>	25084 Lemon Yellow Top End T-Shirt	\$24
<input type="radio"/>	25085 Grey Top End T-Shirt	\$24
<input type="radio"/>	25086 Grey Top End Racer Back Tank Top	\$24
<i>Please select a size for your shirt selection:</i>		
<input type="radio"/>	25087 Extra Small	
<input type="radio"/>	25088 Small	
<input type="radio"/>	25089 Medium	
<input type="radio"/>	25090 Large	
<input type="radio"/>	25091 Extra Large	
<input type="radio"/>	25092 XX Large	
<input type="radio"/>	25093 XXX Large	

FRAME COLOR:				
<input type="radio"/>	55018	Black	Wet Black	\$0
<input type="radio"/>	55019	Black	Black Prism	\$0
<input type="radio"/>	55017	Black	Stealth Black	\$0
<input type="radio"/>	55010	Blue	Midnight Blue	\$0
<input type="radio"/>	55011	Blue	Royal Blue	\$0
<input type="radio"/>	55012	Blue	Electric Blueberry	\$0
<input type="radio"/>	55023	Grey	Stealth Grey	\$0
<input type="radio"/>	55053	Green	Agriculture Green	\$0
<input type="radio"/>	55008	Green	Grasshopper Green	\$0
<input type="radio"/>	55054	Green	Neon Green	\$0
<input type="radio"/>	55006	Orange	Sunkissed Orange	\$0
<input type="radio"/>	55005	Orange	Simple Orange	\$0
<input type="radio"/>	55051	Orange	Neon Orange	\$0
<input type="radio"/>	55015	Pink	Bubblegum Pink	\$0
<input type="radio"/>	55057	Pink	Neon Pink	\$0
<input type="radio"/>	55055	Purple	Lilac Purple	\$0
<input type="radio"/>	55001	Red	Simple Red	\$0
<input type="radio"/>	55002	Red	Tangerine Red	\$0
<input type="radio"/>	55003	Red	Raspberry Red	\$0
<input type="radio"/>	55050	Red	Neon Red	\$0
<input type="radio"/>	55020	Silver	Silver Metallic	\$0
<input type="radio"/>	55016	White	White Snow	\$0
<input type="radio"/>	55007	Yellow	Safety Yellow	\$0
<input type="radio"/>	55052	Yellow	Neon Yellow	\$0

Note: All Special requests are subject to approval, additional charges, and additional lead time.

Special Request 1)	
Special Request 2)	
Special Request 3)	
Special Request 4)	
Special Request 5)	

Invacare Corporation Top End Division

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.

©2015 Invacare Corporation

1181429 REV E KPIPARO 02-2022

Invacare® Top End® XLT™ Handcycle



HANDCYCLE SPECIFICATIONS:	
Weight:	Approximately 35 pounds
Turning Radius:	18 feet for 360 degree circle, 12 feet on a 3-point turn
Transfer Height:	Front 13", Rear 9" from ground with 2" cushion
Frame Style:	Single main tube, recumbent position, 2.5" X .065" ovalized 6061 aluminum aircraft tubing
Fork:	4130 chrome-moly, raked angle
Propulsion / Shifting:	Internal front hub, twist activated shifter mounted on left side of seat frame
Cranks:	175mm X 12" adjustable height
Back:	Fixed height at 20"
Footrest:	Dual adjustable height and angle, includes straps
Hand Pedals:	Ovalized ergonomic with frictionless bearings
Brake:	Reverse Activated Drum and Locking Parking Brake mounted on right side of seat frame
Overall length:	72" overall, tire to tire
Camber:	15 degree
Wheels:	26" Metal spoked wheels
Axles:	Quick release
Tires / Tubes:	26" X 1" high pressure clincher tires (100 PSI)
Positioning	Hook and loop seat restraint included
Cushions:	Seat and back included
Safety / Visibility:	Flag and wheel reflectors included
Chain Guard:	Full length, shock resistant composite material except when Mountain Drive is chosen. With Mountain Drive selection a half moon shaped black composite chainguard is used.

DEFINITIONS:	
Assembly and Shipping:	Arrives factory assembled, rear wheels and some accessories will need to be installed. Shipping rate is \$150 for all shipments in the continental US. Freight quotes for shipments outside the continental US are available upon request.
Quad Twist Shifter Adaptor:	Three pronged adaptor which allows user to tap for shifting
Tri-pin Hand Pedals:	3-Point contact hand pedal similar to hand controls used for vehicles. Allows for the hands to be kept in place for people with limited grip strength.
Quad Grips Hand Pedals:	Combination of Tri-Pin hand pedal and Quad Cuff hand pedal. Allows for the hands to be kept in place for people with limited grip strength.
Quad Cuff Hand Pedals:	Includes Power Plates on the top and the bottom of the pedal to assist through both up and down stroke. Includes an adjustable metal cuff with a hook and loop strap to secure hands.
V - Crankset:	Option for people who need more clearance due to long legs or large midsections. Cranks are 180mm long and 20" wide.
Mountain Drive:	Doubles the gears adding lower gears for help climbing hills. Works similar to a transmission.
Mountain Drive with V-Cranks:	Combines the V Crank upgrade and the Mountain Drive upgrade so the handcycle is equipped with both.

Warranty Info:

Invacare warrants the frames when purchased new and unused to be free from defects in materials and workmanship for a period of three (3) years from the date of purchase from Invacare or a dealer, with a copy of the seller's invoice required for coverage under this warranty. Invacare warrants the upholstered materials (seat and back) and remaining components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of thirteen (13) months from date of purchase from Invacare or a dealer, with a copy of the seller's invoice required for coverage under this warranty