

# ODYSSEY LX



Available  
**1 Pack**

## Odyssey LX 4 Wheel Full Sized Scooter

S45200LX-20

Red, 20" Seat Width, 1/cs

S45200LX-22

Red, 22" Seat Width, 1/cs

### FEATURES:

- Combines attractive styling with rugged construction, featuring a heavy duty drive train.
- The Captains seat has height adjustment in addition to the auto style depth adjustment, allows for the perfect seat positioning of every operator.
- The reclining contour seat and adjustable headrest adds the final touch of comfort to every ride.
- Flip back deluxe arms are height, width and angle adjustable so you can achieve that perfect position.
- Full front suspension insures a comfortable ride.
- Standard anti-tippers are added safety feature.
- Top speed of 6 mph and a cruising range of 20 miles more than covers those extended trips between battery charges
- Heavy duty front bumpers protect the headlight and front shroud.
- Large carry basket is also included for all those personal items
- Front and rear lights with turn signals are additional safety feature
- Flat-free black tires are non-marking and maintenance-free. They provide a secure and comfortable ride
- Easy to adjust infinite tiller range will help find the ideal position for all operators.
- 400 lb. Weight Capacity.

## SPECIFICATIONS

ITEM #	S45200LX-20, S45200LX-22
HCPCS CODE GROUP	K0899
DIMENSION (L X W X H)	55.1" X 26.4" X 39.8"
GROSS WEIGHT (WITH/WITHOUT BATTERY)	210 LBS./162 LBS.
BATTERY WEIGHT	48 LBS.
SEAT WEIGHT	26 LBS.
SEAT DIMENSIONS (D X W)	20" X 20" /20" X 22"
TIRE SIZE (FRONT)	13" SOLID
TIRE SIZE (REAR)	13" SOLID
BATTERY (INCLUDED)	112V, 36AH X2
CHARGER	4 AMP OFF BOARD
MAX. SPEED	6 MPH
CRUISING RANGE	20~25 MILES
CLIMBING ANGLE	8°
TURNING RADIUS	70"
GROUND CLEARANCE	3.5"
CONTROLLER DYNAMIC	75 AMP MULTI MODE
WEIGHT CAPACITY	400 LBS.
WARRANTY	LIFETIME ON FRAME, ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS, SIX MONTHS ON BATTERY

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

powered mobility

To order call toll free: 877.224.0946

drive™