



Company	<input type="text"/>	Name	<input type="text"/>
Address 1	<input type="text"/>	Phone	<input type="text"/>
Address 2	<input type="text"/>	Email *	<input type="text"/>
City	<input type="text"/>	PO Number	<input type="text"/>
State	<input type="text"/>	Reference	<input type="text"/>
Postal Code	<input type="text"/>	Notes	<input type="text"/>

ShowerBuddy SB1

The ShowerBuddy (SB1) transfer system will transform your existing curbed shower into an accessible, barrier-free shower in minutes. This chair allows the user to enjoy showering again without the high cost, delay, and stress of a full bathroom remodel. The SB1 alleviates heavy lifting and multiple transfers for toileting and showering. This chair is separated into three sections: rolling base (shower chair), bridge, and shower base. The rolling base can double as a bedside commode. It also rolls around the house like a wheelchair, allowing a user to easily transfer and roll into the bathroom area and over the toilet. When toileting is complete, the rolling base and shower base are connected using the bridge, allowing the chair to slide over the curb and into the shower stall. Once inside the shower stall, the bridge disconnects, and the shower base can rotate 360 degrees to allow the user to face the shower controls or any other direction. SB1 features a unique secondary safety mechanism that prevents the chair from sliding off the track and causing injury. The SB1 chair comes equipped with a standard front-opening cushion, footplate, flip-back, lockable and removable armrests, 5" lockable casters, a removable commode bucket with lid, chest belt, lap belt and height adjustability. ShowerBuddy (SB1) is a caregiver-assisted product; it is not intended for independent use. Note: SB1 only fits over a limited amount of toilets.

The SB1 weighs 54.8 lbs with a weight capacity of 350 lbs.

1. Primary	Description	MSRP	Quantity
SB1	ShowerBuddy SB1	\$3,395.00	1

2. Accessories		Description	MSRP	Quantity
SBC-5	Custom Padded Back Cushion - SB1, SB2, SB2T, SB3T & SB6		\$685.00	
B2	Lateral Kit (Right)		\$460.00	
B10	Large Headrest Kit (non tilting chairs)		\$395.00	
B8	Leg Adductor Kit		\$383.00	
B6	Padded Back Kit (SB1, SB2, SB2T,SB3T)		\$441.00	
B4	BodyPoint Chest Belt Kit (SB1, SB2, SB2T, SB3T)		\$242.50	
B5	BodyPoint Lap Belt Kit (SB1, SB2, SB2T, SB3T)		\$242.50	
SBLT-B	SolutionBased Belt for Laterals		\$68.00	
SBACD100	ShowerBuddy Channel Armrests (Each) - SB1, SB2, SB2T, SB3T		\$165.00	
TBR	ShowerBuddy Tensioned Backrest Adjustable - SB1, SB2, SB2T, SB3T		\$85.00	
B1	Lateral Kit (Left)		\$460.00	
BP-LH-S	ShowerBuddy BodyPoint Leg Harness - Small		\$110.00	
BP-LH-M	ShowerBuddy BodyPoint Leg Harness - Medium		\$110.00	
B9	Neckrest Assembly		\$230.00	

3. Other		Description	MSRP	Quantity
LCB	Custom Size Bridge (Labor)		\$300.00	
SBR	ShowerBuddy Base Risers 4 inches (set of 4) - SB1		\$100.00	
SBR200	ShowerBuddy Base Risers 8 inches (set of 4) - SB1		\$100.00	
SPH	ShowerBuddy Straight Push Handle		\$96.00	
SBR150	ShowerBuddy Base Risers 6 inches (set of 4) - SB1		\$100.00	
49	ShowerBuddy Suction Cup (Each) - SB1, SB2, SB2T		\$30.00	
SB1E	ShowerBuddy Extension Tubes for SB1		\$180.00	

4. Pediatrics		Description	MSRP	Quantity
RW-M-B		Blue Rifton Wave (Medium) with Chest Strap with Laterals Head Blocks and Calf Rest	\$1,210.00	
P1		ShowerBuddy P1 Petite Converter Kit	\$1,195.00	
P2		ShowerBuddy P2 Petite Converter Kit	\$1,295.00	
RW-S-B		Blue Rifton Wave (Small) with Chest Strap with Laterals Head Blocks and Calf Rest	\$1,165.00	
Rifton Hardware Kit		Rifton Hardware Kit	\$575.00	
B-16		Lateral Kit for Pediatric Converter Kit (P1 & P2)	\$1,200.00	
RW-L-B		Blue Rifton Wave (Large) with Chest Strap with Laterals Head Blocks and Calf Rest	\$1,340.00	
SBC-3		Custom Padded Back Cushion - P1 & P2	\$685.00	

5. Seating options		Description	MSRP	Quantity
B15		22" Wide Bariatric Kit (SB1, SB2, SB2T, SB3T)	\$1,305.00	
B7		Roho Cushion Kit	\$900.00	
SCCF-2		ShowerBuddy Seat Cushion Closed Front (SB1-SB3T)	\$195.00	
SBC-2		ShowerBuddy Custom Visco Seat Cushion	\$745.00	