

SB2 Order Form

For assistance call (877) 769-2833



Company	<input type="text"/>	Name	<input type="text"/>
Address 1	<input type="text"/>	Phone	<input type="text"/>
Address 2	<input type="text"/>	Email *	<input type="text"/>
City	<input type="text"/>	PO Number	<input type="text"/>
State	<input type="text"/>	Reference	<input type="text"/>
Postal Code	<input type="text"/>	Notes	<input type="text"/>

ShowerBuddy SB1

The ShowerBuddy (SB1) transfer system will transform your existing curbed shower into an accessible, barrier-free shower in minutes. This chair allows the user to enjoy showering again without the high cost, delay, and stress of a full bathroom remodel. The SB1 alleviates heavy lifting and multiple transfers for toileting and showering. This chair is separated into three sections: rolling base (shower chair), bridge, and shower base. The rolling base can double as a bedside commode. It also rolls around the house like a wheelchair, allowing a user to easily transfer and roll into the bathroom area and over the toilet. When toileting is complete, the rolling base and shower base are connected using the bridge, allowing the chair to slide over the curb and into the shower stall. Once inside the shower stall, the bridge disconnects, and the shower base can rotate 360 degrees to allow the user to face the shower controls or any other direction. SB1 features a unique secondary safety mechanism that prevents the chair from sliding off the track and causing injury. The SB1 chair comes equipped with a standard front-opening cushion, footplate, flip-back, lockable and removable armrests, 5" lockable casters, a removable commode bucket with lid, chest belt, lap belt and height adjustability. ShowerBuddy (SB1) is a caregiver-assisted product; it is not intended for independent use. Note: SB1 only fits over a limited amount of toilets.

The SB1 weighs 54.8 lbs with a weight capacity of 350 lbs.

1. Primary	Description	MSRP	Quantity
SB1	ShowerBuddy SB1	\$3,395.00	1

2. Accessories	Description	MSRP	Quantity
BP-B	ShowerBuddy Bodypoint Belt	\$181.50	
SB-WST-R	SB Large Headrest Socket-Mounting Tube	\$80.00	
TBR	ShowerBuddy Tensioned Backrest Adjustable - SB1, SB2, SB2T, SB3T	\$85.00	
34-22	ShowerBuddy 22" Arm Rest Right (Each) - SB1, SB2, SB2T & SB3t	\$90.00	
FFP-18	Flat Foot Plate 18" - SB1, SB2, SB2T & SB3T	\$60.00	
33-22	ShowerBuddy 22" Arm Rest Left (Each) - SB1, SB2, SB2T, SB3T	\$90.00	
PBC	ShowerBuddy Solid Padded Backrest	\$195.00	
BP-B-PTC32	ShowerBuddy BP Lap Belt Tube Clamp 32mm (each)	\$16.00	
SB-HL-100-US	SB Large Headrest, Ultra Soft, Adjustable	\$165.00	
SB-WST-TA	ShowerBuddy Large Headrest D-shaped Tube - SB2T, SB3T	\$15.00	
SBLT-L	ShowerBuddy Lateral (LEFT) - SB1, SB2, SB2T, SB3T	\$395.00	
SBLT-R	ShowerBuddy Lateral (Right) - SB1, SB2, SB2T & SB3T	\$395.00	
70A	ShowerBuddy Neckrest - SB2T, SB3T	\$95.00	
SBLAL	ShowerBuddy Leg Adductor (Each) - SB1, SB2, SB2T, SB3T	\$318.00	
BP-B-QR	ShowerBuddy BP Belt Quick Release Knobs (Pair)	\$29.00	
BP-B-PTC25	ShowerBuddy BP Chest Belt Tube Clamp 25mm (each)	\$16.00	
SBACD100	ShowerBuddy Channel Armrests (Each) - SB1, SB2, SB2T, SB3T	\$125.00	
SB-SS-100	SB Lateral/Side Support Pad (each)	\$65.00	

3. Other	Description	MSRP	Quantity
LCB	Custom Size Bridge (Labor)	\$300.00	
SBR200	ShowerBuddy Base Riser 8 inches (set of 4) - SB1	\$100.00	
SPH	ShowerBuddy Straight Push Handle	\$96.00	
SBR	ShowerBuddy Base Risers 4 inches (set of 4) - SB1	\$100.00	
SB1E	ShowerBuddy Extension Tubes for SB1	\$180.00	
SBR150	ShowerBuddy Base Risers 6 inches (set of 4) - SB1	\$100.00	
SBRD	ShowerBuddy Redirectional / Stabilizer	\$395.00	
49	ShowerBuddy Suction Cup - SB1, SB2T	\$30.00	

4. Pediatrics		Description	MSRP	Quantity
P1		ShowerBuddy P1 Petite Converter	\$1,195.00	
P2		ShowerBuddy P2 Petite Converter Kit	\$1,295.00	
PSA		ShowerBuddy Pediatric Seat Adaptor	\$995.00	

5. Seating options		Description	MSRP	Quantity
SCSF-22		ShowerBuddy 22" Soft Cushion Seat Foam - SB1, SB2, SB2T & SB3T	\$90.00	
SCO		ShowerBuddy Seat Cushion Overlay (SB1-SB3T)	\$155.00	
SCP1-500-166		ShowerBuddy Pediatric Seat Cushion 6" Commode Opening (SCP1)	\$195.00	
SCP2-500-167		ShowerBuddy Pediatric Seat Cushion 6.6" Commode Opening (SCP2)	\$195.00	
SCSF-05P		ShowerBuddy Soft Foam Cushion w/ Petite Opening (SB1, SB2, SB2T & SB3T)	\$165.00	
SCSF		ShowerBuddy Soft Cushion (SB3t, SB2t, SB2, SB1)	\$140.00	
30		ShowerBuddy Cushion - SB1, SB2, SB2T, SB3T	\$170.00	
SCCF-2		ShowerBuddy Seat Cushion Closed Front (SB1-SB3T)	\$195.00	
PM-1		ShowerBuddy Roho Shower Commode Cushion	\$675.00	