

SB2 Order Form

For assistance call (877) 769-2833



Company	<input type="text"/>	Name	<input type="text"/>
Address 1	<input type="text"/>	Phone	<input type="text"/>
Address 2	<input type="text"/>	Email *	<input type="text"/>
City	<input type="text"/>	PO Number	<input type="text"/>
State	<input type="text"/>	Reference	<input type="text"/>
Postal Code	<input type="text"/>	Notes	<input type="text"/>

ShowerBuddy SB2 TubBuddy

The TubBuddy (SB2) transfer system will transform your existing bathtub or curbed shower into an accessible, barrier-free shower in minutes. It allows the user to enjoy showering again without the high cost, delay, and stress of a full bathroom remodel. The SB2 alleviates heavy lifting and multiple transfers. This chair is separated into three sections: rolling base (shower chair), bridge, and tub base. The rolling base can double as a bedside commode. It also rolls around the house like a wheelchair, allowing users to transfer comfortably and roll over most toilets to be used as an over-toilet commode. When toileting is complete, the rolling base and tub base are connected using the bridge, allowing the unit to slide over the tub lip or curbed shower. Once over the tub/shower, the rolling base and bridge can be disconnected and removed so that the shower curtain can either be closed or allow the caregiver space to move closer to the patient. A custom bridge may be required if the toilet is too close or too far from the tub. The SB2 features a secondary safety mechanism that prevents the chair from rolling off the tracks and causing injury. The SB2 chair comes equipped with a standard front-opening cushion, footrests, flip-back, lockable and removable armrests, 5" lockable casters, a removable commode bucket with lid, chest belt, lap belt, and height adjustability. The TubBuddy (SB2) is a caregiver-assisted product; it is not intended for independent use.

The SB2 weighs 59.5 lbs with a weight capacity of 350 lbs.

1. Primary	Description	MSRP	Quantity
SB2	ShowerBuddy SB2 TubBuddy	\$3,695.00	1

2. Accessories	Description	MSRP	Quantity
BP-B	ShowerBuddy Bodypoint Belt	\$181.50	
SB-WST-R	SB Large Headrest Socket-Mounting Tube	\$80.00	
TBR	ShowerBuddy Tensioned Backrest Adjustable - SB1, SB2, SB2T, SB3T	\$85.00	
34-22	ShowerBuddy 22" Arm Rest Right (Each) - SB1, SB2, SB2T & SB3t	\$90.00	
FFP-18	Flat Foot Plate 18" - SB1, SB2, SB2T & SB3T	\$60.00	
33-22	ShowerBuddy 22" Arm Rest Left (Each) - SB1, SB2, SB2T, SB3T	\$90.00	
PBC	ShowerBuddy Solid Padded Backrest	\$195.00	
BP-B-PTC32	ShowerBuddy BP Lap Belt Tube Clamp 32mm (each)	\$16.00	
SB-HL-100-US	SB Large Headrest, Ultra Soft, Adjustable	\$165.00	
SB-WST-TA	ShowerBuddy Large Headrest D-shaped Tube - SB2T, SB3T	\$15.00	
SBLT-L	ShowerBuddy Lateral (LEFT) - SB1, SB2, SB2T, SB3T	\$395.00	
SBLT-R	ShowerBuddy Lateral (Right) - SB1, SB2, SB2T & SB3T	\$395.00	
70A	ShowerBuddy Neckrest - SB2T, SB3T	\$95.00	
SBLAL	ShowerBuddy Leg Adductor (Each) - SB1, SB2, SB2T, SB3T	\$318.00	
BP-B-QR	ShowerBuddy BP Belt Quick Release Knobs (Pair)	\$29.00	
TBE24	ShowerBuddy Tub Base Extension Rail Pair - SB2, SB2T	\$220.00	
BP-B-PTC25	ShowerBuddy BP Chest Belt Tube Clamp 25mm (each)	\$16.00	
SBACD100	ShowerBuddy Channel Armrests (Each) - SB1, SB2, SB2T, SB3T	\$125.00	
SB-SS-100	SB Lateral/Side Support Pad (each)	\$65.00	

3. Other	Description	MSRP	Quantity
LCB	Custom Size Bridge (Labor)	\$300.00	
SPH	ShowerBuddy Straight Push Handle	\$96.00	
RBE24	ShowerBuddy Rolling Base Extension Rail Pair	\$220.00	
SBRD	ShowerBuddy Redirectional / Stabilizer	\$395.00	
49	ShowerBuddy Suction Cup - SB1, SB2T	\$30.00	

4. Pediatrics	Description	MSRP	Quantity
P1	ShowerBuddy P1 Petite Converter	\$1,195.00	
P2	ShowerBuddy P2 Petite Converter Kit	\$1,295.00	
PSA	ShowerBuddy Pediatric Seat Adaptor	\$995.00	

5. Seating options	Description	MSRP	Quantity
SCSF-22	ShowerBuddy 22" Soft Cushion Seat Foam - SB1, SB2, SB2T & SB3T	\$90.00	
SCO	ShowerBuddy Seat Cushion Overlay (SB1-SB3T)	\$155.00	
SCP1-500-166	ShowerBuddy Pediatric Seat Cushion 6" Commode Opening (SCP1)	\$195.00	
SCP2-500-167	ShowerBuddy Pediatric Seat Cushion 6.6" Commode Opening (SCP2)	\$195.00	
SCSF-05P	ShowerBuddy Soft Foam Cushion w/ Petite Opening (SB1, SB2, SB2T & SB3T)	\$165.00	
SCSF	ShowerBuddy Soft Cushion (SB3t, SB2t, SB2, SB1)	\$140.00	
30	ShowerBuddy Cushion - SB1, SB2, SB2T, SB3T	\$170.00	
SCCF-2	ShowerBuddy Seat Cushion Closed Front (SB1-SB3T)	\$195.00	
PM-1	ShowerBuddy Roho Shower Commode Cushion	\$675.00	