

ShowerBuddy SB2 TubBuddy Order Form

For assistance call (877) 769-2833



Company	<input type="text"/>	Name	<input type="text"/>
Address 1	<input type="text"/>	Phone	<input type="text"/>
Address 2	<input type="text"/>	Email *	<input type="text"/>
City	<input type="text"/>	PO Number	<input type="text"/>
State	<input type="text"/>	Reference	<input type="text"/>
Postal Code	<input type="text"/>	Notes	<input type="text"/>

ShowerBuddy SB2 TubBuddy

The TubBuddy (SB2) transfer system will transform your existing bathtub or curbed shower into an accessible, barrier-free shower in minutes. It allows the user to enjoy showering again without the high cost, delay, and stress of a full bathroom remodel. The SB2 alleviates heavy lifting and multiple transfers. This chair is separated into three sections: rolling base (shower chair), bridge, and tub base. The rolling base can double as a bedside commode. It also rolls around the house like a wheelchair, allowing users to transfer comfortably and roll over most toilets to be used as an over-toilet commode. When toileting is complete, the rolling base and tub base are connected using the bridge, allowing the unit to slide over the tub lip or curbed shower. Once over the tub/shower, the rolling base and bridge can be disconnected and removed so that the shower curtain can either be closed or allow the caregiver space to move closer to the patient. A custom bridge may be required if the toilet is too close or too far from the tub. The SB2 features a secondary safety mechanism that prevents the chair from rolling off the tracks and causing injury. The SB2 chair comes equipped with a standard front-opening cushion, footrests, flip-back, lockable and removable armrests, 5" lockable casters, a removable commode bucket with lid, chest belt, lap belt, and height adjustability. The TubBuddy (SB2) is a caregiver-assisted product; it is not intended for independent use.

The SB2 weighs 59.5 lbs with a weight capacity of 350 lbs.

1. Primary	Description	MSRP	Quantity
SB2	ShowerBuddy SB2 TubBuddy	\$3,695.00	1

2. Accessories		Description	MSRP	Quantity
SBC-5		Custom Padded Back Cushion - SB1, SB2, SB2T, SB3T & SB6	\$685.00	
B2		Lateral Kit (Right)	\$460.00	
B10		Large Headrest Kit (non tilting chairs)	\$395.00	
B8		Leg Adductor Kit	\$383.00	
B6		Padded Back Kit (SB1, SB2, SB2T,SB3T)	\$441.00	
B4		BodyPoint Chest Belt Kit (SB1, SB2, SB2T, SB3T)	\$242.50	
B5		BodyPoint Lap Belt Kit (SB1, SB2, SB2T, SB3T)	\$242.50	
SBLT-B		SolutionBased Belt for Laterals	\$68.00	
SBACD100		ShowerBuddy Channel Armrests (Each) - SB1, SB2, SB2T, SB3T	\$165.00	
TBR		ShowerBuddy Tensioned Backrest Adjustable - SB1, SB2, SB2T, SB3T	\$85.00	
B1		Lateral Kit (Left)	\$460.00	
BP-LH-S		ShowerBuddy BodyPoint Leg Harness - Small	\$110.00	
BP-LH-M		ShowerBuddy BodyPoint Leg Harness - Medium	\$110.00	
B9		Neckrest Assembly	\$230.00	

3. Other		Description	MSRP	Quantity
LCB		Custom Size Bridge (Labor)	\$300.00	
SPH		ShowerBuddy Straight Push Handle	\$96.00	
RBE24		ShowerBuddy Rolling Base Extension Rail Pair	\$220.00	
49		ShowerBuddy Suction Cup (Each) - SB1, SB2, SB2T	\$30.00	
TBE24		ShowerBuddy Tub Base Extension Rail Pair - SB2, SB2T	\$220.00	

4. Pediatrics		Description	MSRP	Quantity
P1		ShowerBuddy P1 Petite Converter Kit	\$1,195.00	
P2		ShowerBuddy P2 Petite Converter Kit	\$1,295.00	
B-16		Lateral Kit for Pediatric Converter Kit (P1 & P2)	\$1,200.00	
SBC-3		Custom Padded Back Cushion - P1 & P2	\$685.00	

5. Seating options		Description	MSRP	Quantity
B15		22" Wide Bariatric Kit (SB1, SB2, SB2T, SB3T)	\$1,305.00	
B7		Roho Cushion Kit	\$900.00	
SCCF-2		ShowerBuddy Seat Cushion Closed Front (SB1-SB3T)	\$195.00	
SBC-2		ShowerBuddy Custom Visco Seat Cushion	\$745.00	

