

ShowerBuddy SB2T TubBuddy with Tilt Order Form

For assistance call (877) 769-2833



Company	<input type="text"/>	Name	<input type="text"/>
Address 1	<input type="text"/>	Phone	<input type="text"/>
Address 2	<input type="text"/>	Email *	<input type="text"/>
City	<input type="text"/>	PO Number	<input type="text"/>
State	<input type="text"/>	Reference	<input type="text"/>
Postal Code	<input type="text"/>	Notes	<input type="text"/>

ShowerBuddy SB2T TubBuddy with Tilt

The TubBuddy with Tilt (SB2T) transfer system will transform your existing bathtub or curbed shower into an accessible, barrier-free bathroom in minutes. This chair allows the user to enjoy showering again without the high cost, delay, and stress of a complete bathroom remodel. The TubBuddy SB2T alleviates multiple transfers and heavy lifting. The SB2T is composed of three sections: rolling base (shower chair), bridge, and tub base. The rolling base can double as a bedside commode. It also rolls around the house like a wheelchair, allowing users easy transfer. Once in the bathroom, the SB2T can roll over most toilets to be used as an over-toilet commode. When toileting is complete, the rolling base and shower base are connected using the bridge, allowing the chair to slide in a tilted position with legs raised and ready to move over the tub. The 30-degree infinite tilt can also be used for pressure relief, increased comfort, and easy-reach bathing. This feature is unique to SB2T and to the marketplace, and it cannot be found on any other shower transfer system. Once over the tub/shower, the rolling base and bridge can be disconnected and removed so that the shower curtain can be closed for privacy or allow the caregiver space to move closer to the patient. A custom bridge may be required if the toilet is too close or too far from the tub. The SB2T features a secondary safety mechanism to prevent the chair from sliding off the tracks and causing injury. The SB2T chair comes equipped with a standard front-opening cushion, neck rest, swing-away and removable footrests with calf strap, flip-back/lockable and removable armrests, 5" lockable casters, a removable commode bucket with lid, chest belt, lap belt, and height adjustability. The TubBuddy with Tilt (SB2T) is a caregiver-assisted product; it is not intended for independent use.

The SB2T weighs 66.62 lbs with a weight capacity of 350 lbs.

1. Primary	Description	MSRP	Quantity
SB2T	ShowerBuddy SB2T TubBuddy with Tilt	\$4,095.00	1

2. Accessories		Description	MSRP	Quantity
SBC-5		Custom Padded Back Cushion - SB1, SB2, SB2T, SB3T & SB6	\$685.00	
B2		Lateral Kit (Right)	\$460.00	
B3		Large HeadRest Kit (SB2T, SB3T)	\$260.00	
B-17		ShowerBuddy 18" Footplate with Custom Padded Cushion Kit	\$500.00	
B8		Leg Adductor Kit	\$383.00	
B6		Padded Back Kit (SB1, SB2, SB2T,SB3T)	\$441.00	
B4		BodyPoint Chest Belt Kit (SB1, SB2, SB2T, SB3T)	\$242.50	
B5		BodyPoint Lap Belt Kit (SB1, SB2, SB2T, SB3T)	\$242.50	
SBLT-B		SolutionBased Belt for Laterals	\$68.00	
SBACD100		ShowerBuddy Channel Armrests (Each) - SB1, SB2, SB2T, SB3T	\$165.00	
TBR		ShowerBuddy Tensioned Backrest Adjustable - SB1, SB2, SB2T, SB3T	\$85.00	
B1		Lateral Kit (Left)	\$460.00	
BP-LH-S		ShowerBuddy BodyPoint Leg Harness - Small	\$110.00	
BP-LH-M		ShowerBuddy BodyPoint Leg Harness - Medium	\$110.00	
FFP-18		ShowerBuddy Flat Footplate 18" - SB2T, SB3T, SB6	\$150.00	
BPB-CS		Bodypoint Calf Strap	\$125.00	

3. Other		Description	MSRP	Quantity
LCB		Custom Size Bridge (Labor)	\$300.00	
SPH		ShowerBuddy Straight Push Handle	\$96.00	
RBE24		ShowerBuddy Rolling Base Extension Rail Pair	\$220.00	
49		ShowerBuddy Suction Cup (Each) - SB1, SB2, SB2T	\$30.00	
TBE24		ShowerBuddy Tub Base Extension Rail Pair - SB2, SB2T	\$220.00	

4. Pediatrics		Description	MSRP	Quantity
RW-M-B		Blue Rifton Wave (Medium) with Chest Strap with Laterals Head Blocks and Calf Rest	\$1,210.00	
P1		ShowerBuddy P1 Petite Converter Kit	\$1,195.00	
P2		ShowerBuddy P2 Petite Converter Kit	\$1,295.00	
RW-S-B		Blue Rifton Wave (Small) with Chest Strap with Laterals Head Blocks and Calf Rest	\$1,165.00	
Rifton Hardware Kit		Rifton Hardware Kit	\$575.00	
B-16		Lateral Kit for Pediatric Converter Kit (P1 & P2)	\$1,200.00	
RW-L-B		Blue Rifton Wave (Large) with Chest Strap with Laterals Head Blocks and Calf Rest	\$1,340.00	
SBC-3		Custom Padded Back Cushion - P1 & P2	\$685.00	

5. Seating options		Description	MSRP	Quantity
B15		22" Wide Bariatric Kit (SB1, SB2, SB2T, SB3T)	\$1,305.00	
B7		Roho Cushion Kit	\$900.00	
SCCF-2		ShowerBuddy Seat Cushion Closed Front (SB1-SB3T)	\$195.00	
SBC-2		ShowerBuddy Custom Visco Seat Cushion	\$745.00	