

SB2 Order Form

For assistance call (877) 769-2833



Company	<input type="text"/>	Name	<input type="text"/>
Address 1	<input type="text"/>	Phone	<input type="text"/>
Address 2	<input type="text"/>	Email *	<input type="text"/>
City	<input type="text"/>	PO Number	<input type="text"/>
State	<input type="text"/>	Reference	<input type="text"/>
Postal Code	<input type="text"/>	Notes	<input type="text"/>

ShowerBuddy SB2T TubBuddy with Tilt

The TubBuddy with Tilt (SB2T) transfer system will transform your existing bathtub or curbed shower into an accessible, barrier-free bathroom in minutes. This chair allows the user to enjoy showering again without the high cost, delay, and stress of a complete bathroom remodel. The TubBuddy SB2T alleviates multiple transfers and heavy lifting. The SB2T is composed of three sections: rolling base (shower chair), bridge, and tub base. The rolling base can double as a bedside commode. It also rolls around the house like a wheelchair, allowing users easy transfer. Once in the bathroom, the SB2T can roll over most toilets to be used as an over-toilet commode. When toileting is complete, the rolling base and shower base are connected using the bridge, allowing the chair to slide in a tilted position with legs raised and ready to move over the tub. The 30-degree infinite tilt can also be used for pressure relief, increased comfort, and easy-reach bathing. This feature is unique to SB2T and to the marketplace, and it cannot be found on any other shower transfer system. Once over the tub/shower, the rolling base and bridge can be disconnected and removed so that the shower curtain can be closed for privacy or allow the caregiver space to move closer to the patient. A custom bridge may be required if the toilet is too close or too far from the tub. The SB2T features a secondary safety mechanism to prevent the chair from sliding off the tracks and causing injury. The SB2T chair comes equipped with a standard front-opening cushion, neck rest, swing-away and removable footrests with calf strap, flip-back/lockable and removable armrests, 5" lockable casters, a removable commode bucket with lid, chest belt, lap belt, and height adjustability. The TubBuddy with Tilt (SB2T) is a caregiver-assisted product; it is not intended for independent use.

The SB2T weighs 66.62 lbs with a weight capacity of 350 lbs.

1. Primary	Description	MSRP	Quantity
SB2T	ShowerBuddy SB2T TubBuddy with Tilt	\$4,095.00	1

2. Accessories	Description	MSRP	Quantity
BP-B	ShowerBuddy Bodypoint Belt	\$181.50	
TBR	ShowerBuddy Tensioned Backrest Adjustable - SB1, SB2, SB2T, SB3T	\$85.00	
BPB-CS	ShowerBuddy Bodypoint Calf Strap	\$125.00	
34-22	ShowerBuddy 22" Arm Rest Right (Each) - SB1, SB2, SB2T & SB3t	\$90.00	
FFP-18	Flat Foot Plate 18" - SB1, SB2, SB2T & SB3T	\$60.00	
33-22	ShowerBuddy 22" Arm Rest Left (Each) - SB1, SB2, SB2T, SB3T	\$90.00	
PBC	ShowerBuddy Solid Padded Backrest	\$195.00	
BP-B-PTC32	ShowerBuddy BP Lap Belt Tube Clamp 32mm (each)	\$16.00	
SBLT-L	ShowerBuddy Lateral (LEFT) - SB1, SB2, SB2T, SB3T	\$395.00	
SBLT-R	ShowerBuddy Lateral (Right) - SB1, SB2, SB2T & SB3T	\$395.00	
70A	ShowerBuddy Neckrest - SB2T, SB3T	\$95.00	
large-headrest-kit	Large headrest kit	\$260.00	
SBLAL	ShowerBuddy Leg Adductor (Each) - SB1, SB2, SB2T, SB3T	\$318.00	
BP-B-QR	ShowerBuddy BP Belt Quick Release Knobs (Pair)	\$29.00	
TBE24	ShowerBuddy Tub Base Extension Rail Pair - SB2, SB2T	\$220.00	
BP-B-PTC25	ShowerBuddy BP Chest Belt Tube Clamp 25mm (each)	\$16.00	
SBACD100	ShowerBuddy Channel Armrests (Each) - SB1, SB2, SB2T, SB3T	\$125.00	
SB-SS-100	SB Lateral/Side Support Pad (each)	\$65.00	

3. Other	Description	MSRP	Quantity
LCB	Custom Size Bridge (Labor)	\$300.00	
RBE24	ShowerBuddy Rolling Base Extension Rail Pair	\$220.00	
SBRD	ShowerBuddy Redirectional / Stabilizer	\$395.00	
49	ShowerBuddy Suction Cup - SB1, SB2T	\$30.00	

4. Pediatrics	Description	MSRP	Quantity
P1	ShowerBuddy P1 Petite Converter	\$1,195.00	
P2	ShowerBuddy P2 Petite Converter Kit	\$1,295.00	
PSA	ShowerBuddy Pediatric Seat Adaptor	\$995.00	

5. Seating options		Description	MSRP	Quantity
SCSF-22		ShowerBuddy 22" Soft Cushion Seat Foam - SB1, SB2, SB2T & SB3T	\$90.00	
SCO		ShowerBuddy Seat Cushion Overlay (SB1-SB3T)	\$155.00	
SCP1-500-166		ShowerBuddy Pediatric Seat Cushion 6" Commode Opening (SCP1)	\$195.00	
SCP2-500-167		ShowerBuddy Pediatric Seat Cushion 6.6" Commode Opening (SCP2)	\$195.00	
SCSF-05P		ShowerBuddy Soft Foam Cushion w/ Petite Opening (SB1, SB2, SB2T & SB3T)	\$165.00	
SCSF		ShowerBuddy Soft Cushion (SB3t, SB2t, SB2, SB1)	\$140.00	
30		ShowerBuddy Cushion - SB1, SB2, SB2T, SB3T	\$170.00	
SCCF-2		ShowerBuddy Seat Cushion Closed Front (SB1-SB3T)	\$195.00	
PM-1		ShowerBuddy Roho Shower Commode Cushion	\$675.00	