

Nextkid Reha Booster US- Jan. 2026 US Retail prices

Date	Buyer	PO#	Mark for
Address			
Contact	Phone	Fax	

Booster Mode Seat + Accessories - Children weighing between 40 and 110.2 lbs. and whose height is between 43.4 and 59 inches



Dimensions & weights	NEXTKID Reha US
Seat depth (behind crotch strap)	23 cm / 9.1"
Seat depth	37 cm / 14.5"
Seat width	22 - 29 cm / 8.7 - 11.4"
Back height (seat inside)	55 - 75 cm - 21.7 - 29.5"
Shoulder height (harness guide)	32 - 52 cm - 12.6 - 20.4"
Shoulder width	36 - 42 cm / 14.2 - 16.5"
Overall dimensions (w x h x d):	44 x 62 x 48 cm / 17.3 x 24.4 x 18.9"
Weight of seat	8.5 kg / 18.7 lbs
User height (min. - max.)	110.1 - 150 cm / 43.4 - 59"
User weight (min. - max.)	18.2 - 50 kg / 40 - 110.2 lbs



	Model	Price	HCPCS	Description
Nextkid Reha	£ 523	\$3020	E1399	Seat + back w. height-adjustable headrest, 5-point positioning belt w. shoulder pads, chest clip, small crotch pad, adjustable crotch strap, sets of removable head/ trunk & hip support, adjustable side wings, side impact protection, base plate, color black

	Option 1 - Swivel & tilt base w.lockout/right or left door HCPCS	Unit	Item no.	Price
Right	Swivel & tilt base(10°) incl. footrest adapter & LATCH*	E1399	£ 535-0	\$1210
Or	Left turning swivel & tilt base(10°)**	E1399	£ 535-1	\$1210
	Stabilizing bow (if vehicle has no latch)	E1399	£ 536	\$187
Either	Footrest, foldable (depth 6.2") (angle adjust.)	E1399	£ 537	\$742
Or	Front Stand (stability option if no footrest is used)	E1399	£ 538	\$480
Either	Abduction block, rigid, adjustable	E1399	£ 539	\$340
Or	Seat depth extension, 5"	E1399	£ 540	\$321



535 Swivel & tilt base w. lockout

**** IMPORTANT: Swivel direction cannot be changed!**

Option 2 - Footrest adapter only - no swivel & tilt

	Adapter for footrest incl. LATCH*	E1399	Each	£ 532	\$745
	Stabilizing bow (if vehicle has no latch)	E1399	Each	£ 536	\$187
Either	Footrest, foldable (depth 6.2") (angle adjust.)	E1399	Each	£ 537	\$742
Or	Front Stand (stability option if no footrest is used)	E1399	Each	£ 538	\$480
Either	Abduction block, rigid, adjustable	E1399	Each	£ 539	\$340
Or	Seat depth extension, 5"	E1399	Each	£ 540	\$321

*** 535 & 532 latch connectors are for stability options only and not a safety device. Nextkid booster must always be secured with vehicle's lap and shoulder belt, even when unoccupied.**



532 Adapter incl. LATCH + 537 footrest



532 Adapter incl. LATCH + 538 front stand

Nextkid Reha Booster US – 2026 US Retail prices

Option 3 - no footrest adapter 532 - no swivel & tilt 535

Abduction block, soft padded*no combination with 539	E1399	Each	£ 358	\$154
Either Abduction block, rigid, adjustable	E1399	Each	£ 539	\$340
Or Seat depth extension, 5"	E1399	Each	£ 540	\$321



358 abduction block, soft padded
+ 540 seat depth extension 5"



539 abduction block, rigid, adjustable



540 seat depth extension, 5"

This is a Universal Restraint System for disabled children. It is approved in accordance with Federal Motor Vehicle Safety Standard 213 for general use in vehicles and it will fit most but not all cars. A correct fit is likely if the vehicle manufacturer has declared in the vehicle handbook that the vehicle is capable of accepting a Universal Restraint System for this age group. Only suitable for use in **vehicles with 3-point** (lap & shoulder) seat belts. This is a seat for children weighing between 40 – 110 lbs. and whose height is between 43.4 and 57 inches in **Booster Mode**. Suitable only for **forward-facing** installation. **THIS booster** seat is **NOT CERTIFIED** for use in an aircraft.

General Accessories

Anterior support, table, black	E1399	Each	£ 528	\$616
Seat wedge, below (+10°), black **	E1399	Each	£ 810	\$205
Cushion, leveling & protection	E1399	Each	£ 513	\$158
Abduction block, soft padded*no combination with 539	E1399	Each	£ 358	\$154



528 anterior support, table



810 seat wedge below (+10°)



513 Cushion, leveling & protection

**** Seat with lower anchors 523 & seat wedge below 810 cannot be used in some vehicles. Seat wedge below 810 cannot be used with swivle and tilt base 535/ footrest adapter with latch 532.**

Nextkid Reha **AND** all accessories exceed US FMVSS 213 standards

Specs & pricing are subject to change without notice. HCPCS are suggestions only with no guarantee of reimbursement.
ExoMotion, LLC dba Thomashilfen | www.thomashilfen.us Effective Date January 1st 2026