

Defender Reha US Car Seat – Jan. 2026 US Retail prices

Date	Buyer	PO#	Mark for
Address			
Contact	Phone	Fax	

Step 1 – Seat



Model	Price	HCPCS	Description
Seat £ 886	\$993	E1399	Seat and back with height-adjustable headrest, LATCH equipped, 5-point safety/positioning belt with shoulder pads, chest clip, small crotch pad, set of removable head, trunk and hip supports, cup holder, color black-charcoal, 2 years warranty



Harness Mode only



Step 2 – General Options

General Accessories	HCPCS Unit	Item no.	Price
Footrest, short (depth 4.7")* (angle adjust.)	E1399 Each	£ 806	\$356
Footrest, long (depth 7.1")* (angle adjust.)	E1399 Each	£ 805	\$356
Footrest adapter + Base plate 89036*	E1399 Each	£ 821-0 + 89036	\$388
Seat wedge, below, black	E1399 Each	£ 810	\$205
Crotch pad, large, black**	E1399 Each	£ 832	\$83
Anterior support, table, black	E1399 Each	£ 831	\$325
Abduction block, padded, black	E1399 Each	£ 833	\$126
Footrest pad (for 806), height 4", black	E1399 Each	☐ 834	\$104
Front stand	E1399 Each	☐ 817	\$209

Technical Specs Defender Reha

Child weight	26.5 - 120lbs
Seat weight	22 lbs
Seat depth	13.5"
Seat width	8.5-12" ***
Seat to head	23.5-30.5" ****
Shoulder height	10-19.5" *****
Lower leg length	(7")# 11-17"

#7"minimum with footrest pad 834

This is a Universal Restraint System for disabled children. It is approved in accordance with Federal Motor Vehicle Safety Standard 213 for general use in vehicles and it will fit most but not all cars. A correct fit is likely if the vehicle manufacturer has declared in the vehicle handbook that the vehicle is capable of accepting a Universal Restraint System for this age group. Only suitable for use in **vehicles with 3-point** (lap & shoulder) seat belts. This is a seat for children weighing between 26.5 and 65 lbs., whose height is between 27 and 57 inches and who are over 12 month old in **Harness Mode** or 40 – 120 lbs., and whose height is between 43.4 and 57 inches in **Booster Mode**. Suitable only for **forward-facing** installation.

* The footrest adapter 821 & Base plate 89036 must be ordered to use footrest

** We recommend using larger crotch pad only in booster mode

*** To reach max. seat width remove additional hip guide pads from additional seat pad

**** To reach max. back height remove additional seat pad and foam pad under seat cushion

***** To reach max. shoulder height remove additional seat pad and foam pad under seat cushion



Footrest + Footrest pad 4"



Table



Seat wedge below



Abduction block



Crotch pad

Specs & pricing are subject to change without notice. HCPCS are suggestions only with no guarantee of reimbursement. ExoMotion, LLC dba Thomashilfen | www.thomashilfen.us Effective Date January 1st 2026