

EASyS Advantage – 2026 US Retail prices

info@thomashilfen.com | Phone 866.870.2122 | Fax 866.870.0801

Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	

Step 1 – Seat Sizes (Required options are indicated below by † symbol)

6204 EASyS 6404 EASyS



	seat size 1		seat size 2		HCPCS	Description
†Seat	660961	\$441	670961	\$441	E2292	Solid contoured, integrated seat/back, 180° recline, growth adjustable; seat depth/ width adjustable from underneath, articulated foot plate height and angle adjustable; knee angfront/rear facing seat, easy to change; ventilated/washable upholstery ; knee angle adjustment with release lever; cozy comfortable seat cushion, BSS-Balance Safety System; calf strap; 5 different color options
†Back	660962	\$544	670962	\$544	E2293	
†Recline	660963	\$457	670963	\$457	E1014	
†Hip sup.	660964	\$252	670964	\$252	E0956	
†Calf panel	6635	\$69	6735	\$69	E0995	
†Foot plate, flip up	660965	\$105	670965	\$105	K0040	
†Legrest, elevating	660966	\$210	670966	\$210	E0990	

Seat Size Chart	Lower Leg	Seat Depth	Seat width	Seat to Head	Max Height	Hip Angle
Size 1	7 – 11½"	7½ – 11¾"	7½ – 11¾"	16½ – 26¾"	49"	90° - 180°
Size 2	8½ – 14½"	10½ – 15"	10 – 14½"	21½ – 29"	60"	90° - 180°

†Upholstery colors	Red/ Grey	Blue/ Grey	Green/ Grey	Magenta/ Grey	Grey/ Grey
Size 1	<input type="checkbox"/> 6247	<input type="checkbox"/> 6248	<input type="checkbox"/> 6249	<input type="checkbox"/> 6246	<input type="checkbox"/> 6245
Size 2	<input type="checkbox"/> 6447	<input type="checkbox"/> 6448	<input type="checkbox"/> 6449	<input type="checkbox"/> 6446	<input type="checkbox"/> 6445

Step 2a – Outdoor Frame Option (Required options are indicated below by † symbol)

Advantage-chassis - Folding aluminum frame, tilt-in-space, back frame suspension, face-to-face option, 0° safety lock out, cargo basket, unfold lock



	HCPCS	SIZE 1	SIZE 2
Advantage-chassis	E1234	<input type="checkbox"/> 62023 \$2830	<input type="checkbox"/> 64023 \$2830
†Castor pin lock, pair	K0073	67830 \$101	67830 \$101
†Bus tie-down loop kit (*)	K0108	67839 \$221	67839 \$221
Pneumatic wheels, full set	K0070	<input type="checkbox"/> 66461 \$325	<input type="checkbox"/> 66461 \$325
Frame padding, black	K0108	<input type="checkbox"/> 6622-6 \$132	<input type="checkbox"/> 6622-6 \$132
Frame padding, red/ black	K0108	<input type="checkbox"/> 6622-7 \$132	<input type="checkbox"/> 6622-7 \$132
Tip assist/ curb tipper	K0108	<input type="checkbox"/> 67502 \$111	<input type="checkbox"/> 67502 \$111
Under seat storage basket	K0108	NC	NC

EASyS Advantage Technical Specs	SIZE 1	SIZE 2
Seat capacity	77.2 lbs	88.2 lbs
Frame capacity	110.2 lbs	121.3 lbs
Unit weight	41 lbs	48 lbs
Max seat tilt	-10° to +35°	
Frame suspension	Yes	
Upholstery	Trevira-CS	
Foot angle	78° - 100°	
Size of foot rest	11" x 7.1"	12.6" x 8.3"
Floor to handle	28.2 – 46.1"	28.2 – 46.1"
Floor to seat	19"	19.5"
Max width	23.4"	27.2"
Fold length	33.5"	34"
Fold height	23.6"	27.2"
Brakes	foot operated	
Seat direction	front & rear facing	
Front wheels	7.5" solid	
Rear wheels	11.4" solid	



* Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

Step 2b – Indoor Frame Option (Required options are indicated below by † symbol)

Q-chassis - hi-low indoor base/feeder/therapy frame, height/ tilt adjustable, pedal controlled gas springs, lockable rear wheels



	HCPCS	SIZE 1 & 2	
Q-chassis indoor hi/low, pedal controlled	E1399	<input type="checkbox"/> 6666	\$2308
†EASyS Adapter for Q-chassis	K0108	<input type="checkbox"/> 62951	\$137
Push bar for Q-chassis	K0108	<input type="checkbox"/> 6660	\$346

Q-chassis Technical Specs **

Capacity	88 lbs
Seat tilt	-7 to +31°
Floor to seat	11 - 28"
Seat direction	front facing

Step 3 – Positioning & General Options (Required options are indicated below by † symbol)

Lower Body Positioning	HCPCS	Unit	Small	Medium	Large	Price
Foot straps with cam buckle	E0952	Pair		<input type="checkbox"/> 6648		\$37
Ankle Huggers, Bodypoint®	K0108	Pair		<input type="checkbox"/> 213151		\$221
Footrest pad, height 1¾"	K0108	Each	<input type="checkbox"/> 10622		<input type="checkbox"/> 10023	\$104

Mid Body Positioning

Abduction block, padded	E0957	Each	<input type="checkbox"/> 6620		<input type="checkbox"/> 6720-6	\$179	
Hip supports	2 point, padded, SR buckle	E0978	Each	<input type="checkbox"/> 8420		<input type="checkbox"/> 8470	\$116
	4 point, padded, SR buckle	E0978	Each	<input type="checkbox"/> 8421		<input type="checkbox"/> 8471	\$158
	Soft pelvic harness/ leg harness	E0978	Each	<input type="checkbox"/> 6614	<input type="checkbox"/> 6714	<input type="checkbox"/> 213152	\$221
Crotch strap (fits 6629,8420,8470,8421,8471)	E0978	Each		<input type="checkbox"/> 6341		\$37	
Tray, height adjustable, clear with picture window	E0950	Each	<input type="checkbox"/> 66281		<input type="checkbox"/> 67281	\$237	
Mounting hardware tray - †required when ordering 6628/6728	E1034	Each	<input type="checkbox"/> 66282		<input type="checkbox"/> 67282	\$258	
Arm rest, height, depth and angle adjustable	E0973	Pair			<input type="checkbox"/> 8478	\$394	

Upper Body Positioning

Chest harness supports	Butterfly style, neoprene	E0960	Each	<input type="checkbox"/> 6385		<input type="checkbox"/> 6585	\$221
	Vest style	E0960	Each	<input type="checkbox"/> 6631	<input type="checkbox"/> 6632	<input type="checkbox"/> 6731-0	\$221
	H-style, posture support, small, Bodypoint®	E0960	Each		<input type="checkbox"/> 214806		\$221
	PivotFit shoulder harness, dynamic, Bodypoint®	E0960	Each	<input type="checkbox"/> 214270	<input type="checkbox"/> 214250		\$242
	5-point harness	E0960	Each		<input type="checkbox"/> 8480		\$153
Trunk supports	Rigid lateral, left & right, height adjustable	E0956	Pair	<input type="checkbox"/> 6634		<input type="checkbox"/> 6734	\$283
	Flexible lateral, l&r w/ chest strap, height adj.	E0956	Pair		<input type="checkbox"/> 6633		\$310
Head supports	Contoured pillow	E0955	Each		<input type="checkbox"/> 6824	<input type="checkbox"/> 68370	\$221
	Rigid lateral, left & right, height adjustable	E0955	Pair	<input type="checkbox"/> 6624		<input type="checkbox"/> 6724	\$263

General Accessories

Grab bar, padded, black	K0108	Each	<input type="checkbox"/> 6621		<input type="checkbox"/> 6721-5	\$168	
Footrest tube upholstery, grey	K0108	Pair		<input type="checkbox"/> 6230		\$37	
Vent tray	E1029	Each	<input type="checkbox"/> 6294		<input type="checkbox"/> 6494	\$504	
Medical necessities pack, black	K0108	Each		<input type="checkbox"/> 6639		\$95	
Oxygen tank holder, for D/E-tanks	E2208	Each		<input type="checkbox"/> 6297		\$221	
Straps	Secure folded seat/frame	K0108	Each		<input type="checkbox"/> 667011	\$37	
	Wrist "runaway" tether	K0108	Each		<input type="checkbox"/> 667111	\$37	
Weather protection	Canopy, sun/rain with window, black	K0108	Each		<input type="checkbox"/> 663761	\$342	
	Lower body windbreaker, black for 663761	K0108	Each		<input type="checkbox"/> 661852	\$104	
	Canopy, rain only, clear	K0108	Each		<input type="checkbox"/> 6638	\$258	
	Rain shield, clear for canopy 663761/ 637811	K0108	Each		<input type="checkbox"/> 6739	\$132	
	Rain cape, grey	K0108	Each	<input type="checkbox"/> 63404		<input type="checkbox"/> 6540	\$172
	Sleeping bag, summer weight, black	K0108	Each	<input type="checkbox"/> 6619		<input type="checkbox"/> 6719	\$314
	Sleeping bag, winter weight, black	K0108	Each	<input type="checkbox"/> 6618-5		<input type="checkbox"/> 6718	\$314
Umbrella, angle adjustable	K0108	Each		<input type="checkbox"/> 6594		\$100	

*Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

Specs & pricing are subject to change without notice. HCPCS are suggestions only with no guarantee of reimbursement. ExoMotion LLC dba Thomashilfen | www.thomashilfen.us Effective Date January 1st 2026