

# Hercules Small –2026 US Retail prices

info@thomashilfen.com | Phone 866.870.2122 www.Thomashilfen.us

Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	

## Step 1 – Booster Seat



Booster seat

Model	Price	HCPCS	Description
<b>Right Swivel</b> £ 471-0 \$4185		E1399	Booster seat with tilt function (up to 10° adjustable), height adjustable backrest, height & width adj. lateral head supports, back angle adjustment, side impact protectors, width adjustable hip & leg guides, seat depth adjustment, 5-point positioning harness with shoulder pads, crotch pad, chest clip, mounting options for accessories, removable & washable covers - swivel function to <b>RIGHT</b> vehicle door, incl. swivel lock†
<b>Left Swivel</b> £ 471-1 \$4185		E1399	Booster Seat as above - swivel function to <b>LEFT</b> vehicle door† <b>By Special order only!</b>

Swivel function\* with lock out

Seat tilt up to 10°



\*Swivel direction **CANNOT** be changed!

Hip & leg guide adjustment



This is a Universal Restraint System for disabled children. It is approved in accordance with Federal Motor Vehicle Safety Standard 213 for general use in vehicles and it will fit most but not all cars. A correct fit is likely if the vehicle manufacturer has declared in the vehicle handbook that the vehicle is capable of accepting a Universal Restraint System. Only suitable for use in vehicles with 3-point (lap & shoulder) seat belts.

**This is a booster seat for children, young persons and adults weighing between 40 and 110 lbs. (18.1 and 50 kg) AND whose height is between 43.4 and 59 inches (110.1 cm to 150cm). Suitable only for forward-facing installation. THIS booster seat is NOT CERTIFIED for use in an aircraft.**



## Step 2 – † Required Options: **ALWAYS** use LATCH connectors (479) or stabilizing bow (421) for stability purposes!

	HCPCS	Unit	Item no.	Price
† Lower anchors connection	E1399	Each	£ 479	\$178
† Stabilizing bow (if car has no latch)	E1399	Each	£ 421	\$132



479



421

### Accessories for additional stability purposes

Front stand (min. 15.75" footwell height)	E1399	Each	£ 480	\$415
Cushion, leveling (for shaped seats)	E1399	Each	£ 513	\$158



480



513

### Accessories for i ddYf VcXm

Head pad, width reducer	E1399	Each	£ 494	\$103
Shoulder pads, width reducer	E1399	Pair	£ 483	\$238
Upper belt guide, additional*	E1399	Each	£ 492	\$84
Chest clip**	E1399	Each		NC



494 Head pad, width reducer



483 Shoulder pads



492 Add. upper belt guide\*

Technical Specs Hercules Small	
Seat depth	31 - 39 cm / 12.2 - 15.35"
Depth behind crotch strap	20 - 33 cm / 7.9 - 13"
Seat width	21 - 33 cm / 8.7 - 13"
Seat to head	57 - 80 cm / 22.4 - 31.5"
Back width (seat)	45 cm / 17.7"
Shoulder height	36 - 60 cm / 14.2 - 23.6"
Shoulder width	45 cm / 17.7"
Shoulder width with shoulder pad	37 cm / 14.6"
Back angle (max.)	90° - 145°
Seat tilt-in-space (with tilting base 10°)	0° / 2.5° / 5° / 7.5° / 10°
Lower leg length (with foot rest)	16 - 42 cm / 6.3 - 16.5"
Height lateral trunk supports (upper edge)	28 - 34 cm / 11 - 13.4"
Distance lateral trunk supports	20 - 33 cm / 7.9 - 13"
Max. overall dimensions (w x h x d)	55 x 85 x 34 cm / 21.7 x 33.5 x 13.4"
Min. folded dimensions (w x h x d)	55 x 23 x 60 cm / 21.7 x 9.1 x 23.6"
Weight of seat	11.6 kg / 25.6 lbs
Body size user	110.1 - 150 cm / 43.4 - 59"
User weight	18.2 - 50 kg / 40 - 110.2 lbs

**Accessories for trunk support & mid body**

Lateral trunk supports, swing-away, w. brackets	E1399	Pair	£ 487	\$835
Lateral trunk supports, rigid, w. brackets	E1399	Pair	£ 486	\$283
Armrests ( <b>table 491 can only be used with 489 armrests</b> )	E1399	Pair	£ 489	\$394
Anterior support, table ( <b>required for usage: 498 armrests</b> )	E1399	Pair	£ 491	\$657



487 Lateral trunk supports, swing-away



486 Lateral trunk supports, rigid



489 Armrests



491 Anterior support, table

**Accessories for lower body and legs**

Abduction block	E1399	Each	£ 488	\$235
Footrest, foldable, angle adjust. (min. 15.75" footwell height)	E1399	Each	£ 481	\$588
Incontinence seat pad	E1399	Each	£ 497	\$94
Front cover	E1399	Each	£ 433	\$58
Crotch pad with magnetic lock	E1399	Each	£ 498	\$247
Extra key for magnetic lock	E1399	Each	£ E66331	\$25



488 Abduction block




481 Footrest, foldable, angle adjustable



497 Incontinence seat pad



433 Front cover



498 Crotch pad with magnetic lock

\* Swivel direction **CANNOT** be changed!

**\*\*CAUTION**  
**ALWAYS** keep clip on chest height level.

Hercules Small **AND** accessories exceed US FMVSS 213 standards

Specs & pricing are subject to change without notice. HCPCS are suggestions only with no guarantee of reimbursement.  
ExoMotion, LLC dba Thomashilfen| www.thomashilfen.us Effective Date January 1st 2026