

## Hercules Prime –2026 US Retail prices

Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	

### Step 1 – Booster Seat



Booster seat

Model	Price	HCPCS	Description
<b>Right Swivel</b> £ 461-0 \$4017		E1399	Booster seat with tilt function (up to 10° adjustable), height adjustable backrest, height adj. head cushion, back angle adjustment for proper installation, seat depth adjustment, 5-point positioning harness with shoulder pads, crotch pad, chest clip, mounting options for accessories, removable covers, machine washable, seat & accessories tested acc. FMVSS 213 - Swivel function to <b>RIGHT</b> vehicle door, incl. swivel lock out†
<b>Left Swivel</b> £ 461-1 \$4017		E1399	- Swivel function to <b>LEFT</b> vehicle door† <b>By Special order only!</b>



Swivel function\* with lock out

Seat tilt up to 10°



\*Swivel direction **CANNOT** be changed!

This is a Universal Restraint System for disabled children, young persons and adults. It is approved in accordance with Federal Motor Vehicle Safety Standard 213 for general use in vehicles and it will fit most but not all cars. A correct fit is likely if the vehicle manufacturer has declared in the vehicle handbook that the vehicle is capable of accepting a Universal Restraint System. Only suitable for use in vehicles with 3-point (lap & shoulder) seat belts.

**This is a booster seat for children, young persons and adults weighing between 40 and 176.4 lbs. (18.1 and 80 kg) AND whose height is between 51.1 and 70.9 inches (130 cm to 180cm). Suitable only for forward-facing installation. THIS booster seat is NOT CERTIFIED for use in an aircraft.**



### Step 2 – † Required Options: **ALWAYS** use LATCH connectors (419) **or** stabilizing bow (421) for stability purposes!

	HCPCS	Unit	Item no.	Price
†Lower anchors connection	E1399	Each	£ 419	\$178
†Stabilizing bow (if car has no latch)	E1399	Each	£ 421	\$132



419



421

#### Accessories for additional stability purposes

Front stand (min. 15.75" footwell height)	E1399	Each	£ 422	\$415
Cushion, leveling (for shaped seats)	E1399	Each	£ 513	\$158



422



513

#### Accessories for i ddYf VcXm

<YUXfYgřž \Y][ \hUX1 gtable, w. bracket	E1399	Each	£ 404	\$326
Shoulder guides, left & right, w. brackets	E1399	Pair	£ 427	\$300
Upper belt guide, additional*(see pic below)	E1399	Each	£ 425	\$84
Headrest, adjustable wings	E1399	Each	£ 436	\$431



404 <YUXfYgřž \Y][ \hUX1 gtable



427 Shoulder guides



425 Add. upper belt guide\*



436 <YUXfYgřž aX1 gtable wings

Technical Specs Hercules Prime	
Seat depth	40 - 46 cm / 15.7 - 18.1"
Depth behind crotch strap	24 - 38 cm / 9.5 - 15"
Seat width	40 cm / 15.7"
Seat to head	63 - 76 cm / 24.8 - 29.9"
Back width (seat)	45 cm / 17.7"
Shoulder height	38 - 60 cm / 15 - 23.6"
Shoulder width with shoulder guide	42 - 51 cm / 16.5 - 20.1"
Back angle (max.)	120°
Seat tilt-in-space (with tilting base 10°)	0° / 2.5° / 5° / 7.5° / 10°
Lower leg length (with foot rest)	17 - 43 cm / 6.7 - 16.9"
Height of lateral trunk supports (upper edge)	26 - 36 cm / 10.2 - 14.2"
Distance lateral trunk supports (max.)	35 cm / 13.8"
Max. overall dimensions (w x h x d)	52 x 88 x 43 cm / 20.5 x 34.6 x 16.9"
Min. folded dimensions (w x h x d)	52 x 21 x 73 cm / 20.5 x 8.3 x 28.7"
Weight of seat	13.4 kg / 29.5 lbs
Body height user	130 - 180 cm / 51.1 - 70.9"
User weight	18.1 - 80 kg / 40 - 176.4 lbs

## Hercules Prime –2026 US Retail prices

info@thomashilfen.com | Phone 866.870.2122 | www.Thomashilfen.us

### Accessories for trunk support

Lateral trunk supports, swing-away, w. brackets*	E1399	Pair	£ 405	\$835
Lateral trunk supports, rigid, w. brackets*	E1399	Pair	£ 429	\$304
Upper body guides, swing-away, w. brackets*	E1399	Pair	£ 406	\$909
Upper body guides, rigid, w. brackets*	E1399	Pair	£ 428	\$314

\*Upper body guides cannot be used together with lateral trunk supports



405 Lateral trunk supports, swing-away



429 Lateral trunk supports, rigid



406 Upper body guides, swing-away



428 Upper body guides, rigid

### Accessories for mid and lower body

Abduction block	E1399	Each	£ 412	\$293
Leg guides, left + right, w. brackets	E1399	Pair	£ 415	\$283
Footrest, foldable, angle adjust. (min. 15.75" footwell height)	E1399	Each	£ 417	\$588
Incontinence seat pad	E1399	Each	£ 497	\$94
Front cover	E1399	Each	£ 433	\$58
Crotch pad with magnetic lock	E1399	Each	£ 434	\$247
Extra key for magnetic lock	E1399	Each	£ E66331	\$25



412 Abduction block



415 Leg guides, left + right, w. brackets




417 Footrest, foldable, angle adjustable



497 Incontinence seat pad



433 Front cover



434 Crotch pad with magnetic lock

\*Swivel direction **CANNOT** be changed!

#### \*\*CAUTION

**ALWAYS** keep clip on chest height level.

Hercules Prime **AND** accessories exceed US FMVSS 213 standards

Specs & pricing are subject to change without notice. HCPCS are suggestions only with no guarantee of reimbursement.  
ExoMotion, LLC dba Thomashilfen | www.thomashilfen.us Effective Date January 1st 2026