

Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	

Step 1 – Seat & Frame (Required options are indicated below by † symbol)



£ 6808 Swiftly 2	HCPCS	Description
†Seat 680821 \$2967 E1236		Folding frame with integrated/ contoured seat, fixed tilt, adjustable seat depth from underneath, footrest height adjustable, back recline, suitable for Early Intervention use
†Back 680822 \$132 E1014		
†Castor locks 6810 \$101 K0108		Castor pin locks, pair
†Bus tie-down 6830 \$221 K0108		Bus tie-down loop kit *
†5-pt. harness 6629 \$153 E0960		5-point harness
Seat Size Chart		
	Lower leg length	Seat depth
Swiftly 2	8.3" – 14.6"	14.8"
	Seat width	Back height
	15.6"	27.6"
	Min. shoulder height	Hip angle
	11.8"	90° to 100°

Step 2 – Positioning & General Options

Lower Body Positioning	HCPCS	Unit	Price
Foot straps with cam buckle	E0952	Pair	£ 6648 \$37
Ankle Huggers**	K0108	Pair	£ 213151 \$221
Footrest pad, height 1¾"	K0108	Each	£ 10023 \$104

Mid Body Positioning	HCPCS	Unit	Price
Abduction block, padded	E0957	Each	£ 6720-6 \$179
Pelvic positioning supports	2-point, padded, SR buckle	E0978	Each £ 8470 \$116
	4-point, padded, SR buckle	E0978	Each £ 8471 \$158
	Soft pelvic harness, medium	E0978	Each £ 6714 \$221
	Leg harness, large**	E0978	Pair £ 213152 \$221
Crotch strap (fits 8470)	E0978	Each	£ 6341 \$37
Tray, height adjustable, clear with picture window	E0950	Each	£ 67281 \$237
Mounting hardware tray †required ordering 6728	E1034	Each	£ 67282 \$258

Upper Body Positioning	HCPCS	Unit	Price
Chest harness supports	Butterfly style, neoprene	E0960	Each £ 6585 \$221
	Vest style	E0960	Each £ 6632 \$221
	H-style, posture support, small**	E0960	Each £ 214806 \$221
	PivotFit shoulder harness, small**	E0960	Each £ 214270 \$242
	PivotFit shoulder harness, medium**	E0960	Each £ 214250 \$242
	5-point harness	E0960	Each £ 8480 \$153
Trunk supports	Rigid lateral	E0956	Pair £ 6734 \$283
	Flexible lateral with chest strap	E0956	Pair £ 6633 \$310
Head supports	Contoured pillow, black	E0955	Each £ 6837 \$221
	Contoured pillow, red	E0955	Each £ 68371 \$221

General Accessories	HCPCS	Unit	Price
Crotch pad with magnetic lock	E1399	Each	£ 434 \$247
Extra key for magnetic crotch pad	E1399	Each	£ E66331 \$25
Seat minimizer kit for Early Intervention use	K0108	Set	£ 6836 \$305
Seat width reducer, black (seat width 11.6")	K0108	Pair	£ 6841-0 \$147
Seat width reducer, red (seat width 11.6")	K0108	Pair	£ 6841-1 \$147
Pneumatic wheels, full set, black	K0070	Set	£ 6671 \$325
Under seat storage basket	K0108	Each	£ 6835 \$79
Frame padding, black	K0108	Pair	£ 6622-6 \$132
Frame padding, red/black, reversible	K0108	Pair	£ 6622-7 \$132
Grab bar, padded, black	K0108	Each	£ 6721-5 \$168
Tip assist	K0108	Each	£ 6829 \$111
Straps	Secure folded seat/frame	K0108	Each £ 667011 \$37
	Wrist "runaway" tether	K0108	Each £ 667111 \$37
Weather protection	Canopy, sun/rain with window, black	K0108	Each £ 663761 \$342
	Lower body windbreaker for 663761	K0108	Each £ 661852 \$104
	Canopy, rain only, clear	K0108	Each £ 6638 \$258
	Rain shield, clear for canopy 663761	K0108	Each £ 6739 \$132
	Rain cape, grey	K0108	Each £ 6540 \$172
	Sleeping bag, summer weight, black	K0108	Each £ 6719 \$314
	Sleeping bag, winter weight, black	K0108	Each £ 6718 \$314
Umbrella, angle adjustable, grey	K0108	Each	£ 6594 \$100

Swiftly 2	Technical Specs #
Capacity	110 lbs
Unit weight	41 lbs
Seat tilt	20°
Back recline	90-100°
Upholstery	Trevira-CS
Floor to handle	32.3"-48.8"
Floor to seat	22.8"
Max width	26.4"
Fold length	31.1"
Fold height	19.3"
Brakes	foot operated
Seat direction	front facing
Front wheels	7.1" solid
Rear wheels	9.8" solid



* Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

** Bodypoint®