





Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	

**Step 1 – Seat & Frame (Required options are indicated below by † symbol)**

	£ 6809 Swifty 3	HCPCS	Description
†Seat	680901 \$3471	E1236	Folding frame with integrated/ contoured seat, fixed tilt, adjustable seat depth from underneath, footrest height adjustable, reclining back rest
†Back	680902 \$150	E1014	
†Tie-down	680903 \$221	K0108	Bus tie-down loop kit *
†Castor locks	680904 \$101	K0073	Castor pin locks, pair
†Tip assist	680905 \$111	E0971	Tip assist, anti-curb tipper
†5-pt. harness	680906 \$230	E0960	5-point positioning harness
†Foot plate	680907 \$105	K0040	Foot plate, flip up
†Leg rest	680908 \$210	E0990	Legrest, elevating

**Step 2 – Positioning & General Options**

Lower Body Positioning	HCPCS	Unit	Price
Foot straps with cam buckle	E0952	Pair	£ 6648 \$37
Calf strap	E0978	Each	£ coming soon



6648 Footstraps w. cam buckle

**Mid Body Positioning**

Abduction block, adjustable	E0957	Each	£ 412	\$293
2-point pelvic belt, padded	E0978	Each	£ coming soon	



412 Abduction block, adjustable

**Upper Body Positioning Lateral**

Head supports, adjustable	E0955	Pair	£ 6955	\$263
Rigid lateral, left & right, height adj.	E0960	Pair	£ 6956	\$283
Swing away lateral, l&r, height adj.	E0960	Pair	£ 6957	\$462
Chin support	K0108	Each	£ 562	\$321



6955 Head supports, adj.



6956 Rigid lateral, left & right



6957 Swing away lateral, l & r



562 Chin support

Swifty 3	Technical Specs
Thigh length / seat depth (with seat minimizer)	37 - 46 cm / 14.6 - 18.1" (31) / (12.2)
Seat width (with seat minimizer)	41 cm (33) / 16.1" (13)
Vertex height (back height) (min. - max.)	68 - 83 cm / 26.8 - 32.7"
Back height 1 (height to armpit) (min. - max.)	26 - 53 cm / 10.2 - 20.9"
Back height 2 (shoulder height) (min. - max.)	49 - 65 cm / 19.3 - 25.6"
Lower leg length	25 - 45 cm / 9.8 - 17.7"
Back angle	90° - 115°
Knee angle	90° / 108° / 126° / 144° / 162° / 180°
Foot angle	90°
Seat tilt-in-space, fixed	7.2°
Size of footrest (W x D)	29.5 x 23 cm / 11.6 x 9.1"
Overall dimension (L x W x H) (settings min. / max.)	113 x 69 x 119 cm / 44.5 x 27.2 x 46.9" / 129 x 69 x 134 cm / 50.8 x 27.2 x 52.8"
Folded dimensions (L x W x H)	60 x 69 x 90 cm / 23.6 x 27.2 x 35.4" (standing)
Push bar height	79 - 116 cm / 31.1 - 45.7"
Front wheel, swiveling (size / width)	19 cm (7.5") / 5 cm (2")
Rear wheel, plug-in (size / width)	31 cm (12.2") / 5.5 cm (2.2")
Weight	24.5 kg / 54.0 lbs
Max. user size	175 cm / 68.9"
Max. user weight	82 kg / 180.8 lbs

\* Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness or 5-point harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

**General Accessories**


Crotch pad with magnetic lock	E1399	Each	£ 434	\$247
Extra key for crotch pad	E1399	Each	£ E66331	\$25
Seat minimizer kit (6959+6958)	E0190	Set	£ 6960	\$273
Seat depth reducer (seat depth 12.2")	E0190	Pair	£ 6959	\$80
Seat width reducer (seat width 13")	E0190	Pair	£ 6958	\$178
Under seat storage basket	K0108	Each	£ 6970	\$154
Under seat storage carry-on bag	K0108	Each	£ 6971	\$231
Canopy	K0108	Pair	£ 6963	\$339
Rain cape, poncho	K0108	Pair	£ 6964	\$250
Frame padding, black	K0108	Each	£ 6965	\$132
Incontinence pad	T4540	Each	£ 6961	\$94
Pneumatic wheels, full set	K0070	Set	£ 6966	\$337
Drum brake kit	E2228	Each	£ 6967	\$290



434 Crotch pad with magnetic lock



6959 Seat depth reducer 11.6"



6958 Seat width reducer 11.6"



6970 Under seat storage basket



6971 Under seat storage carry-on bag



6964 Rain cape, poncho



6963 Canopy



6965 Frame padding, black



6961 Incontinence pad



6966 Pneumatic wheels, full set



6967 Drum brake kit

\*Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt, 5-point and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.  
 Specs & pricing are subject to change without notice. HCPCS are suggestions only with no guarantee of reimbursement.  
 ExoMotion LLC dba Thomashilfen | [www.thomashilfen.us](http://www.thomashilfen.us) Effective Date June 1<sup>st</sup> 2026