



Yes, you can.®

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Invacare® IVC™ 9000 XT Wheelchair

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

Existing information (serial number) will need to be supplied to Customer Service to ensure correct configuration.

PRE-CONFIGURED STOCK WHEELCHAIRS

STOCK CHAIRS

Quantity	Model Number	Description	List Price
<input type="checkbox"/>	9XT_PTO_29154	9XT/WD66/88 9000XT 16" x 16" Adult, Desk Length Fixed Ht Space Saver Arm (9153629154)	\$1,070.00
<input type="checkbox"/>	9XT_PTO_29153	9XT/WD86/88 9000XT 18" x 16" Adult, Desk Length Fixed Ht Space Saver Arm (9153629153)	\$1,070.00
<input type="checkbox"/>	9XT_PTO_29155	9XT/WD86/78 9000XT 18" x 16" Adult, Full Length Fixed Ht Space Saver Arm (9153629155)	\$1,070.00
<input type="checkbox"/>	9XT_PTO_33131	9XT/WD06/88 9000XT 20" x 16" Adult, Desk Length Fixed Ht Space Saver Arm (9153633131).....	\$1,212.00

STOCK FRONT RIGGINGS

Quantity	Model Number	Description	List Price
<input type="checkbox"/>	T93HA_PTO_37801	T93HA Hemi Sw. Footrests Aluminum Footplates with Heel Loops (9153637801)	\$93.00 E0951
<input type="checkbox"/>	T93HC_PTO_37786	T93HC Hemi Sw. Footrests Composite Footplates with Heel Loops (9153637786)	\$93.00 E0951
<input type="checkbox"/>	T94HA_PTO_37804	T94HA/U67 Hemi Sw. Elev Legrests Aluminum Footplates (9153637804)	\$180.00 K0195 E0990
<input type="checkbox"/>	T94HC_PTO_37788	T94HC/U67 Hemi Sw. Elev Legrests Composite Footplates (9153637788)	\$180.00 K0195 E0990
<input type="checkbox"/>	AHL4A_PTO_19911	AHL4A/U67 Hemi Smartleg Art. Legrests Aluminum (9153619911)	\$320.00 K0053

STOCK OPTIONS

Quantity	Model Number	Description	List Price
<input type="checkbox"/>	9153632671	1497 Oxygen Holder Pkg Assy	\$151.00 E2208
<input type="checkbox"/>	9153632659	1496 IV / Oxygen Hldr Pkg Assy	\$158.00 K0108
<input type="checkbox"/>	9153629710	1360 Adj Rear Anti - Tip 8" Caster	\$82.00 E0971
<input type="checkbox"/>	9153629726	1360A Adj Rear Anti - Tip 6" Caster	\$82.00 E0971
<input type="checkbox"/>	9153629713	A1576 Wheel Lock Extensions	\$39.00 E0961

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® IVC™ 9000 XT Wheelchair

FRAME AND SEAT SELECTIONS

FRAME TYPE

- 9XT 9000XT \$1,070.00 K0004

FRAME WIDTH AND DEPTH

- WD46 14" x 16" No Charge A9900
 WD48 14" x 18" \$118.00 A9900
 WD56 15" x 16" \$118.00 A9900
 WD58 15" x 18" \$118.00 A9900
 WD66 16" x 16" No Charge
 WD67 16" x 17" \$118.00 A9900
 WD68 16" x 18" \$118.00 A9900
 WD76 17" x 16" \$118.00 A9900
 WD77 17" x 17" \$118.00 A9900
 WD78 17" x 18" \$118.00 A9900
 WD86 18" x 16" No Charge
 WD87 18" x 17" \$118.00 A9900
 WD88 18" x 18" \$118.00 A9900
 WD96 19" x 16" \$142.00 A9900
Not available with Super Low seat-to-floor.
 WD97 19" x 17" \$254.00 A9900
Not available with Super Low seat-to-floor.
 WD98 19" x 18" \$254.00 A9900
Not available with Super Low seat-to-floor.
 WD06 20" x 16" \$142.00 E2201
Not available with Super Low seat-to-floor.
 WD07 20" x 17" \$254.00 E2201
Not available with Super Low seat-to-floor.
 WD08 20" x 18" \$254.00 E2201
Not available with Super Low seat-to-floor.
 WD26 22" x 16" \$142.00 E2201
Not available with Super Low seat-to-floor.
 WD27 22" x 17" \$254.00 E2201
Not available with Super Low seat-to-floor.
 WD28 22" x 18" \$254.00 E2201
Not available with Super Low seat-to-floor.

SEAT-TO-FLOOR HEIGHT

- A Adult (19 1/2" STF)..... No Charge
Based on 24" wheels and 8" casters. See seat-to-floor matrix.
 H Hemi (17 1/2" STF) No Charge
Based on 24" wheels and 8" casters. See seat-to-floor matrix.
 SH Super Low (14 1/2" STF) \$85.00
 SHE Super Hemi (15 1/2" STF) No Charge
Based on 20" wheels and 6" casters. See seat-to-floor matrix. Requires AT902 (6" casters), U1763BK (20" composite wheels, urethane tires) AHDRM (aluminum handrims) and A1576 (wheel lock extensions). Seat-to-floor height is approximately 15 1/2"

ARM AND BACK OPTIONS

ARM TYPE

Arm types are removable, flip-back.

- 18 Full Length Fixed Ht Conv (10") No Charge A9900
 19 Full Length Adj. Ht. Conv (10" - 14") \$126.00 E0973
 28 Desk Length Fixed Ht. Conv (10") No Charge A9900
 29 Desk Length Adj. Ht. Conv (10" - 14") \$126.00 E0973
 78 Full Length Fixed Ht. SS (10") No Charge
Not available with adjustable angle back type.
 79 Full Length Adj. Ht. SS (10" - 14") \$126.00 E0973
Not available with adjustable angle back type.
 88 Desk Length Fixed Ht. SS (10") No Charge
Not available with adjustable angle back type.
 89 Desk Length Adj. Ht. SS (10" - 14") \$126.00 E0973
Not available with adjustable angle back type.
 79JY Full Length Adj. Ht. SS (8" - 10") \$130.00 E0973
Not available with adjustable angle back type.
 89JY Desk Length Adj. Ht. SS (8" - 10") \$130.00 E0973
Not available with adjustable angle back type.

BACK TYPE

- ADJ Adj. Back with 10° Bend No Charge
 S1418 Straight Back w/ Push Handles (15" - 18") ..No Charge K0108
 ANGLE Adj. Angle Back with 10° Bend \$100.00 K0108
Not available with Super Low seat-to-floor. Standard with rear anti-tippers.
 ANGLE-S Adj. Angle Back with Straight \$100.00
Not available with Super Low seat-to-floor. Standard with rear anti-tippers.

COLOR OPTIONS

FRAME FINISH

- 08P Deep Purple No Charge
 17P Cobalt Blue No Charge
 24P Wet Black No Charge
 34P Satin (Flat) Black No Charge
 71P Silver Vein No Charge
 88P Wet Black with Twilight Sparkle No Charge
 105P Midnight Blue No Charge

CANDY COLORS

- 61P Electric Red No Charge
 62P Electric Blue No Charge
 66P Electric Purple No Charge
 104P Electric Teal No Charge

UPHOLSTERY COLORS

- U67 Black Vinyl No Charge
Vinyl upholstery must choose WRAP or SOUPH.
 U240 Black Nylon Upholstery Standard
 U245 Black Sail Cloth Adj. Tension Seat \$10.00 A9900
U245 sail cloth adjustable tension seat upholstery available with U240 black nylon back upholstery only.

NAUGAHYDE

- WRAP Wrap Around Upholstery \$77.00 A9900
 SOUPH Screwed On Back Upholstery \$89.00 A9900
Not available with Adjustable Angle or Adjustable Angle-S Back Type.

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® IVC™ 9000 XT Wheelchair

CASTERS

CASTER OPTIONS

- AT902 6" Composite Urethane 3 SpokeNo Charge A9900
- 1228 8" x 1 3/4" Semi. Pneu Casters Lt Grey \$56.00
Not available with Super Low seat-to-floor.
- AT903 8" x 1 1/4" Composite Urethane Standard
Not available with Super Low seat-to-floor.
- 1299 8" x 1" Solid Rubber CastersNo Charge
Not available with Super Low seat-to-floor.
- 1292 8" x 1 1/4" Pneu Casters \$87.00 A9900
Not available with Super Low seat-to-floor.
- 1492 8" x 1 1/4" Pneu. with Flat Free Inserts \$104.00 E2217
Not available with Super Low seat-to-floor.

CASTER JOURNAL OPTION

Not available with Adjustable Angle or Adjustable Angle-S Back Type.
Not available with Super Low seat-to-floor.
Fixed caster journals are standard.

- ADJC Adjustable Caster Journal \$29.00 A9900
- ADJCPB Adjustable Precision Bearing Caster Journal \$41.00 A9900

WHEELS

WHEEL OPTIONS

- 1763BK 20" Composite Pneu. Tires \$117.00 A9900
- 4763BK 20" Composite Pneu. Flat Free Inserts \$232.00 E2213
- U1763BK 20" Composite Urethane Tires \$113.00 A9900
- 1764BK 22" Composite Pneu. Tires \$117.00 A9900
- 4764BK 22" Composite Pneu Flat Free Inserts \$232.00 E2213
- 5244BK 24" HD Comp. Pneu. with FF Inserts \$298.00 E2213
Not available with Super Low seat-to-floor.
- U1764BK 22" Composite Urethane Tires \$113.00 A9900
- 2244BK 24" Composite Pneu. Tires \$140.00 A9900
Not available with Super Low seat-to-floor.
- 4244BK 24" Composite Pneu. with FF Inserts \$225.00 E2213
Not available with Super Low seat-to-floor.
- U2244BK 24" Composite Urethane Treaded Tires \$120.00 A9900
Not available with Super Low seat-to-floor.
- U2222C 24" Composite Urethane Corded Tires Standard
Not available with Super Low seat-to-floor.

WHEEL OPTIONS (ONE ARM DRIVE - OAD)

The Cyclical Lever drive requires treaded rear wheels and is only available on 16", 17", and 18" deep. Hemi seat-to-floor must use space saver arms. Weight limit is 250 lbs. The above options are not available with amputee attachment or offset axle.

- 99457R RH OAD with 4244BK / Flat Free \$697.00 E0958
- 99458L LH OAD with 4244BK / Flat Free \$697.00 E0958
- 9456L LH One Arm Drive with U2222C \$606.00 E0958
- 9456R RH One Arm Drive with U2222C \$606.00 E0958
- LDR96R1 6" Right Hand Cyclical Lever One Arm Drive .. \$797.00 E0958
Not available with Super Low seat-to-floor.
Not available with adjustable caster journal.
- LDR96L1 6" Left Hand Cyclical Lever One Arm Drive \$797.00 E0958
Not available with Super Low seat-to-floor.
Not available with adjustable caster journal.
- LDR98R1 8" Right Hand Cyclical Lever One Arm Drive .. \$797.00 E0958
Not available with Super Low seat-to-floor.
Not available with adjustable caster journal.
- LDR98L1 8" Left Hand Cyclical Lever One Arm Drive \$797.00 E0958
Not available with Super Low seat-to-floor.
Not available with adjustable caster journal.

HANDRIM OPTIONS

- AHDRM Aluminum HandrimsNo Charge
Not available with U2222C rear wheels.
- CHDRM Chrome Plated HandrimsNo Charge
- COM Composite Handrims Standard
Not available on 20" and 22" composite rear wheels.
- 1512 Plstc Coated Aluminum Handrims \$85.00
- 1224-8 8 Alum Oblique Proj. Hdrm, Plstc Coat \$252.00 E0967
Not available on 20" and 22" composite rear wheels.
- 1225-8 8 Alum Vert. Proj. Hdrm, Plstc Coat \$252.00 E0967
Not available on 20" and 22" composite rear wheels.

REAR AXLE

- OFFAXL Offset Axle \$40.00 K0108
Not available with Adjustable Angle or Adjustable Angle-S Back Type. Not available with Super Low seat-to-floor.

WHEEL LOCKS

The above selections are not available with Adjustable Angle or Adjustable Angle-S Back Type.

- 9690B Pull To Lock Whl Locks (Std on Hemi) No Charge
- 9691B Push To Lock Whl Locks (Std On Adult) No Charge
- A1567 Pull To lock Clamp On Whl Locks \$76.00 A9900
- 1721 High Push To Lock Clamp On Whl Locks \$71.00 A9900

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® IVC™ 9000 XT Wheelchair

POSITIONING STRAP OPTIONS

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap \$50.00 E0978
- 1515 Push Button Style Seat Positioning Strap \$50.00 A9900 E0978
- 1321BK Hook and Loop Seat Positioning Strap \$40.00 E0978

CHEST POSITIONING STRAPS

- 9515 Push Button Style Chest Positioning Strap \$50.00 E0960

ACCESSORIES AND STOCK OPTIONS

BACK ACCESSORIES

- L8601 Omit Back Upholstery Less \$10.00
- ADJTENS Adjustable Tension Back Upholstery \$112.00
- 1160 Carrying Pocket \$54.00 A9270
- STC Stencil Back Upholstery \$29.00 A9270
Only available with vinyl upholstery. Block letters are .75" in height.
- RB Dual Back Liner \$49.00 A9270
- BP02 Back Pack \$50.00 A9270

SEAT ACCESSORIES

- L8651 Omit Seat Upholstery Less \$10.00
- 9120 Solid Seat Insert \$98.00 A9900
Not available on 22" wide chairs.
- RS Dual Seat Liner \$49.00 A9270

ARM ACCESSORIES

- NRAS Non-Removable Arm With Screws \$10.00 A9900
- 1300 Utility Tray \$157.00 E0950
Not available with 24" rear wheel set in Hemi position.

MISCELLANEOUS ACCESSORIES

- 1305L Crutch and Cane Carrier \$43.00 E2207
Not available with traditional one arm drive, but is available with CLD one arm drive.
- 1577QR Quick Release Axles \$86.00 K0108
**Not available with offset axle.
Not available with any one arm drive.**
- 2218 IV Rod, Telescopic \$150.00 K0105
Not available with traditional one arm drive, but is available with CLD one arm drive.
- 2226 Front Anti - Tippers \$83.00 E0971
- 9357SH Super Low Amputee Attachment \$95.00 E0959
Only available on Super Low Seat-to-Floor.
- 1360SH Super Low Adj Rear Anti - Tipper \$82.00 E0971
Only available on Super Low Seat-to-Floor.
- 9357L Amputee Attachment \$95.00 E0959
**Not available with offset axle.
Not available with any one arm drive.**

STOCK OPTIONS

- 1497 Oxygen Holder Package Assembly \$151.00 E2208
**Not available with traditional one arm drive, but is available with CLD one arm drive.
Not available with offset axle.
For 0° and 6° seat angles. See service manual for details.**
- 1496 IV / Oxygen Holder Package Assembly \$158.00 K0108
**Not available with offset axle.
Not available with any one arm drive.**
- 1360 Adjustable Rear Anti - Tipper (8" Casters) \$82.00 E0971
**Only available with 8" casters.
Standard with Adjustable Angle Back.**
- 1360A Adjustable Rear Anti - Tipper (6" Casters) \$82.00 E0971
Only available with 6" casters.
- A1576 Wheel Lock Extensions \$39.00 E0961
Elastic cord wheel lock extensions. Hemi seat to floor come standard with non-cord wheel lock extensions.
- 1100 Tool Kit No Charge

FRONT RIGGING OPTIONS

FRONT RIGGING

- T93HA Hemi Sw. Footrests Aluminum Footplates with Heel Loops .. \$93.00
- T93HC Hemi Sw. Footrests Composite Footplates with Heel Loops .. \$93.00
- T94HA Hemi Sw. Elev. Legrests Aluminum Footplates \$180.00 K0195 E0990
- T94HC Hemi Sw. Elev. Legrests Composite Footplates .. \$180.00 K0195 E0990
- SH904M Super Low Elev. Legrests 11" - 17" \$266.00 K0195
- SH904MS Super Low Elev. Legrests 14" Wide \$266.00 K0195
- SH93E Super Low Ext. Sw. Footrests 11" - 17" \$185.00 A9900
- SH93ES Super Low Ext. Sw. Footrests 14" Wide \$185.00 A9900
- H93M Hemi Minimum Sw. Footrests \$185.00 A9900
- H93E Hemi Ext Sw. Footrests \$185.00 A9900
- AHL4A Hemi Smartleg Art. Legrests 16" and Above \$320.00 K0053
- AHL4C Hemi Smartleg Art. Legrests 16" and Above \$320.00 K0053
- AHS4A Hemi Smartleg Art. Legrests 14" - 15" Wide ... \$320.00 K0053
- T93HS Hemi Sw. Footrests Aluminum Footplates 14" Wide \$90.00 A9900
- T94HS Hemi Sw. Elev. Legrests Aluminum Footplates 14" Wide .. \$175.00 K0195
- T963A 60° Footrests Aluminum Footplate \$90.00 A9900
- T963C 60° Footrests Composite Footplate \$90.00 A9900
- T963S 60° Footrests Composite Footplate 14" Wide \$90.00 A9900

FOOTPLATE OPTIONS

- 1350 Extra Large Aluminum Footplates \$79.00 K0041
- AT5543 Adjustable Angle Flip - up Footplates \$67.00 K0040

FOOTREST / LEGREST ACCESSORIES


- 1820 Non - Removable Device \$12.00 A9900
- 1337 Calf Strap \$36.00 K0038
- 1346 Heel Loops \$39.00 E0951
- IMPACTG Impact Guards \$51.00

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® IVC™ 9000 XT Wheelchair

Seat-to-Floor 6" Casters			
Seat-to-Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
15 1/2"	Top	20"	Top
16 1/2"	Middle	22"	Top
17 1/2"	Bottom	24"	Top
17 1/2"	Bottom	20"	Bottom

Seat-to-Floor 8" Casters			
Seat-to-Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
17 1/2"	Top	20"	Bottom
17 1/2"	Top	24"	Top
18 1/2"	Middle	22"	Bottom
19 1/2"	Bottom	24"	Bottom



Weight Limitation 250 Pounds
 Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

Seat-to-Floor 6" Casters and Offset Axle Bracket					
Seat-to-Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position	Axle Plate Mounting on Frame	
				Top Lug	Bottom Lug
15 1/2"	Top Hole	22"	#5 Hole	#4 Hole	#1 Hole
15 1/2"	Top Hole	20"	#3 Hole	#4 Hole	#1 Hole
16 1/2"	Middle Hole	24"	#5 Hole	#4 Hole	#1 Hole
16 1/2"	Middle Hole	22"	#3 Hole	#4 Hole	#1 Hole
16 1/2"	Middle Hole	20"	#2 Hole	#4 Hole	#1 Hole
17 1/2"	Bottom Hole	24"	#3 Hole	#4 Hole	#1 Hole
17 1/2"	Bottom Hole	22"	#2 Hole	#4 Hole	#1 Hole
17 1/2"	Bottom Hole	20"	#3 Hole	#4 Hole	#1 Hole

For a complete list of seat-to-floor heights and angles that are configurable by the dealer refer to the service manual.

Seat-to-Floor 8" Casters and Offset Axle Bracket					
Seat-to-Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position	Axle Plate Mounting on Frame	
				Top Lug	Bottom Lug
17 1/2"	Top Hole	20"	#3 Hole	#1 Hole	#4 Hole
17 1/2"	Top Hole	22"	#2 Hole	#4 Hole	#1 Hole
17 1/2"	Top Hole	24"	#3 Hole	#4 Hole	#1 Hole
18 1/2"	Middle Hole	22"	#3 Hole	#1 Hole	#4 Hole
18 1/2"	Middle Hole	24"	#2 Hole	#4 Hole	#1 Hole
19 1/2"	Bottom Hole	20"	#3 Hole	#1 Hole	#4 Hole

For a complete list of seat-to-floor heights and angles that are configurable by the dealer refer to the service manual.

Comments

NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.



Invacare Corporation
 USA
 One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

www.invacare.com

©2012 Invacare Corporation
 All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.