

The newly redesigned Outlander Full Platform Exterior Lift features an on/off key switch, a toggle switch for up and down, a puddle light and an integrated license plate bracket. Easy drive-on loading from either end of the lift takes the hassle out of loading your scooter or power chair in a busy parking area. Simple hitched-based installation requires little to no modification to your vehicle and allows for easy removal.

FEATURES AND BENEFITS

- Puddle light for ease of use at night
- Toggle switch for ease of up and down use
- Drive on/drive off
- Integrated license plate bracket
- Compatible with Class II, III, or IV hitches
- Manual override
- LED indicator for On/Off
- Phosphate undercoating for weather resistance

Outlander

BASE WEIGHT	81.4 lbs.
MAXIMUM SCOOTER & PC OVERALL WEIGHT	350 lbs.
MAXIMUM PRODUCT WIDTH	28"
OVERALL PLATFORM LENGTH	42"
OVERALL PLATFORM WIDTH	30"
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	20"
MOTOR	Fully sealed 12V DC
WARRANTY	3-Year transferable
HITCH CLASS	II, III, IV
OPTIONS	Swing-away mount, Scooter lockdown arm, Power chair hold-down bar, Strap kit
VEHICLE PLATFORM	Full-size van, Minivan, Truck, Sedan, SUV and Crossover
LIFTS CURRENT PRIDE® PRODUCTS	Accepts all Pride® Scooters and Power Chairs except the Pursuit® and Pursuit® XL

