

Hill-Rom®

P310 Wound Surface Prevention and Treatment Surface



At a glance

- Multi-zoned, alternating pressure with low air loss
- Designed for low to moderate-risk patients to prevent and treat Stage I through uncomplicated Stage IV pressure ulcers
- Accommodates patients weighing up to 400 lbs
- Unique cell-on-cell technology provides added support and minimizes bottoming out
- Mattress length options of 80" and 84"
- Standard quick inflate device for faster mattress inflation
- nano Ag+ Technology that protects the surface against stain and odor-causing bacteria is standard*
- Supports Hill-Rom® 1000 bed, CareAssist® ES bed, Advanta™ 2 bed and other flat deck frames

A cost-effective option for the prevention and treatment of pressure ulcers

The Hill-Rom® P310 Wound Surface is a compact, lightweight prevention and treatment surface designed with the value-conscious customer in mind. This surface combines both alternating pressure and low air loss modalities in one mattress, resulting in an easy to use, affordable option with added safety components that both patients and caregivers will appreciate.

In addition to its performance features, the surface has also been designed to fit the CareAssist® ES bed and other flat deck frames to maximize the level of therapy and safety to the patient and caregiver.

* This product is not intended to protect users or others against bacteria, viruses, or other disease organisms. Always disinfect Hill-Rom® surfaces thoroughly after each use, according to facility protocol.

Hill-Rom®

Enhancing Outcomes for Patients and Their Caregivers.™

Hill-Rom®

P310 Wound Surface Prevention and Treatment Surface

Standard antimicrobial ticking
nano Ag+ technology

Bi-directional stretch ticking
relieves shear



Alternating pressure bladders
with low air loss



Lightweight control box
Choice of four alternating pressure
cycles with ten comfort control settings



Optional bolsters
Integrated cradle design



Anti-skid bottom
Added resistance to movement

Technical Specification

Surface and Controller – 80" version

P310MAT

Internal Dimension (LxWxH).....80" x 35.4" x 8"
Weight.....16.5 lbs

Surface and Controller (with Bolsters) – 80" version

P310MAT

Internal Dimension (LxWxH).....80" x 35.4" x 8"/10"
Weight.....17.5 lbs
Meets flammability requirements: 16CFR1632, 16CFR1633,
CAL129, CAL603

Surface and Controller – 84" version

P310MAT

Internal Dimension (LxWxH).....84" x 35.4" x 8"
Weight.....16.5 lbs

Surface and Controller (with Bolsters) – 84" version

P310MAT

Internal Dimension (LxWxH).....84" x 35.4" x 8"/10"
Weight.....17.5 lbs
Meets flammability requirements: 16CFR1632, 16CFR1633,
CAL129, CAL603

Control Unit

P310ACU

Power Supply .. AC 120V 50/60 Hz, 0.17A (for 120V system)
Fuse Rating..... T1A, 250V
Cycle Time..... Four settings (10,15,20 minute and OptiRest)
Dimension (LxWxH) 11.4" x7.9" x4.6"
Weight.....5 lbs

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2010 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

162239 rev 1 1/22/10



Enhancing Outcomes for Patients and Their Caregivers.™

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com