




PRECISION
COMFORT

GEL



Features and Benefits

- High-density, high-resilience foam for long lasting support
- Coccyx and sacral relief for longer sitting comfort and skin protection at these vulnerable points
- Super soft silicone gel insert in overlay which forms to the IT and coccyx for pressure reduction and prevents shearing
- Gel-infused viscoelastic upper foam in the overlay to help keep the skin cool and comfortable
- HD high-resilience base foam around the gel insert in the overlay to help support the greater trochanter and femurs for enhanced pressure reduction
- Removable gel overlay designed to stay in place to allow the clinician to modify the base for a custom fit
- Ergonomic contour which helps maintain proper posture and reduces fatigue
- Moisture wicking, four-way stretch, and fluid resistant cover which keeps the person dry and comfortable while allowing the user to form into the layers of visco foam and gel properly



- Custom sizes
- Replaceable and billable base and overlay
- 2 year warranty

PC Size Chart												
WIDTH												
DEPTH		14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"
	14"	SD	CS	SD	CS	CS	CS	CS	CS	CS	CS	CS
	15"	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
	16"	SD	CS	SD	CS	SD	CS	SD	CS	CS	CS	CS
	17"	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
	18"	SD	CS	SD	CS	SD	CS	SD	CS	SD	CS	SD
	19"	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
	20"	CS	CS	SD	CS	SD	CS	SD	CS	SD	CS	SD
	21"	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
	22"	CS	CS	CS	CS	CS	CS	CS	CS	SD	CS	SD
	23"	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
	24"	CS	CS	CS	CS	CS	CS	CS	CS	SD	CS	SD

SD Standard Size CS Custom Size

HCPCS: E2603: 14"x14" - 20"x20"

HCPCS: E2604: 22"x18" - 24"x24"