

# BESTLIFT™

NEW!

## PL500HE | FULL BODY PATIENT LIFT

The BestLift™ PL500HE is the perfect blend of form and function for patient transfers. With an overall base length of 44.9", maximum width of 42.9", and base height of 4.5", it is the perfect solution for transferring bariatric patients in small rooms with low beds.



- **500lb** capacity
- Dual front castors
- New electronics with 24/7 capability
- New hand-operated base opening

# | PL500HE




## | DIMENSIONS

Key dimensions of PL500HE	
maximum capacity	500Lb / 230Kg
base legs length	44.9"
base closed external width	24.9"
based open external width	42.9"
mast height	48.3"
base height	4.5"
min. lifting height	25.8"
max. lifting height	69.1"
net weight	94Lb/43Kg
shipping weight	106Lb/48Kg


\* SPECIFICATIONS PROVIDED ARE A GENERAL GUIDELINE  
AND ARE SUBJECT TO CHANGE WITHOUT NOTICE



2920 Pacific Dr., Norcross, GA 30071  
www.bestcarellc.com  
678-679-6690



TOP VIEW OF PL500HE



SIDE VIEW OF PL500HE