

PR458-MED

Model Number	PR458-MED
Size	Medium
Weight Capacity	375 LBS
Recline Positions	3
Standard Fabrics	3
Fabric Options	All
Chaise Seat	YES
Coil Springs	NO
Back Type	Button
Head Rest Cover	YES
Arm Covers	YES
Battery Back Up	YES
Overall Width	32"
Overall Height	42"
Width Between Arms	21.5"
Floor to Top of Seat	20"
Seat Depth	21"
Seat to Top of Back	28"
Seat to Top of Arm	6.5"
Seat Angle Raised	15°
Floor to Top of Foot Rest Reclined	23"
Floor to Top of Back Reclined	33.5"
Back Angle Closed	107°
Back Angle Reclined	145°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	175°
Distance Required from Wall Reclined	17"
Extension reclined from Back of Seat to end of Foot Rest	42.5"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 Yrs.
Pro Rated Years 4 Thru 7	NO
Warranty on Mechanical Labor	3 Yrs.
Warranty on Electrical Labor	1 Yr.
Removable Back Option Available	YES
Smart Tek System Standard	YES
Customizable	NO
Left Hand Control Option Available	NO
Extra Pocket Option Available	NO
Heat & Massage Option Available	YES
Gel Memory Foam Option Available	NO
Power Pillow Option Available	NO
AutoFlex Option Available	NO
Foot Rest Extension Option Available	YES
Straight Lift Option Available	YES

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.



Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with
UL 73 and
CSA C22.2 No. 68

