



PRIDE[®]
M O B I L I T Y
L I V E Y O U R B E S T[®]

TRAVEL/SCOOTER COLLECTION



3
Wrangler®

4-5
Zero Turn 8¹ / Zero Turn 10¹

6-7
iRIDE™ /
Go-Go Elite Traveller® Plus¹

8-9
Go-Go Elite Traveller®¹ /
Go-Go® Sport¹

10-11
Go-Go® Folding Scooter¹ /
Go-Go® LX with CTS Suspension¹

12-13
Go-Go® Ultra X¹ / Go-Go® ES2¹

14-15
Revo® 2.0¹ / Victory® LX Sport¹

16-17
Raptor / Victory® 10.2¹

18-19
Victory® 10¹ / Victory® 9¹

20-21
Maxima¹ / Accessories

Up To
29.4
Miles Per Charge at 200 lbs.

94.75"
Turning Radius

Up To
11.4mph
Maximum Speed

350lbs.
Weight Capacity



Non-Medical Device

WRANGLER®

Feel the crisp, fresh air and enjoy the peaceful serenity of the outdoors. Forge a new path with the **Wrangler®**, an aggressive, outdoor scooter designed for extreme and rugged conditions. With powerful suspension, bright lighting and dual hydraulic brakes, this beast of a scooter is a powerful force.



Stone Grey

Desert Sand

ZERO TURN8¹

Your mobility has turned a corner with the Zero Turn 8, an innovative scooter equipped with *iTurn Technology™*, enabling you to maneuver this scooter through small spaces and around tight corners! Plus, enjoy full lighting, plenty of storage and a USB charger.



Up To
13.1
Miles Per Charge at 200 lbs.

37.25"
Turning Radius

Up To
5.9mph
Maximum Speed

325lbs.
Weight Capacity

FDA Class II Medical Device¹

Up To
24
Miles Per Charge at 200 lbs.

43"
Turning Radius

Up To
7.2mph
Maximum Speed

400lbs.
Weight Capacity



FDA Class II Medical Device¹

ZERO TURN10¹

The Zero Turn 10 can traverse the most beautiful, out-of-reach places and move seamlessly through doors and hallways in your home. This scooter combines powerful dual motors and two-wheel drive with revolutionary turning technology. Let the Zero Turn 10 carry you to places unknown.

Sapphire Blue

Pewter

Pearl White

Black Onyx

Green Moonstone

Fire Opal

Candy Apple Red

Ocean Blue
(Matte)

Pearl White

Black¹¹
(Matte)

Travel Mobility 5

iRIDE™

Journey far and wide with the iRIDE! This ultra-lightweight scooter is super portable and perfect for traveling. Equipped with both a manual hand brake and supplemental foot brake, the iRIDE delivers a safe and comfortable ride!



Up To
12.7
Miles Per Charge at 200 lbs.

33.5"
Turning Radius

Up To
3.7mph
Maximum Speed

250lbs.
Weight Capacity

Non-Medical Device

Up To
14.5
Miles Per Charge

45.75"
Turning Radius

Up To
4.5mph
Maximum Speed

300lbs.
Weight Capacity



Also available in 3-wheel. Please see back two pages for specs.

FDA Class II Medical Device¹

GO GO¹

ELITE TRAVELLER® PLUS

It's perfect for traveling, PLUS so much more! The Go-Go Elite Traveller® Plus offers easy disassembly and includes an ergonomic tiller for one-handed operation. Today's fast-paced world is no match for this stylish and portable scooter.

Raspberry

Mango

Arctic Ice

Lime

Blackberry



Red

Blue

Interchangeable Shrouds

Travel Mobility 7

GOGO[®] ELITE TRAVELLER[®]

Every day is a new adventure with the Go-Go Elite Traveller[®], an exceptionally stable and maneuverable travel scooter. With easy disassembly and a convenient storage basket, the Go-Go Elite Traveller[®] keeps up with your active life!



Up To
12
Miles Per Charge

44"
Turning Radius

Up To
4mph
Maximum Speed

300lbs.
Weight Capacity

FDA Class II Medical Device¹

Also available in 3-wheel. Please see back two pages for specs.

Up To
14.5
Miles Per Charge

50.5"
Turning Radius

Up To
4.7mph
Maximum Speed

325lbs.
Weight Capacity



Also available in 3-wheel. Please see back two pages for specs.

FDA Class II Medical Device¹

GOGO[®] SPORT

It's time to make your mark on this world! The Go-Go[®] Sport offers simple disassembly and an ergonomic delta tiller, perfect for one-handed operation. Enjoy the versatility of this sporty scooter, as well as bright LED lighting.

Red

Blue

Interchangeable Shrouds



Red

Blue

Interchangeable Shrouds

GOGO^{®1} FOLDING SCOOTER

With the Go-Go[®] Folding Scooter, the world is at your fingertips! With a unique folding design, this scooter folds in just three simple steps and fits easily in a closet, corner or vehicle trunk! Enjoy full lighting and speeds up to 4 mph.



Up To
13
Miles Per Charge

32"
Turning Radius

Up To
4mph
Maximum Speed

250lbs.
Weight Capacity

FDA Class II Medical Device¹

Up To
8.6
Miles Per Charge

44.75"
Turning Radius

Up To
4.4mph
Maximum Speed

300lbs.
Weight Capacity



Also available in 3-wheel. Please see back two pages for specs.

FDA Class II Medical Device¹

GOGO^{®1} LX WITH CTS SUSPENSION

There's smooth, and then there's the Go-Go[®] LX with CTS! This scooter has revolutionary suspension for the smoothest, most comfortable ride. From the park to your backyard, this feature-rich scooter offers best-in-class ride quality.

White



Red

Blue

Interchangeable Shrouds

Orion Yellow

Onyx Black

Pearl Pink

Sunburst Orange

Optional

Travel Mobility II

GOGO^{®1} ULTRAX

Enjoy worry-free travel with the Go-Go[®] Ultra X, an ultra-compact scooter that is simple to disassemble. Ride in comfort and style with the foldable swivel seat and keep your personal items handy with the front basket.



Up To
7.2
Miles Per Charge

44"
Turning Radius

Up To
4mph
Maximum Speed

260lbs.
Weight Capacity

FDA Class II Medical Device¹

Also available in 3-wheel. Please see back two pages for specs.

Red

GOGO^{®1} ES2

Get ready to go with the Go-Go[®] ES 2 Scooter! Compact and lightweight, this scooter has two speed settings for ultimate versatility, whether you are indoors or outdoors. It's perfect for errands and quick trips! Plus, a front basket provides plenty of storage space.



Up To
6.8
Miles Per Charge

31"
Turning Radius

Up To
4mph
Maximum Speed

250lbs.
Weight Capacity

FDA Class II Medical Device¹

Black



REVO® 2.0¹

Robustly rugged meets brilliant convenience. With the Revo® 2.0, there is nothing you can't handle! This tough and dependable scooter has powerful suspension for a smooth ride, yet easily disassembles for full portability.



Up To
17.8
Miles Per Charge at 200 lbs.

52"
Turning Radius

Up To
5mph
Maximum Speed

400lbs.
Weight Capacity

FDA Class II Medical Device¹

Also available in 3-wheel. Please see back two pages for specs.

Up To
19.7
Miles Per Charge at 200 lbs.

61.25"
Turning Radius

Up To
8mph
Maximum Speed

400lbs.
Weight Capacity



FDA Class II Medical Device¹

VICTORY® LX^{SPORT}¹

You're active and always on the go. You need a scooter that is up to the task! Appeal to your sporty side with the Victory® LX Sport. Enjoy premium suspension, bright lighting and other sophisticated features, all in one sporty package.

Grey Street

True Blue

Travel And
Scooter Mobility



Candy Apple Red

True Blue

Scooter Mobility 15





Do you have the need for speed? The **Raptor** drives at speeds up to 14 mph. Equipped with a regenerative braking system, ample storage and complete LED lighting, the Raptor appeals to the adventurer within you!



Up To
31
Miles Per Charge at 220 lbs.

60"
Turning Radius

Up To
14mph
Maximum Speed

400lbs.
Weight Capacity

Non-Medical Device

Matte Black with Optional Inserts

True Blue

Candy Apple Red

Green Machine

Lemon Crush

Pink Topaz

Evolution Orange

Matte Black

Optional Inserts



Up To
13.5
Miles Per Charge

52.25"
Turning Radius

Up To
5.2mph
Maximum Speed

400lbs.
Weight Capacity



Also available in 3-wheel. Please see back two pages for specs.

FDA Class II Medical Device¹

VICTORY® 10.2¹

Embrace your potential! The **Victory® 10.2** is fully loaded with style and advanced features, including frosted LED headlights and an integrated storage space with cup holder. Drive at speeds up to 5.2 mph and see what life has in store for you!

Candy Apple Red

Ocean Blue
(Matte)



Scooter Mobility 17

VICTORY 10¹

Versatility knows no bounds with the Victory® 10. Enjoy simple disassembly with one hand and drive up to 16 miles on a single charge. You won't miss a beat with this scooter, which includes a bright LED headlight. Anything is possible with the Victory® 10!



Up To
16
Miles Per Charge

54"
Turning Radius

Up To
5.3mph
Maximum Speed

400lbs.
Weight Capacity

FDA Class II Medical Device¹

Also available in 3-wheel. Please see back two pages for specs.

Candy Apple Red

Viper Blue



Up To
13
Miles Per Charge

51.5"
Turning Radius

Up To
5.3mph
Maximum Speed

300lbs.
Weight Capacity

Also available in 3-wheel. Please see back two pages for specs.

FDA Class II Medical Device¹

VICTORY 9¹

Sophisticated and stylish, the Victory® 9 is everything you could want in a scooter. The delta tiller lets you operate the scooter with a single hand. Enjoy easy disassembly and plenty of lighting with the LED headlight.



Mobilify Scooter 19

Candy Apple Red

MAXIMA¹

Heavy duty and robust, the **Maxima** scooter maximizes your potential! Drive at speeds up to 5.3 mph and experience a fully adjustable tiller. Enjoy the convenience of the front basket, perfect for storing your personal items.



Up To
18
Miles Per Charge

50.5"
Turning Radius

Up To
5.3mph
Maximum Speed


500lbs.
Weight Capacity

FDA Class II Medical Device¹

Also available in 3-wheel. Please see back two pages for specs.

Candy Apple Red

Viper Blue



1. Cane Clip Attachment | 2. Cup Holder | 3. Under Seat Storage | 4. Single OR Double Cane Holder | 5. XLR USB CHarger
6. Rear Basket | 7. Oxygen Holder | 8. Ram® X-Grip Cell Phone Holder | 9. Wishbone Crutch Holder | 10. Lap Belt | 11. Saddlebag
12. Walker Holder | 13. Rearview Mirror | 14. Weather Cover | 15. Military Patches | 16. Safety Flag

ACCESSORIES





MODEL NAME / NUMBER	ZERO TURN 8 ¹ / JAZZY ZTS	ZERO TURN 10 ¹ / S7102T	WRANGLER® / MV600	iRIDE™ / S25	GO-GO ELITE TRAVELLER® ¹ PLUS / SC53 (3W), SC54 (4W)	GO-GO ELITE TRAVELLER® ¹ / SC40E (3W), SC44E (4W)	GO-GO® SPORT ¹ / S73 (3W), S74 (4W)	GO-GO® FOLDING ¹ / S19WH1001 (SLA) S19WH1002 (Lithium Ion)	GO-GO® LX WITH CTS ¹ / SC50LX (3W), SC54LX (4W)	GO-GO® ULTRA X ¹ / SC40X (3W), SC44X (4W)
WEIGHT CAPACITY	325 lbs.	400 lbs.	350 lbs.	250 lbs.	300 lbs.	300 lbs.	325 lbs.	250 lbs.	300 lbs.	260 lbs.
MAXIMUM SPEED ²	Up to 5.9 mph	Up to 7.2 mph	Up to 11 mph at 350 lbs. Up to 11.4 mph at 200 lbs.	Up to 3.7 mph	Up to 4.5 mph	Up to 4 mph	Up to 4.7 mph	Up to 4 mph	Up to 4.4 mph	Up to 4 mph
GROUND CLEARANCE ^{3,8,9}	1.3" at motor, 1.5" at anti-tip, 3.5" (side deck)	2" at motor	5.5" at motor	1.3" at motor	1.7" at transaxle, 0.83" at anti-tip, 2.6" at side deck	1.38" at motor, 0.75" at anti-tip, 2.25" at mid frame	1.5" at motor, 0.75" at anti-tip	2.3"	1.5" at motor / 0.75" at anti-tip	3W: 3.1" at deck 4W: 3" at deck
TURNING RADIUS ⁵	37.25"	43"	94.75"	33.5"	3W: 34.5" 4W: 45.75"	3W: 33" 4W: 44"	3W: 36" 4W: 50.5"	32"	3W: 33" 4W: 44.75"	3W: 33" 4W: 44"
OVERALL LENGTH / WIDTH ^{3,7}	45.75" / 20.75"	48" / 24.25"	61.25" / 32.5"	38" / 19.25"	3W: 40", 4W: 40.25" / 21.25"	3W: 37" 4-W: 39.5" / 19.5"	3W: 42" 4W: 43" / 21.25"	37" / 19"	3W: 39.75" 4W: 41" / 20.5"	3W: 37" 4W: 39.5" / 19.5"
SEAT-TO-GROUND HEIGHT	21.5" - 23.5"	23" - 25"	30" - 31.5"	18" – 20"	20.5" - 22"	18.25"-20.25"	20.5" - 22"	20.5"	20.75" - 22.25"	18.25" - 20.25"
SEAT-TO-DECK HEIGHT	17" - 18.5"	18"-20"	18.75" - 20.25"	8.75" – 10.75"	16.25"- 17.75"	15.75"-17.25"	17" - 18"	13"	16.25" - 17.25"	15.75" - 17.25"
FRONT WHEELS (TIRES)	2"x7" solid	9" solid	14.5" tubeless pneumatic	6.5" solid	3W: 2" x 8" 4W: 2" x 8"	3W: 2" x 8" solid 4W: 2" x 7" solid	3W: 3" x 8" solid 4W: 3" x 9" solid	1.5" x 7" solid	3W: 2" x 8" solid 4W: 2" x 7" solid	3W: 2" x 8" solid 4W: 2" x 7" solid
REAR WHEELS (TIRES)	2.5"x9" solid	10.75" solid	14.5" tubeless pneumatic	6.75" solid	3" x 9" solid	2.5" x 8" solid	3" x 9" solid	2" x 8" solid	2.5" x 8" solid	2.5" x 8" solid
SUSPENSION	Front and rear (CTS)	Front and rear (CTS)	Front and rear (CTS)	N/A	N/A	N/A	N/A	N/A	Front and rear (CTS)	N/A
RANGE PER CHARGE ⁵	Up to 11 miles (at 325 lbs.) Up to 13.1 miles (at 200 lbs.)	Up to 18 miles (at 400 lbs.) Up to 24 miles (at 200 lbs.)	Up to 18.6 miles (75 AH @ 350 lbs.) Up to 22.6 miles (75AH @ 200 lbs.) Up to 24.2 miles (100 AH/ travel battery @350 lbs.) Up to 29.4 miles (100 AH/ travel battery @200 lbs.)	Up to 10.3 miles at 250 lbs. Up to 12.7 miles at 200 lbs.	3W: Up to 14.5 miles 4W: Up to 14.5 miles	3W: Up to 8.5 miles (12 AH) 4W: Up to 8.2 miles (12 AH) 3W: Up to 12.4 miles (18 AH) 4W: Up to 12 miles (18 AH)	3W: Up to 15.9 miles 4W: Up to 14.5 miles	Up to 13 miles (lithium ion) Up to 9.3 miles (sealed lead acid)	3W: Up to 8 miles 4W: Up to 8.6 miles	3W: Up to 6.9 miles 4W: Up to 7.2 miles
HEAVIEST PIECE ³	48 lbs. (rear section)	75 lbs. (front section)	219 lbs. (unit does not disassemble)	32 lbs. (Base only)	3W: 32 lbs. (rear section) 4W: 34 lbs. (front section)	3W: 29 lbs. (front section) 4W: 34 lbs. (front section)	3W: 35 lbs.(front section) 4W: 44 lbs. (front section)	51 lbs. (unit does not disassemble, unit weight without battery pack)	3W: 34 lbs. (front section) 4W: 37 lbs. (front section)	3W: 29 lbs. (front section) 4W: 34 lbs. (front section)
STANDARD SEATING	Type: Compact, foldable Weight: 27 lbs. Material: Black vinyl Dimensions: Width: 17", Depth: 17", Back height: 14"	Type: Foldable, deluxe contour Weight: 34 lbs. or Foldable, compact Weight: 30 lbs. Material: Black vinyl Dimensions: Width: 18", Depth: 17"	Type: Pillow-top, limited recline Weight: 48.6 lbs. Material: Black vinyl Dimensions: Width: 18.5", Depth: 16.5"	Type: Compact, removable Weight: 10 lbs. Material: Molded plastic Dimensions: 14.75" x 13.5"	Type: Compact, foldable Weight: 23 lbs. Material: Black vinyl Dimensions: Width: 17", Depth: 17"	Type: Compact, foldable Weight: 23 lbs. Material: Black vinyl Dimensions: Width: 17", Depth: 17"	Type: Compact, foldable Weight: 23 lbs. Material: Black vinyl Dimensions: Width: 17", Depth: 17"	Type: Compact, foldable Material: Gray vinyl Dimensions: Width: 16", Depth: 13"	Type: Compact, foldable Weight: 23 lbs. Material: Black vinyl Width: 17", Depth: 17"	Type: Compact, foldable Weight: 23 lbs. Material: Black vinyl Dimensions: Width: 17", Depth: 17"
DRIVE SYSTEM	24-volt DC motor, rear-wheel drive, dual in-line	24-volt DC motor, rear-wheel drive, dual in-line	24-volt DC motor, rear-wheel drive, dual in-line	36 volt brushless hub motor	24-volt DC motor, sealed transaxle	24-volt DC motor, sealed transaxle	24-volt DC motor, sealed transaxle	4-Pole 270w	24-volt DC motor, sealed transaxle	24-volt DC motor, sealed transaxle
BRAKING SYSTEM	Regenerative and electromechanical	Regenerative and electromechanical	Regenerative and electromechanical	Manual, dual disc	Regenerative and electromechanical	Regenerative and electromechanical	Regenerative and electromechanical	Regenerative and electromechanical	Regenerative and electromechanical	Regenerative and electromechanical
BATTERIES ⁴	(2) 22 AH	(2) 40 AH	(2) 75 AH (standard) (2) 100 AH (optional)	(1) 36 volt, 5.2 AH lithium ion	(2) 18 AH	(2) 12 AH (standard) or (2) 18 AH (optional)	(2) 18 AH	(2) 12 AH (standard, sealed lead acid) or (1) 11.5 AH (optional, lithium ion)	(2) 12 AH	(2) 12 AH
BATTERY WEIGHT ^{5,6}	Battery box 32 lbs., 13.5 lbs. each	30.25 lbs. each	75 AH – 48 lbs. each; 100 AH – 60 lbs. each	4.5 lbs.	14 lbs. each ⁷ 29.5 lbs. per battery pack	9 lbs. each (12 AH), 20 lbs. per battery pack 14 lbs. each (18 AH), 29.5 lbs. per battery pack	14 lbs. each 29.5 lbs. per battery pack	9 lbs. each (12 AH), 5 lbs. (11.5 AH)	9 lbs. each (12 AH), 20 lbs. per battery pack	9 lbs. each 20 lbs. per battery pack
BATTERY CHARGER	Off-board, 3.5 amp	Off-board, 3.5 amp	Off-board, 8 amp	36 volt lithium, off-board	Off-board, 2 amp	Off-board, 2 amp	Off-board, 2 amp	Off-board, 2-amp (sealed lead acid) or Off-board, 2.5-amp (lithium ion)	Off-board, 2 amp (12 AH)	Off-board, 2 amp
HCPCS CODES	K0899*	K0899*	N/M (non-medical)	N/M (non-medical)	K0806 (Group 2, Standard)	K0800 (Group 1, Standard)	K0801 (Group 1, Heavy Duty)	K0899*	K0899*	K0899*
WARRANTY	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Batteries: 6-months	Frame: lifetime limited Drive motors: 2-year limited Electronics: 2-year limited Batteries: 6-months	Frame: lifetime limited Drive motors: 2-year limited Electronics: 2-year limited Batteries: 6-months	Frame: 2-year limited Drivetrain: 1-year limited Electronics: 1-year limited Battery: 6-months	Frame: lifetime limited 2-year limited Electronics: 2-year limited Batteries: 6-months	Frame: lifetime limited 2-year limited Electronics: 2-year limited Batteries: 6-months	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Batteries: 6-months	Frame: 2-year limited Controllers: 18-month limited Electronics: 1-year limited Batteries: 6-months	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Batteries: 6-months	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Batteries: 6-months

[illegible]

PRIDE[®]

MOBILITY

LIVE YOUR BEST[®]

(US) 800-800-8586 | (CAN) 888-570-1113 | pridemobility.com |    



Note: Transport of the iRIDE™ item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your iRIDE™ and its lithium-ion battery are permitted.

1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2. Varies with user weight, terrain type, battery amp- hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
3. Varies with user weight, terrain type, battery amp- hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
4. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
5. Deep cycle (AGM or gel cell type recommended).
6. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
7. Battery weight may vary based on manufacturer.
8. Width of base (not including seat).
9. Does not account for anti-tip clearance.
10. At midframe.
11. Product is available exclusively to Authorized Pride Brick & Mortar Providers. No internet sales are permitted. Available in store only. Other conditions apply.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

N/M: This product is not a medical device and is not intended to assist, treat, diagnose, or alleviate any medical condition or disability.

*K0899 – power mobility device, not coded by dme pdac or does not meet criteria

**Not all accessories are available for all products. Please speak with your Provider for details and availability.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

