

Phoenix HD 4

HCPCS: K0801-1,K0801



FIGURE A




FIGURE B

- Quick and easy disassembly
- Includes interchangeable color panels in Red and Blue (Figure A)
- Armrests are padded and adjustable
- Quick connect batterie
- Ergonomic throttle control
- Height-adjustable swivel seat with two-tone upholstery
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust tiller
- Anti-tip wheels
- Headlight
- Swivel seat with fold-down backrest
- Easy free wheel operation

PRODUCTS

Description

PHOENIXHD4	Phoenix HD 4Whl Travel POV 17.5" Seat
PHOENIXHD4-20	Phoenix HD 4Whl Travel POV 20" Seat

SPECIFICATIONS

Description

Base Weight	77 lbs.
Batteries (pair)	12V x 20AH
Battery Weight (pair)	31 lbs.
Brakes	Electromagnetic
Charger	1.5A Offboard
Climbing Angle	6°
Dimensions	41.5" (L) x 22.5" (W)
Floor to Seat Height	20" - 22"
Front Wheels	9" x 3" Flat-free
Ground Clearance	3.75"
HCPCS	K0801
Heaviest Piece	46 lbs. (Front Section)
Maximum Range	15 miles
Motor	24V x 350W x 3600 rpm
Rear Wheels	9" x 3" Flat-free
Seat Dimensions	17.5" (W) x 17" (D)
Seat Type	Folding

SPECIFICATIONS	Description
Seat Weight	21 lbs.
Top Speed	4 mph
Total Weight (excl. hand control)	129 lbs.
Turning Radius	54"
Warranty on Batteries	12 Months
Warranty on Drive Train Components	24 Months
Warranty on Electronic Controller	24 Months
Warranty on Frame	Lifetime
Weight Capacity	350 lbs
Weight Capacity (13046 13046-HD)	350 lbs.
Wheel Base Dimensions (L x W)	30.7" x 19.1" 30.7" x 19.1"

ACCESSORIES	Description
SF8010	Oxygen Holder Set, L Type,1/ea, EA
SF8030	Crutch Holder Set, SType, 1/ea, EA
SF8020	Rear Basket Set, LType, 1/ea, EA
SF8050	Sun Shade for Scooter,1/ea, EA

PARTS	Description
S31004GT	HandleBarCover, Digital Batt Gauge,1/ea, EA
S3KEY	Key for Phoenix,Daytona3,4,1/ea, EA
C20-007-00200	PlasticHandGrip,BcatX/PhoenixHD,Mav,1/ea, PR
S31130-4HD	Rear Frame, Phoenix 4 HD, 1/ea, EA
LRM412123	21AH Battery,POV, 1/ea, EA
BC35051	BasketHook Set,Phoenix,Scout,Ventura,1/e, EA
S31052	Armrest Knob POV 1/ea, EA
S31026	Anti Tip Wheel,Phoenix/Bobcat,1/ea, EA
S350016910	Rear Axle Bolt,Phoenix, Bobcat X, 1/ea, EA
LRM412103-11	Speed Control Knob, EA
S35048MFR	Mud Flap, Right, Phoenix HD, 1/ea, EA
S350169706-HD	BatteryCaseOnly,Phoenix-3HD,4HD,1/ea, EA
S350169706-HDCP	Battery Box Charge Port, Phoenix HD,1/ea, EA
S35048MFL	Mud Flap, Left, Phoenix HD, 1/ea, EA
S31031K	FreeWheel Lever Knob,Phoenix,1/ea, EA
S350104	Front Bumper,Phnx4,1/ea, EA
S350103	Connecting Rod,Phnx4/HD4,1/ea, EA
S350108	Front Frame,Phnx4/HD4,1/ea, EA
S350107	Front Shroud,Phoenix4/HD4,1/ea, EA
S3501BL-1	FrontDecor,Blue,Left,Phoenix4/HD4,1/ea, EA