



Pilot 2410

Features



Pilot 2310

- Available in **3- and 4-wheel** models.
- **22 mile maximum range** excellent for indoor and outdoor use.
- Standard comfortable and durable **Captain's Seats** feature deep-contoured cushions and backs for maximum comfort and support.
- Captain's Seats swivel and feature a **fold-down backrest** with 25 degrees of recline and an adjustable height headrest.
- Seats have four **height adjustments**, from 20" to 23".
- Seat can be moved forward or backward using the auto-style seat slides.
- **Padded, flip-up armrests** allow for easy transfers and user access. Armrests are also removable and height and width adjustable.
- Adjustable angle, **delta tiller** allows the user to position the controls at the most comfortable angle and depth.
- **Ergonomic controls** make it easier for those with limited strength or dexterity to handle the scooter.
- **Anti-tip wheels** for greater safety.
- Standard **safety package** includes headlight, horn, turn signals and dual rear view mirrors.



PilotTM Series

Standard Mobility Scooters

ActiveCare
driveTM Mobility

Pilot™ Series Standard Mobility Scooters

Specifications

Pilot 2310

HCPCS	K0807
Weight Capacity	350 lbs.
Maximum Speed	6 mph
Maximum Range	22 miles
Turning Radius	53"
Ground Clearance	3"
Overall Length	48"
Overall Width	24.5"

Seat Dimensions

18"(w) x 18"(d) x 21"(h) Captain's Seat
20"(w) x 19"(d) x 21"(h) Captain's Seat

Floor to Seat Height	20"-23"
Motor	24V x 350W x 4200rpm
Batteries	Two 12V x 35AH (U1)
Battery Charger	5 Amp, Offboard
Brakes	Electromagnetic
Wheels	10" x 3" Flat Free
Base Weight	84 lbs.
Battery Weight (pair)	48 lbs.
Seat Weights	
18" Captain's Seat	49 lbs.
20" Captain's Seat	55 lbs.

Pilot 2410

HCPCS	K0807
Weight Capacity	350 lbs.
Maximum Speed	6 mph
Maximum Range	22 miles
Turning Radius	54"
Ground Clearance	3"
Overall Length	48"
Overall Width	24.5"

Seat Dimensions

18"(w) x 18"(d) x 21"(h) Captain's Seat
20"(w) x 19"(d) x 21"(h) Captain's Seat

Floor to Seat Height	20"-23"
Motor	24V x 350W x 4200rpm
Batteries	Two 12V x 35AH (U1)
Battery Charger	5 Amp, Offboard
Brakes	Electromagnetic
Wheels	10" x 3" Flat Free
Base Weight	90 lbs.
Battery Weight (pair)	48 lbs.
Seat Weights	
18" Captain's Seat	49 lbs.
20" Captain's Seat	55 lbs.

Inventory Numbers

PILOT2310BL18CS	3 Wheel, 18" Captain's Seat, Blue
PILOT2310BG18CS	3 Wheel, 18" Capt. Seat, Burgundy
PILOT2310BL20CS	3 Wheel, 20" Captain's Seat, Blue
PILOT2310BG20CS	3 Wheel, 20" Capt. Seat, Burgundy
PILOT2410BL18CS	4 Wheel, 18" Captain's Seat, Blue
PILOT2410BG18CS	4 Wheel, 18" Capt. Seat, Burgundy
PILOT2410BL20CS	4 Wheel, 20" Captain's Seat, Blue
PILOT2410BG20CS	4 Wheel, 20" Capt. Seat, Burgundy

Accessories (sold separately)

Carriers & Holders

- ☐ Crutch/Cane Holder
- ☐ Cylinder Tank Carrier
- ☐ Rear Basket

Carry Bags

- ☐ Armrest Bag
- ☐ Backpack

Other Accessories

- ☐ Scooter Trailer
- ☐ Scooter Cover
- ☐ Cup Holder

Warranty

Lifetime on frame, 14 months on electronic controller and drive train components, six months on batteries.

Color Options



Burgundy



Midnight Blue