



# XtendaCare®

Providing Pressure Reduction and Comfort  
with options tailored to your needs.



## Standard XtendaCare®



## XtendaCare® TE ( transfer edging )

Provides ultra-firm support along each side of the mattress for edge of bed sitting and supported transferring.



## XtendaCare® SE ( secure edging )

Raised edging directs residents away from the edge of the bed, helping prevent side-rail entrapment and accidental falls. Designed to enhance resident safety without compromising dignity.



## Stretch Fabric...Why we don't

When subjected to the real-world conditions of weight, heat and moisture, stretchable fabrics can allow moisture vapor to be drawn into the cushion, causing infection control issues. Comfortex's non-shearing Stratos nylon **DOES NOT** stretch and has an unmatched 20-year history of protection and performance.



### Low Shear – High Performance Top Cover

Patients and Nursing Home Residents are protected by Comfortex's exclusive anti-shear nylon, treated to be both fluid & vapor proof. Sewn with a slight billowiness, allows each person's unique shape to descend safely and comfortably into their sleep surface.

### High Resiliency Support

HR foams have a stronger, more stable cell structure to provide years of cushioning comfort. The XtendaCare's soft top, over supportive layers and strong laminations, maximize comfort, pressure redistribution and durability.

**Comfortex**  
**Innovative Healthcare Surfaces**  
**1-800-445-4007 [www.comfortexinc.com](http://www.comfortexinc.com)**