

The Pool Access Chair and Heavy Duty Pool Access Chair both provide for safe, convenient and easy access for people with disabilities to pools with ramps or beach entries.

The Pool Access Chair features a weight capacity of 250 lbs while our Heavy Duty Pool Access Chair features an extra wide seat and a 400 lb weight capacity.

Both provide an inexpensive way to meet ADA compliance standards at your commercial facility and make your facility more user friendly to all users. These chairs can also be used effectively in shower and locker rooms.

Features:

- 250 lb. (113 kg) weight capacity - Pool Access Chair
- 400 lb. (181 kg) weight capacity - Heavy Duty Chair
- Durable non-corrosive PVC construction
- Stainless steel hardware
- All plastic rear wheels with rubber tread
- Anti-tip design
- Swing out armrest for lateral transfers
- One year warranty

Options:

- *Pool Access Chair with Optional Mesh Seat*
Item #F-11PACSMS
- *Chest strap* - Item #F-426CS
- *Seat belt* - Item #F-11PACSB



PROUDLY SOLD BY: