

**Item# F-03PEZ**



The Power EZ™ Lift is a great economical choice for individuals needing to gain access to their pool or spa.

The Power EZ™ Lift provides a fully automatic actuator to raise and lower with 360° manual rotation. The lift can accommodate spa walls up to 26 inches (66cm) high and can be set back up to 20 inches (50cm) from the pool edge.

The Power EZ™ Lift now comes with a rotation bar to allow the attendant to rotate the unit with greater ease.

**FEATURES:**

- 400 Pound (181 kg) weight capacity
- Rechargeable battery operation with charger
- Clears walls up to 26" (66 cm) high
  - Durable mesh sling- standard size: 38"x20"
  - Adjustable stainless steel chains
  - Limited lifetime structural warranty
  - 5 year pro-rated electronics warranty

**OPTIONS:**

- *Hard seat option - Item #F-0320HSO*
- *Sling, youth - Item #F-34EZCS / 34"x18" (86x46cm)*
- *Sling, extra long - Item #F-037EZLS / 50"x28" (127x71cm)*
- *Sling, standard with rigid insert & seat belt - Item #F-036RSS*
- *Anchor, standard in-ground (additional) - Item #F-033EZA*
- *Anchor, wood deck - Item #F-30600*
- *Anchor, surface mount - Item #F-035EZSM*
- *Upgrade Pack - Item #F-03PEZAP - Includes battery, rigid sling with seat belt, sling bars and cover*
- *Protective Cover - Item #F-440EZC*
- *Hard seat cover - Item #F-440EZHSC*



**PROUDLY SOLD BY:**