

REFER TO PAGE 182-190
FOR LIFT ACCESSORIES
AND SLINGS.

RECOMMENDED
FOR USE WITH 2 &
4 POINT SLINGS!

Battery Powered Patient Lift with Six Point Cradle

13240 1/cs

BARIATRIC
450 LB
WEIGHT CAPACITY

FEATURES:

- Fast, easy charge from AC outlet.
- **To charge, simply plug in unit. Battery does not need to be removed to be charged. Lift will not operate while plugged in.**
- Heavy gauge construction.
- Jumbo actuator pump with emergency button can switch to manual mode lowers patient down safely.
- Audible warning when battery power falls below 50%.
- The "CHARGE" light will automatically shut off when battery is fully charged.
- **Ideal for Homecare use.**
- For use with 4 point and 2 point slings.
- Wide base provides additional stability.
- Lifting capacity: 450 lbs.
- Limited Lifetime Warranty.



NEW & IMPROVED
CRADLE

13240
with Sling
(Sold Separately)

CAN FIT UNDER
A LOW BED!

Battery Powered Patient Lift with Rechargeable, Removable Battery

13242 1/cs

BARIATRIC
450 LB
WEIGHT CAPACITY

FEATURES:

- Same features as 13240.
- Includes removable, rechargeable battery, (item #13258) external charger (item #13255) and wall mounted bracket (item #13260).
- **Ideal for institutional use.**

To avoid lift downtime - order an additional re-chargeable battery, item # 13258. One for the lift and one for the charging station.

REPLACEMENT PARTS:



**External Charger
with Power Cord**

13255 1/ea



**Wall Mounted
Bracket**

13260 1/ea



**Rechargeable
Battery**

13258 1/ea



NEW & IMPROVED
CRADLE

13242
with Sling (Sold Separately)

CAN FIT UNDER A LOW BED!

SPECIFICATIONS

| Item # | Boom Maximum Height | Boom Minimum Height | Base Width (open) | Base Width (closed) | Base Height | Base Clearance | Casters | Product Weight | Carton Shipping Weight | Weight Limit |
|--------|---------------------------|---------------------------|-------------------------|---------------------------|----------------|-------------------|---------|-------------------|------------------------------|-----------------|
| 13240 | 76" | 27.5" | 41" | 15" | 4.5" | 2" | 5" Rear | 106 lbs. | 116 lbs. | 450 lbs. |
| 13242 | 76" | 27.5" | 41" | 15" | 4.5" | 2" | 5" Rear | 106 lbs. | 116 lbs. | 450 lbs. |

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

patient room

drive
WWW.DRIVEMEDICAL.COM

To order call toll free: 877.224.0946