

Jazzy® 1122



JAZZY 1122 shown in Candy Apple Red with deluxe contoured high-back seat, PG Remote Plus controller and foot platform



At first glance, the Jazzy® 1122 catches your eye with its sleek styling and sporty edge. Yet, the 1122 is so much more. The 1122 features standard Active-Trac® and rear caster suspension and a high-speed/high-torque motor package that allows an optimal speed of 6 mph and a weight capacity of 300 lbs. In addition to these impressive performance-enhancing features, the 1122 also offers easy battery access through the front of the chair without the removal of the shroud or seating system and a height and depth-adjustable solid foot platform. The Jazzy 1122 is truly the most advanced power chair available today!

SPECIFICATIONS

Weight capacity:	300 lbs.		
Maximum speed:	Up to 6 mph		
Ground clearance:	3.25"		
Turning radius:	21"		
Overall length:	40.5"		
Overall width:	25.4"		
Wheels:	Drive: 14" pneumatic Casters: 8" solid Anti-tips: 6" solid		
Suspension:	Active-Trac® and rear caster		
Drive train:	Two-motor, mid-wheel drive; high-speed motor package (std.)		
Braking system:	Intelligent braking (electronic, regenerative disc brakes)		
Std. Electronics:	70 amp, PG VSI controller		
Specialty controls:	N/A		
Battery charger:	Off-board, 8 amp (standard)		
Per-charge range:	Up to 25 miles		
Battery requirements:	(2) 12 volt, deep cycle		
Battery size:	70AH, Group 24		
Transport weight			
Battery weight:	53.5 lbs. each		
Base weight:	129 lbs.		
Std. seat weight:	42.5 lbs. (high-back)		
Std. Seat dimensions:	18" x 18"		
Warranty:	Lifetime limited warranty on frame 2-year limited warranty on electronics 18-months limited warranty on drive motors		
Colors:			
Candy Apple Red	Viper Blue	Onyx Black	Champagne

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

FEATURES

- High-speed motor package—top speed of 6 mph
- Easy front-access batteries
- Powerful, patented, two-motor, mid-wheel drive design provides superb power and acceleration
- Active-Trac® and rear caster suspension for outstanding outdoor performance
- Adjustable front anti-tips to meet your preferences and terrain challenges

OPTIONS

- Solid seat pan
- Power elevating seat
- Solid tires
- Lap belt
- Cup holder
- Heavy-duty dust cover
- Swing-away joystick
- Swing-away footrests
- Elevating leg rests
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Rear basket
- Specialty controls
- Angle-adjustable footplates
- 70° and 80° heavy-duty drop-in footrests
- Saddlebag
- 5 amp, on-board/8 amp, off-board charger combo