

Jazzy® 1170 XL Plus  
shown in Candy Apple Red



# Jazzy®

## BARIATRIC SERIES

Jazzy® bariatric power chairs deliver the high-performance maneuverability of mid-wheel drive design with heavy-duty construction for maximum durability. With features like enhanced seating sizes, powerful drivetrains, and weight capacities up to 650 lbs., Jazzy bariatric models are the perfect choice for the larger individual.

Jazzy® 1420  
shown in Sunburst Orange



Jazzy® 1650  
shown in Onyx Black



(US) 800-800-8586 • (Canada) 888-570-1113  
Visit us on the web: [www.pridemobility.com](http://www.pridemobility.com)

**Pride®**  
Mobility Products Corp.



Specifications	Jazzy® 1170 XL Plus	Jazzy® 1420	Jazzy® 1650
<b>Weight capacity</b>	500 lbs.	500 lbs.	650 lbs.
<b>Maximum speed</b>	Up to 5 mph	Up to 4.5 mph	Up to 5 mph
<b>Ground clearance</b>	4"	4.5"	4.25"
<b>Turning radius</b>	26.75"	25.75"	26.25"
<b>Overall length</b>	45.25"	46"	47"
<b>Overall width</b>	30"	27.25"	30.375"
<b>Drive Wheels</b>	16" solid or pneumatic	16" solid or pneumatic	16" solid or pneumatic
<b>Front anti-tips</b>	8" solid	8" solid	8" solid
<b>Rear casters</b>	9" solid or pneumatic	9" solid or pneumatic	9" solid or pneumatic
<b>Suspension</b>	Active-Trac® suspension	Active-Trac® w/ rear caster suspension	Limited
<b>Drive train</b>	Two-motor, Mid-Wheel Drive	Two-motor, Mid-Wheel Drive	Two-motor, 4-pole, in-line, Mid-Wheel Drive
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)	Intelligent braking (electronic, regenerative disc brakes)	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	70 amp, PG VSI controller	70 amp, VR2 controller	100 amp, PG Remote Plus controller
<b>Std. Battery charger</b>	Off-board, 8 amp	Off-board, 8 amp	Off-board, 8 amp
<b>Per-charge range</b>	Up to 15 miles	Up to 20 miles	Up to 15 miles
<b>Battery requirements</b>	(2) 12 volt, deep cycle (gel cell or AGM)	(2) 12 volt, deep cycle (gel cell or AGM)	(2) 12 volt, deep cycle (gel cell or AGM)
<b>Battery size</b>	Group 24	Group 24	Group 24
<b>Battery weight</b>	53.5 lbs. each	53.5 lbs. each	53.5 lbs. each
<b>Base weight</b>	149.5 lbs.	141.5 lbs.	176 lbs.
<b>Standard seat weight</b>	42.5 lbs. (high-back)	42.5 lbs. (high-back)	75 lbs. (high-back)
<b>Max. seat dimensions</b>	28" x 26"	24" x 24"	32" x 28"
<b>Warranty</b>	Lifetime limited warranty on frame; 2-year limited warranty on electronics; 18-month limited warranty on drive motors	Lifetime limited warranty on frame; 2-year limited warranty on electronics; 18-month limited warranty on drive motors	Lifetime limited warranty on frame; 2-year limited warranty on electronics; 18-month limited warranty on drive motors
<b>Colors</b>	Candy Apple Red, Viper Blue, Onyx Black, Champagne	Candy Apple Red, Viper Blue, Onyx Black, Champagne, Sunburst Orange	Candy Apple Red, Viper Blue, Onyx Black, Champagne
<b>Available SADMERC packages</b>	2VHD-C, 2VHD-SS, 3VHD-C, 3VHD-SS	2HD-C, 2HD-SS, 2VHD-C, 2VHD-SS, 3VHD-C, 3VHD-SS	K0899, 3EHD-C, 3EHD-SS

800-800-8586 (US) ■ 888-570-1113 (Canada)  
Visit us on the web: [www.pridemobility.com](http://www.pridemobility.com)

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %.

**PRIDE**®  
PERFORMANCE MOBILITY  
— LIVE YOUR BEST —