





The new Legend Scooter combines the latest style and performance technology at a great value. Standard features include a top speed of up to 5.75 mph, 350 lbs. maximum weight capacity and an easy lock-up system for quick disassembly. Plus, the new lower dual-bulb headlight offers excellent pathway illumination. The Legend's sleek styling and advanced features will have you turning heads every time you turn a corner.

Experience the new Legend,® and life, like never before!







## **Specifications**

•	
Model number(s)	3-wheel: SC3000
	4-wheel: SC3400
Weight capacity	350 lbs.
Maximum speed	Up to 5.75 mph
Ground clearance	3.5"
Turning radius	3-wheel: 41"
	4-wheel: 57"
Overall length	3-wheel: 45.5"
	4-wheel: 47"
Overall width	3-wheel: 24"
	4-wheel: 24"
Tires	
Front	3.5" x 10"
Rear	4" x 10"

111 05	
Front	3.5" x 10"
Rear	4" x 10"
Anti-tips	Rear
Standard seating	CRS deluxe reclining high-back
	with headrest and sliders
Suspension	Rear
Drive train	Sealed transaxle,
	24-volt DC motor
Dual braking system	Electronic, regenerative
	and electromechanical
Battery charger	On-board, 3 amp
Battery requirements	(2) 12V deep-cycle (AGM or gel cell
	type recommended)

Weight	
Battery	24.5 lbs. each
Total w/o batteries	3-wheel: 155 lbs.
	4-wheel: 166 lbs.
Heaviest piece when disassembled	61 lbs. (rear)

## **Colors**







Champagne Viper Blue Candy Apple

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

Red

## **Benefits**

- State-of-the-art rear suspension
- Top speed of up to 5.75 mph
- Maximum weight capacity of 350 lbs.
- Ultra-simple lock-up for quick disassembly
- Advanced new reclining high-back seat with height-adjustable headrest
- Easily accessible tie down points (for transport of unoccupied scooter)
- Reduced length:
  - 3-wheel: 45.5" / 4-wheel: 47"
- New dual-bulb, adjustable lower headlight for optimal pathway illumination
- Easy access to rear electronics
- Tiller-mounted, recessed "marker" headlight
- Tight turning radius:
  - 3-wheel: 41" / 4-wheel: 57"
- Reinforced front nose section for durability
- 3-year limited warranty

## **Options**

- Power elevating seat
- Forearm crutch holder
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Rear basket
- Safety flag
- Walker holder
- Heavy-duty dust cover
- Saddlebag
- Lap belt



