

Model number:	SC151	SC155
Weight capacity:	300 lbs.	300 lbs.
Maximum speed:	Up to 4.25 mph	Up to 4.25 mph
Ground clearance:	3"	3"
Turning radius:	39"	54"
Overall length:	41"	44"
Overall width:	21"	21"
Wheels: Front: Rear: Anti-tips:	2" x 8" pneumatic 2.5" x 9" pneumatic Rear	2" x 8" pneumatic 2.5" x 9" pneumatic Rear
Suspension:	N/A	N/A
Transport weight Total w/o batteries:	93 lbs.	116 lbs.
Heaviest piece when disassembled:	45 lbs. (rear section)	45 lbs. (front or rear section)
Drive train:	Sealed transaxle, 24 volt DC motor	
Dual braking system:	Electronic, regenerative and electromechanical (electronic, regenerative disc brakes)	
Power elevating seat:	N/A	
Std. Seating:	Foldable molded plastic	
Battery charger:	On-board, 3 amp (standard)	
Per-charge range:	Up to 25 miles	
Battery requirements:	(2) 12V deep cycle (AGM or gel cell type recommended)	
Battery size:	U-1 or optional 17AH	
Warranty:	3-year limited warranty	

Benefits

- Controls at your fingertips
- Easy-to-use external freewheel lever
- Aluminum alloy mags for great style
- Front basket to hold all your belongings
- Seat offers height adjustment and swivel with flip-up armrests for comfortable positioning
- Adjustable tiller for maximum comfort
- On-board accessory bracket
- Easily accessible charger port and circuit breaker
- Easily disassembles for transport

Options

- Forearm crutch holder
- Lap belt
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- · Heavy-duty dust cover
- Saddlebag
- Quad cane holder
- Clip, cane attachment
- Flat-free tires



- · Backpacker
- · Scooter double entry
- · Scooter single entry
- · Full Platform



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

Viper Blue

Colors:



Candy Apple

Red





STYLISH DESIGN & PREMIUM PERFORMANCE AT AN EXCELLENT VALUE.

RALLY®

The Rally® is sure to cause a stir with its futuristic styling, technologically advanced performance and overall great value—a combination irresistible and unmatched in the small scooter class. With the Rally, you'll stake your claim on a lifestyle of fun and freedom.



