



Victory® 3-wheel Victory® 4-wheel

Model number:	SC1600	SC1700
Weight capacity:	350 lbs.	350 lbs.
Maximum speed:	Up to 5 mph	Up to 5 mph
Ground clearance:	4"	4"
Turning radius:	41"	58"
Overall length:	47"	49"
Overall width:	24"	24"
Wheels:		
Front:	3.5" x 10" pneumatic	3.5" x 10" pneumatic
Rear:	3.5" x 10" pneumatic	3.5" x 10" pneumatic
Anti-tips:	Rear	Rear
Suspension:	N/A	N/A
Transport weight		
Total w/o batteries:	110.5 lbs.	116.5 lbs.
Heaviest piece when disassembled:	50 lbs. (rear section)	50 lbs. (rear section)
Drive train:	Sealed transaxle, 24 volt DC motor	
Dual braking system:	Electronic, regenerative and electromechanical (electronic, regenerative disc brakes)	
Power elevating seat:	N/A	
Std. Seating:	Foldable molded plastic	
Battery charger:	On-board, 3 amp (standard)	
Per-charge range:	Up to 25 miles	
Battery requirements:	(2) 12V deep cycle (AGM or gel cell type recommended)	
Battery size:	U-1	
Warranty:	3-year limited warranty	

Colors:



Benefits

- Single-point lock-up for easy scooter disassembly
- Ultra-sleek, flush-mount oval headlight with mirror finish
- Restyled sporty alloy mag wheels
- Heavy-duty front (4-wheel) and rear polyurethane foam rubber bumpers
- Water and shock-resistant Guest battery charger
- Convenient, externally mounted charger ammeter
- Standard front basket
- On-board accessory bracket
- Easily accessible charger port and circuit breaker

Options

- Forearm crutch holder
- Lap belt
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Heavy-duty dust cover
- Saddlebag
- Solid tires
- Quad cane holder
- Clip, cane attachment

SilverStar Lifts A Division of Pride Mobility Products® Corp.

- Backpacker
- Scooter double entry
- Scooter single entry
- Full Platform



Backpacker
These lifts are compatible with the Victory®.

Double entry

Pride®
Mobility Products Corp.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

Victory®

**STYLISH DESIGN & PREMIUM
PERFORMANCE AT AN
EXCELLENT VALUE.**

VICTORY®

Determined focus led to the ideal balance of hot style, solid performance and exceptional value found in the Victory® Scooter—a focus that provides you with a truly fun-loving scooter equipped standard with loads of luxury features at an affordable price. Live large even on a tight budget with a Pride® Victory.



Pride®
Mobility Products Corp.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

VICTORYSS9-21-05