

Top End® Pro™ 2 Wheelchair



PROVIDER / CUSTOMER INFORMATION:			
Provider Information:		End User Information:	
Date:	Name:	Height:	
OMS Acct #:	Ship to Address:	Weight (lbs.):	
Contact:	Address Line 2:	Leg Length (in.):	
Phone #:	City:	Hip Width (in.):	
Fax/Email:	State:	Zip:	Phone #:
<input type="radio"/> Ship to Provider	<input type="radio"/> Ship to End User		Email:

ORDER REQUEST TYPE:		
<input type="radio"/> Quote	<input type="radio"/> Order	PO # <input type="text"/>

BASE MODEL:		
<input checked="" type="radio"/> 10014	Top End Pro 2 - All Sport	\$2,714

Assembly: <i>(see definitions)</i>		
<input checked="" type="radio"/> 10048	Assembled	\$0
Shipping:		
Shipping Fee <i>(Continental US Only - See definitions)</i>		\$325

SIZING:		
SEAT WIDTH:		
<input type="radio"/> 15013	14" Seat Width - White Snow	\$0
<input type="radio"/> 15017	16" Seat Width - Simple Red	\$0
<input type="radio"/> 15021	18" Seat Width - Royal Blue	\$0
SEAT DEPTH: <i>(See Definitions)</i>		
<input checked="" type="radio"/> 15357	Adjustable Seat Depth	\$0
SEAT HEIGHT:		
<input checked="" type="radio"/> 15096	Adjustable Front Seat Height: 18" - 20" Adjustable Rear Seat Height 15" - 20"	\$0
BACK HEIGHT:		
<input checked="" type="radio"/> 15136	Adjustable Back Height 11" to 14"	\$0

BACK TYPE:		
<input checked="" type="radio"/> 15159	Angle Adjustable, Non-Folding Backrest	\$0

FOOTREST:		
<input checked="" type="radio"/> 60003	Clamp-On Tubular, Height Adjustable Footrest	\$0

CAMBER ANGLE:		
<input checked="" type="radio"/> 45008	15 Degree Camber	\$0

WHEELBASE: <i>(Front to Rear Axle)</i>		
<input checked="" type="radio"/> 75068	17" Wheelbase	\$0

CENTER OF GRAVITY: <i>(Front of Backrest to Center of Axle)</i>		
<input checked="" type="radio"/> 45073	Adjustable 2.5" to 5" Center of Gravity	\$0



WHEEL SETS:		
WHEEL SIZE:		
<input checked="" type="radio"/> 50191	25" (559)	\$0
WHEEL TYPE:		
<input checked="" type="radio"/> 50196	High Performance	\$0
HANDRIM TYPE:		
<input checked="" type="radio"/> 50297	6 Tabs, Aluminum, Anodized	\$0
HANDRIM MOUNTING:		
<input checked="" type="radio"/> 50308	Standard Mount	\$0
TIRES:		
<input checked="" type="radio"/> 50216	Primo Grey	\$0

SIDE GUARDS:		
<input checked="" type="radio"/> 20014	Screw-On, Aluminum Side Guards	\$0

ANTI-TIP:		
<input checked="" type="radio"/> 20030	Single Telescoping, Removable, Height Adjustable Anti-Tip	\$0

Top End® Pro™ 2 Wheelchair



UPHOLSTERY:		
MATERIAL & COLOR:		
<input checked="" type="radio"/> 30001	Nylon, Black Upholstery	\$0
SEAT UPHOLSTERY:		
<input checked="" type="radio"/> 30019	Screw-On Seat Upholstery	\$0
BACK UPHOLSTERY:		
<input checked="" type="radio"/> 30012	Adjustable Tension back upholstery	\$0
CUSHIONS:		
<input checked="" type="radio"/> 30005	2" Soft Seat Cushion	\$0

ACCESSORIES:		
STANDARD:		
<input checked="" type="radio"/> 25066	Calf Strap	\$0
<input checked="" type="radio"/> 25065	Rear Frame Protector	\$0
<input checked="" type="radio"/> 25060	Hook and Loop Positioning Strap	\$0
UPGRADABLE:		
POSITIONING:		
<input type="radio"/> 60030	Click Strap, Single Pad for Footrest	\$90
<input type="radio"/> 25001	Click Strap, 8" Small	\$90
<input type="radio"/> 25002	Click Strap, 10" Medium	\$97
<input type="radio"/> 25003	Click Strap, 12" Large	\$109

DEFINITIONS:	
Assembly and Shipping:	Arrives factory assembled, rear wheels and some accessories will need to be installed. Shipping rate is \$325 for all shipments in the continental US. Freight quotes for shipments outside the continental US are available upon request.
Adjustable Seat Depth:	All Top End Pro-2 All Sport Wheelchairs have adjustable seat depths. For seat widths of 10" and 12" the seat depth will be adjustable to depths of 10" and 12". For seat widths of 14", 16", and 18" the seat depth will be adjustable to depths of 13", 15" and 17".
Clamp-on Basketball Wing:	The Clamp-on Basketball Wing set components extend the protective bumper around the front of the chair. This is typically desired for use in Wheelchair Basketball. This option is not assembled on the chair and requires minimal assembly upon arrival.
Accessories:	All Additional Accessories may be shipped in a separate box than the Pro-2. There will be no additional shipping charges when accessories are purchased with a product.
Upgradable Accessories:	All Add-on accessories may be boxed separately and require self installation.

Warranty Info:

Top End warrants the frames when purchased new and unused to be free from defects in materials and workmanship for a period of two (2) years from the date of purchase from Top End or a dealer, with a copy of the seller's invoice required for coverage under this warranty. Top End warrants the upholstered materials (seat and back) and remaining components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of thirteen (13) months from date of purchase from Top End or a dealer, with a copy of the seller's invoice required for coverage under this warranty