

PROTEKT® 600 STAND

34600 | SIT-TO-STAND PATIENT LIFT

The Protekt® 600 lb. Sit-to-Stand has extra wide adjustable kneepads that will help meet a wider range of patient needs. To support the extra weight, the 600 Stand has a heavy duty boom to help lift and support the patient for smoother and safer transfers. The emergency quick release makes the 600 Stand the safest lift in any situation.

- **600 lb. weight capacity with manual emergency lowering.**
- **6 sling hooks allows for multiple sling options, safer patient handling and a wider range of transfers.**
- **Fast speed lift motor reduces time for the transfer.**
- **Battery lasts twice as long as competing brands for less downtime (24 volts).**
- **Ergonomically designed foot pedal offers hands free opening of base.**
- **Dual rear locking casters.**



- **Clip free sling hookups improve safety and reduce downtime.**
- **Full parts breakdowns and modular design on lift and electronics.**
- **Advanced electronics with auxiliary up/down lift controls.**

SPECIFICATION FOR ITEM # 34600

Weight Capacity	600 lbs.
Base Legs Length	40"
Internal Width	19"
External Width	26"
Expandable Internal Width	38"
Base Height	4.7"
Min. Lifting Height	37"
Max. Lifting Height	71"
Net Weight	134 lbs.
Shipping Weight	150 lbs.

WARRANTY:

Electronics 2 years
Frame 3 years