

# PROTEKT® ONYX

## 32475 | HYDRAULIC PATIENT LIFT

HCPCS  
E0630

The innovative Protekt® Onyx home care lift utilizes high performance hydraulics for gradual movement and provides a safe and comfortable patient transfer. This easy to operate hydraulic lift is built to last, and yet extremely affordable.

- 450 lb. weight capacity.
- 6 point spreader bar allows for multiple sling options, safer patient handling and a wider range of transfers.
- High performance hydraulics raise and lower patients gradually and safely from any stationary position.
- User friendly hand operated base opening adjusts and locks easily.
- Dual rear locking casters.
- Includes chain connections.



- Sturdy lift welding ensures safe patient transfers.
- No clips to lose ensures compliance.
- No tool required for assembly.

### SPECIFICATION FOR ITEM # 32475

Weight Capacity	450 lbs.
Base Legs Length	40.5"
Base Closed External Width	22"
Base Open External Width	42.25"
Mast Height	51"
Base Height	5.25"
Min. Lifting Height	17"
Max. Lifting Height	61"

**WARRANTY:**  
Pump 1 year Limited  
Frame Limited Lifetime

[www.proactivemedical.com](http://www.proactivemedical.com)  
[sales@proactivemedical.com](mailto:sales@proactivemedical.com)

270 Washington St., Mount Vernon, NY 10553  
Tel: 855-BE-PROACTIVE (855-237-7622)  
Fax: 845-205-6005