

Protekt™ Aire 2000

Low Air Loss/Alternating Pressure Overlay System



 **Proactive**
medical products

Medicare
Code
EO277

270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com

The Protekt™ Aire 2000 5" mattress overlay provides alternating pressure and low air loss therapies in one pump. It is the perfect overlay system to prevent and treat pressure ulcers. The pump is user friendly, reliable and quiet. The Protekt™ Aire 2000 incorporates a system of 17 individual bladders. The pump offers a 10 minute cycle time and adjustable air flow settings to meet individual patient's needs.

FEATURES AND BENEFITS

MATTRESS & COVER

- 16 laser air holes to keep patient cool and dry.
- Removable, water resistant, non-shear, anti-microbial nylon cover has a non-skid bottom to prevent slipping.
- CPR valve provides rapid deflation of mattress.
- Dimensions: 35"x80"x5".
- Weight: 9 lbs.
- Alternating cell depth: 5".

PROTEKT™ AIRE 2000 PUMP

- Low air loss and alternating pressure therapies in one system.
- Air flow dial to obtain custom patient comfort and pressure levels.
- Low pressure alarm.
- Adjustable hanging brackets.
- UL approved power cord.
- Kink resistant tubing.
- Weight Capacity: 300 lbs.
- Dimension: 12"x6"x3.5".
- Weight: 3.5 lbs.
- 18 month warranty.

80020	Protekt™ Aire 2000	5" Alternating/Low Air Loss System
80021	Protekt™ Aire 2000 Pump	Pump only for Protekt™ Aire 2000
80022	Protekt™ Aire 2000 Mattress Overlay	Mattress overlay only for Protekt™ Aire 2000
80023	Protekt™ Aire 2000 Cover	Cover only for Protekt™ Aire 2000
80024	Protekt™ Aire 2000 Bladder	Individual Bladder for Protekt™ Aire 2000



270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com