

# Protekt™ Aire 3000

Low Air Loss/Alternating Pressure Mattress System



 **Proactive**  
medical products

Medicare  
Code  
EO277

270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com

The value oriented Protekt™ Aire 3000 is an ideal mattress system for a provider with a limited budget that won't sacrifice quality. The Protekt™ Aire 3000 combines dual therapies of low air loss and alternating pressure in one unit. The durable, powerful, yet quiet user friendly pump is designed to treat pressure ulcers stages I-IV. The Protekt™ Aire 3000 mattress has 20 individual bladders that allow for deep submersion and excellent pressure redistribution. The Protekt™ Aire 3000 pump alternates in 10 minute intervals and the adjustable air flow dial provides custom pressure settings for each patients individual requirements.

## FEATURES AND BENEFITS

### MATTRESS & COVER

- 16 laser air holes to keep patient cool and dry.
- Removable, water resistant, non-shear, anti-microbial nylon cover is quilted for extra comfort.
- CPR valve provides rapid deflation of mattress.
- Non-skid cover prevents the mattress from slipping.
- Dimensions: 35"x80"x8".
- Weight: 9 lbs.
- Alternating cell depth: 8".

### PROTEKT™ AIRE 3000 PUMP

- Low air loss and alternating pressure therapies in one system.
- Air flow dial to obtain custom patient comfort and pressure levels.
- Low pressure alarm.
- Adjustable hanging brackets.
- UL approved power cord.
- Kink resistant tubing.
- Weight Capacity: 350 lbs.
- Dimension: 12"x6"x3.5".
- Weight: 3.5 lbs.
- 18 month warranty.

80030	Protekt™ Aire 3000	8" Alternating/Low Air Loss Mattress System
80031	Protekt™ Aire 3000 Pump	Pump only for Protekt™ Aire 3000
80032	Protekt™ Aire 3000 Mattress	Mattress only for Protekt™ Aire 3000
80033	Protekt™ Aire 3000 Cover	Cover only for Protekt™ Aire 3000
80034	Protekt™ Aire 3000 Bladder	Individual Bladder for Protekt™ Aire 3000



270 Washington Street Mount Vernon, NY 10553

**Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005**

[www.proactivemedical.com](http://www.proactivemedical.com) • [sales@proactivemedical.com](mailto:sales@proactivemedical.com)